



Table of Contents

Everything you need for a well stocked laboratory
or educational environment



4	Laboratory Equipment	68	Tube Racks & Storage
22	Workstation Essentials	87	Cold Products & Storage
33	Sample Collection & Transport	103	Slide Products & Accessories
41	Liquid Sample Handling	110	Safety Products & Dispensers
53	Sample Handling	119	Sealing & Labeling
65	Tubes & Vials	127	Charts & References

Laboratory Equipment

High-quality equipment and instruments for a vast range of protocols



Centrifuges & Stirrers

Mixers & Shakers

Sampling & Testing Instruments

It is our objective to exceed customer expectations in service and product quality

Our strict standards include
3rd party plant reviews
and **100% multi-level**
product inspections.

These actions result in a
99.7% reliability on
our equipment and
99.99% reliability on
our laboratory supplies.

We offer an industry leading
FIVE-YEAR limited warranty
for our Benchtop Laboratory
Equipment Product line.



MagFuge®

A high-speed mini centrifuge and magnetic stirrer in one unit



Holds
5.0 mL
Tubes



High-Speed Mini Centrifuge

Two Units
in One!

Magnetic Stirrer

High-speed Centrifuge is easy to use and designed to deliver maximum speed (12,500rpm/9,783x g) and flexibility with rapid acceleration and fast controlled breaking

Ideal for most protocols requiring fast spins including PCR, micro- filter cell separation and HPLC. Versatility to handle 0.2, 0.5, 1.5, 2.0 and 5.0 mL tube sizes

Intuitive One-Dial Programming

- Set speed and time with ease:
- Adjustable speed from 500 to 12,500 rpm
- Speed setting for RPM or RCF
- Adjustable time from 30 seconds to 60 minutes

Safe to Use

- Safety features help ensure it will not operate when the lid is open, a rotor is jammed, or the unit is out of balance

High-performance Magnetic Stirrer is designed using high quality rare earth magnets creating a robust magnetic coupling that significantly reduce the chance of spin-outs

Perfect for work with sensitive samples and protocols up to 3L requiring aggressive mixing for full dilution

Easy to Use

- Set stir mode, speed and time with easy one-dial programming:
- Select Mode: clockwise, c. clockwise and oscillate
- Adjustable speed from 50 to 2,500 rpm
- Continuous and adjustable time settings from 30 seconds to 60 minutes

Durable Accessories Included

- Re-usable, easy to clean silicone mats included to keep vessel stable and prevent spills and two NdFeB stir bars (40 x 12 mm) with chemical resistant PTFE coating included

Item No.	Color	L x W x H in	L x W x H cm
120581	Grey/Purple	9.2 x 7.3 x 4.7	23.5 x 18.5 x 12
120582	Grey/Blue		

Replacement parts

120622	MagFuge® Rotor, 12 x 1.5/2.0 mL tube
120623	MagFuge® Rotor, 6 x 5.0 mL tube
120624	MagFuge® Stirrer Rotor
120682	MagFuge® Silicone Mat
120527	Power Adapter
120015	Tube Adapters, 0.2 mL, 0.5 mL, Pack of 12 (6 ea.)
HS100502	Universal Power Cord Set

Patent Information

Community Design No. 003515337-0001-0002


Euro. Pat. App. 3 246 088

US Design Patent No. D814,045

US Patent No. 10471439



Package Includes		
Centrifuge/Stirrer	Stirrer Rotor	12 Tube Adapters 0.2 mL & 0.5 mL
5.0 mL Tube Rotor	2 Stir Bars, 35 mm x 12 mm, cylindrical shape	
1.5/2.0 mL Tube Rotor	2 Silicone Mats for Stirrer Lid	
1 Low Voltage, Double Insulated Power Adapter (100-240 VAC, 50/60 Hz) with 4 Interchangeable Plugs		



Specifications: Centrifuge

Speed Range	500 – 999 rpm (increments of 1); 1,000 – 5,000 rpm (increments of 5); 5,000 – 12,500 rpm (increments of 25)
RCF	16 – 9,783 x g (max RCF calculated)
Run Time Ranges	30 seconds – 60 minutes (5 second increments)

Specifications: Magnetic Stirrer

Speed Range	50 – 999 rpm (increments of 1); 1,000 – 2,500 rpm (increments of 5)
Run Time Range	Stirrer Mode 30 seconds to 60 minutes or continuous (5 second increments)
	Oscillating Stirrer Mode 1 second – 60 minutes or continuous
	Oscillating Period 30 seconds – 15 minutes
Capacity	3 L plus, Diameter of vessel not to exceed 130 mm



Gusto® High-Speed Mini Centrifuge

Maximum speed with minimal space

This compact centrifuge is ideal for protocols calling for fast spins with a higher RCF requirement

Made of durable ABS with UV-Resistant coating

User Friendly

- Intuitive Push button programming and bright LCD screen make it easy to use

Quiet Operation

- Brush-less motor enables quiet operation and low vibration

Safety Features

- Lid safety feature prevents opening while rotor is turning and clear polycarbonate allows user to see rotor movement for added visual safety
- Includes automatic shut-off feature for sample imbalance and tilt

Item No.	Color	L x W x H in	L x W x H cm
HS10050	Blue ■	8 x 6.7 x 4.5	20.3 x 17.1 x 11.4

Replacement parts

HS100502	Universal Power Cord Set
HS100503	Gusto® Rotor 12 x 1.5/2.0 mL Tube (with Cover and Knob)
HS100504	Gusto® Rotor 4 x 8 Strip 0.2 mL Tube (with Knob)
120015	Tube Adapters, 0.2, 0.5 mL, Pack of 12 (6 ea.)
HS100506	Gusto® Rotor Cover (for 12 x 1.5/2.0 mL Tube Rotor Only)
120527	Power Adapter


Specifications

Speed Range	500–12,500 rpm
Run Time Ranges	1 to 30 minutes
RCF	15 - 9,783 x g (max RCF calculated)



Package Includes

- Gusto® High-Speed Mini Centrifuge
- 12-Place 1.5/2.0 mL Tube Rotor
- 32-Place 0.2 mL Tube Rotor (singles, or 4 x 8 Strip)
- 12 Tube Adapters 0.2 mL
- 12 Tube Adapters 0.5 mL
- 1 Low Voltage, Double Insulated Power Adapter (100-240 VAC, 50/60 Hz) with 4 Interchangeable Plugs



The lid can be unlocked and opened with one hand

5
Year
Warranty

EASY SWITCH OVER
NO
TOOLS REQUIRED

sprout
plus

Sprout® Plus Mini Centrifuge

Mini but Mighty - ideal for small runs

Easy to Use

- Tool-free SnapSpin rotors and quick-spin options make the centrifuges fast and easy to use

Ergonomic Design

- The lid was functionally designed to fit in the palm of the hand for one handed lock and unlock allowing for quick tube extraction
- Large lid size accommodates longer and skirted tubes and tall tube lids
- The housing is made from ABS plastic and lid is polycarbonate plastic enabling the unit to be easily cleaned

Safety Measures

- Fitted with a braking function, the unit quickly stops spinning when powered off or when lid is opened
- A weighted, rubberized base eliminates potential movement and dampens vibration and noise



Accessory Kit Included



Specifications

Speed Range	Approx. 6,000 rpm (fully loaded)
RCF	2,000 x g



Sprout® Plus

Item No.	Color	L x W x H in	L x W x H cm
120610	Green	6 x 5 x 4.3	15.3 x 12.8 x 10.9
120611	Blue		

Replacement parts

HS100011	Sprout® Plus Rotor SnapSpin 6 x 1.5/2.0 mL Tube
HS100012	Sprout® Plus Rotor SnapSpin 2 x 8 Strip, 0.2 mL
120148	Sprout® Plus Lid O-Ring, Pack of 10
120491	Sprout® Plus Power Adapter
120015	Tube Adapters, 0.2 mL, 0.5 mL, Pack of 12 (6 ea.)

Package Includes

Sprout® Plus Mini Centrifuge
6-Place 1.5/2.0 mL Tube Rotor
16-Place 0.2 mL Tube Rotor (singles, or 2 x 8 strip)
6 Tube Adapters 0.2 mL
6 Tube Adapters 0.5 mL
4 Lid O-Rings
Storage Case for Rotors and Adapters
Tube Rack
1 Low Voltage, Double Insulated Power Adapter (100-240 VAC, 50/60 Hz) with 4 Interchangeable Plugs



Need tubes?

See our portfolio of 0.2 mL, 0.5 mL, 1.5 mL Microcentrifuge Tubes





Specially Designed Inserts
seat the centrifuge, all its accessories,
and includes a tube rack

12-Volt Car Adapter
for on-the-go testing

Portable Centrifuge Kit

Spin samples *on-the-go*

The Portable Centrifuge Kit offers a convenient way to spin down and prepare samples outside of a conventional laboratory environment

Immediate Field Processing

- Process samples in remote areas before deterioration can take place
- Allows for immediate on-site testing for faster results
- Made of durable ABS plastic for easy clean in field environments

Easy to Use

- Centrifuge is equipped with tool-less SnapSpin rotors for easy change out
- Simple operation with lock/unlock activation

Portable Convenience

- All items in kit are housed in a handled DuraPorter® Transport Box for portability
- The handled box has an air-tight seal for optimal sample protection

*See page 36 to learn more about the
DuraPorter® Transport Box*

Package Includes

Sprout® Plus Mini-Centrifuge

- 6-Place 1.5/2.0 mL Tube Rotor
- 16-Place 0.2 mL Tube Rotor (singles, or 2 x 8 strip)
- 6 Tube Adapters 0.2 mL
- 6 Tube Adapters 0.5 mL

4 Lid O-Rings

Storage Case for Rotors and Adapters

Tube Rack

- 1 Low Voltage, Double Insulated Power Adapter (100-240 VAC, 50/60 Hz) with 4 Interchangeable Plugs



DuraPorter® Transport Box

- Tube Rack Insert (12-Place 1.5/2.0 mL, 14-Place 0.2 mL)

Transport Box Dividers

Sprout® Plus Cradle

Portable 12V Sprout® Plus Centrifuge Car Adapter

Item No.

120568	Portable Centrifugation Kit, 100-240 VAC, 50/60Hz Universal Plug + 12V Car Adapter
120569	Portable Centrifuge Transport Case with Inserts - Only
120570	12V Sprout® Plus Centrifuge Car Adapter*

*Adapter can be used with Sprout® Plus Mini Centrifuge

See page 8 for functional details on the Sprout® Plus

Magnetic Induction Stirrer

Strong, motor-less and maintenance free

Superior Functional Design

- Advanced magnetic coil technology drives the spinning of the stir bar maintaining continuous, constant speed for reproducible results
- Magnetic coupling is optimized by using alternating current (AC) as its driving force strong enough to stir viscous solutions
- Motor-less operation generates minimal heat ideal for stirring temperature sensitive samples and environments for long periods
- With no moving parts, produces less vibration and quieter operation over traditional magnetic stirrers

Intuitive Operation

- Multiple stirring modes are easy to program
- Bright digital display precisely controls time and speed
- Quickly meets stirring speeds with rapid acceleration and braking (less than 10 seconds)



Protective silicone anti-skid plate cover for improved safety



Durable ABS housing is corrosion resistant and easy to clean



Specifications

Speed Range	50–2,000 rpm variable, 1 rpm increments
Run Time Range	30 seconds to 60 minutes, variable, 1 second increments or continuous
Stirrer Modes	Single clockwise Single counter clockwise Auto reverse rotation 1 to 99 cycles 50% duty
Capacity	3 L plus
Mixing Surface	165 mm diameter, silicone mat 160 mm diameter



Package Includes

- Rare Earth, NdFeB Magnetic Induction Stirrer
- 2 Stir Bars 35 mm x 12 mm Cylindrical Shape
- 1 Silicone Anti-Skid Plate Cover
- 1 Low Voltage, Double Insulated Power Adapter (100-240 VAC, 50/60 Hz) with 4 Interchangeable Plugs



Item No.	Color	L x W x H in	L x W x H cm
120584	Grey/Purple	9.7 x 7.7 x 3	24.5 x 19.5 x 7.5
120585	Grey/Blue		

Replacement parts

120738	MI Stirrer Silicone Mat
120739	MI Stirrer Universal Power Adapter



Mini Magnetic Stirrer

Affordable solution, ultimate function

The Mini Magnetic Stirrer is an affordable, dependable solution ideal for general research, academic, life sciences, industrial and clinical lab settings

Durable Design

- Constructed of chemical resistant ABS and PC blend
- Compact footprint saves valuable space on benchtop
- Lightweight for easy transport

Ultimate Function

- Powerful magnetic coupling ensures accurate alignment, reducing spin-outs during acceleration or deceleration
- Performs routine stirring applications with low vibration and quiet operation

U.S. and E.U. Utility Patent Pending
U.S. Design Patent No. D814,045

Item No.	Color	L x W x H in	L x W x H cm
120155	Grey/Purple	7.2 x 5.2 x 1.85	18.3 x 13.2 x 4.7
120595	Grey/Blue		

Replacement Part

120815	Mini Stirrer Universal Power Adapter
--------	--------------------------------------

Specifications

Speed Range 200–2,200 rpm (Infinitely adjustable)

Capacity Up to 3 L




Package Includes

Mini Magnetic Stirrer

2 Stir Bars 35 mm x 9 mm, Octagonal Shape

1 Low Voltage, Double Insulated Power Adapter (100-240 VAC, 50/60 Hz) with 4 Interchangeable Plugs





Stir Bars & Stir Bar Retriever

- Magnetic bar is made from ALNICO, one of the strongest permanent magnets available
- The sintered manufacturing process reduces the chance of decoupling (spin outs), encapsulation in PTFE protects both magnet and what it is stirring
- Chemically inert and easy to clean
- Can withstand temperatures up to 250°C, making them suitable for use with both unheated and heated stirrers
- Low friction octagon shape with a center pivot ring helps create a strong mixing vortex, even with more viscous solutions
- The retrieval bar can recover stir bars even in hot solutions when stirring and mixing is no longer required

10 x 5 mm


20 x 8 mm

25 x 8 mm

28 x 8 mm

30 x 8 mm

35 x 9 mm



Octagonal Shape

Stir Bars

Pack of 10

Item No.	L x W (in)	L x W (mm)
120547	0.4 x 0.2	10 x 5
120548	0.8 x 0.3	20 x 8
120549	1 x 0.3	25 x 8
120550	1.1 x 0.3	28 x 8
120551	1.2 x 0.3	30 x 8
120552	1.4 x 0.4	35 x 9

Stir Bar Retriever

Pack of 5

Item No.	L (in)	L (cm)
120553	11.8	30

Lab Mats

Two Lab Mat options to keep your workstation safe, quiet, clean and organized



Workstation Lab Mat

gridded surface for tube staging
& preventing movement of small
objects

Silicone Lab Mat

large noise dampening work surface
with lipped edge to contain spills



Vortexer

Multi-purpose, maximum convenience

Innovative Design

- One mixing head does it all with this sleek, modern and compact designed vortex mixer
- Hold slot configurations for multiple sizes of tubes, vessels and assay plates
- Center pad for quick hand-held mixing

Easy to Use

- Control switch with three operation selections (intermittent, off, continuous)
- Variable Speed Rotary Dial provides pre-set speed control selection during shaking mode

Durable Construction

- ABS material resist most commonly used chemicals and is easy to clean
- Stable elastomer feet dampen vibration and prevent movement during operation



Center pad area for hand held mixing



Assay plate x 1
(to be used with retention cords)



25 mL, 50 mL conical tubes x 2
(horizontally mounted)



15 mL conical tubes x 2
(horizontally mounted)



5 mL tubes x 2
(horizontally mounted in the same retention space as 15 mL)



0.5 mL microfuge tubes x 6
(vertically mounted)



1.5/2.0 mL microfuge tubes x 4
(vertically mounted)



0.2 mL microfuge tubes
(strip tubes) x 24
(vertically mounted)

US Design Patent No. D787,085

US Utility Patent No. 9895670

Item No.	Speed Control	Power Input	L x W x H in	L x W x H cm
120209	Infinity Dial	110-120 VAC, 50/60 Hz, 55 W, US/NA		
120212	Infinity Dial	220-240 VAC, 50 Hz, 44 W, CE		6.7 x 7.5 x 7.8
120214	Infinity Dial	220-240 VAC, 50 Hz, 44 W, UK		17.1 x 19 x 19.7
120210	Preset Dial	220-240 VAC, 50 Hz, 44 W, AUS		(with mixing head in place)
120318	Preset Dial	100-110 VAC, 50/60 Hz, 66 W(PSE), Japan		(with mixing head in place)

Replacement parts	
120327	Vortexer Multi-Head
120328	Vortexer Central Pad, 3 pk
120422	Vortexer Retention Cords, 3 pk
120329	100-120 VAC US Cord

Replacement parts	
120333	220-240 VAC AUS Cord
120332	220-240 VAC CE Cord
120331	220-240 VAC UK Cord
120330	100-110 VAC Japan Cord

Specifications:

Speed Range Approx. 1,000; 2,000; 3,000 rpm

Orbit 3.7 mm

Control 3 Position Control — (top) intermittent, (middle) off, (bottom) continuous



Mini Vortexer

Small footprint with maximum power


Simple, touch start operation offers exceptional mixing performance with one hand

Robust Design

- Low profile design and small footprint saves valuable space
- Made of ABS and silicone materials that resist most commonly used chemicals
- Stable platform and elastomer feet ensure unit stays stable during operation

Powerful Function

- Rapid vortex action with low activation force
- Quiet operation and low vibration



U.S. and E.U. Utility Patent Pending
U.S. Design Patent No. D814,045

Item No.	Color	L x W x H in	L x W x H cm
120567	Grey/Purple	5.1 x 4.3 x 3	13 x 11 x 7.5
120598	Grey/Blue		
120827	Black/Black		

Replacement parts

120617	Mini Vortexer Silicone Cover, set of 2
120814	Mini Vortexer Universal Power Adapter

Specifications

Speed	4,500 rpm
Orbit	6 mm
Capacity	Up to 50 mL tube



Package Includes

- Mini Vortexer
- 2 Silicone Caps
- 1 Low Voltage, Double Insulated Power Adapter (100-240 VAC, 50/60 Hz) with 4 Interchangeable Plugs





Digital Orbital Shaker

Shaker for a vast range of protocols

Versatile Design

- Create your own matrix for vessels of all shapes and sizes with Hi-Lo Platform and Bungee Cords
- Instantly adjusts to weight load changes and unbalanced loads with no-stop electronic motor feedback allowing for truly reproducible results
- Durable ABS construction suitable for use in incubators (excluding CO₂) and cold rooms within the stated temperature range

High-Featured Operation

- Suitable for wide range of protocols including very low start speed for sensitive sample agitation to high speeds for vigorous applications
- Soft start and shut down feature reduces spilling from full vessels
- Designed to handle heavy loads up to 10 lbs (4.5 kg)
- Overload protection shut-off prevents overloading



Item No.	L x W x H in	L x W x H cm	Platform Interior in	W x D cm
120460	12.4 x 7.8 x 12.2	31.4 x 19.7 x 31	11 x 11	27.9 x 27.9

Replacement parts

120521	Digital Orbital Shaker Silicone Mat & Bumpers
120524	Digital Orbital Shaker Bungee Cords, Pack of 8
120491	Digital Orbital Shaker Universal Power Adapter

Specifications

Speed Range	20–300 rpm, 1 rpm increments
Orbit Diameter	19 mm (0.8")
Run Time Range	0 to 74 hours and 59 minutes; 1 second increments or continuous
Maximum Load	4.5 kg (10 lbs)



Package Includes

Digital Orbital Shaker

Skid-Free Rubber Platform Mat

8 Bungee Cord Set

Universal Hi-Lo Platform with Soft Rubber Edge Bumpers

1 Low Voltage, Double Insulated Power Adapter (100-240 VAC, 50/60 Hz) with 4 Interchangeable Plugs





Battery and Filter not included in warranty



Soft-Coated Triggers
with direction indicators

Ergonomically Balanced
in the hand with a loaded pipette

Universal Grip
for right- or left- handed users



Adapter
for small pipettes

RF3000®

Packed with features for comfortable and accurate pipetting

The powerful pump offers rapid, yet sensitive aspiration and dispensing, has adjustable valves and is microprocessor controlled for constant reliability and accuracy

- Long-life li-ion battery, with quick charge—up to 20 hours use
- Toggle operation by pressing the large mode button
- An adapter can be fitted to the nose cone to stabilize small volume and thinner pipettes

- New USB connector ensures the appropriate charger is used with the unit
- Low-battery indicator
- Can be recharged when in use



Nose cone assembly only



Storage Stand
doubles as a charging station

Easy Identification

5 different colored nose cones help differentiate department, task or user; as well as a labeling area on the side of the unit



Magnetic & Screw Wall-Mount



3 Modes

Bright LEDs indicate 3 operation modes: high, low, & gravity



Specifications

Power	Micro-USB connector Rechargeable li-Ion battery (replaceable)
Filter	PTFE/25mm dia/0.45 µm filter
Pipettes	1–100 mL, glass and plastic
Pipette Nozzle	Polypropylene



Package Includes

- RF3000® Pipette Controller
- Charging Stand and Wall-Mountable/Magnetic Holder
- Nozzle Adapter for use with Small Pipettes
- Li-Ion Rechargeable Battery
- Five Nose Cones: Black, Yellow, Red, Grey and Blue
- 5VDC Micro USB Plug
- 6 ft (183 cm) Cable
- 1 Low Voltage, Double Insulated Power Adapter (100-240 VAC, 50/60 Hz) with 4 Interchangeable Plugs



Item No.	Volume Range mL	L x W x H in	L x W x H cm
HS3000	1–100	9 x 8.5 x 5	22.9 x 21.6 x 12.7

Replacement parts

HS5000RFFIL	RF3000® Replacement Filters, Pack of 5
120669	RF3000® Universal Power Adaptor 5VDC Micro USB Plug
HS300502	RF3000® Li-Ion Rechargeable Battery, One year warranty on battery included
120487	RF3000® Nose Cones Assorted* ■ ■ ■ ■ ■ Set of 5
HS5000RFRFA	RF3000® Replacement Filter Adapter
HS5000RFRP	RF3000® Replacement Pipette Silicone Adapter

*Assorted colors (Black, Blue, Grey, Red and Yellow)

Pipette Controller

Precise fluid control at your finger tips

Comfortable to hold and easy to use, the manual inline pipette controller was designed with an ergonomic grip and bulb for simple and precise operation

A T-bar thumb lever controls the aspirating and dispensing speeds from serological pipettes from 1 to 100 mL in size

- Conveniently accommodates either glass or plastic serological or volumetric pipettes ranging from 1–100 mL
- Includes a replaceable .45 µm hydrophobic filter to protect your samples and internal components



Nose cone assembly only

Item No.	Color	Volume Range mL	L x W x H in	L x W x H cm
120479	Blue	1–100	9 x 2.1 x 2.1	23 x 5.4 x 5.4

Replacement Parts

Item No.	Description	Color
120604	Three Noses Cone with Adapter	Green, White, Blue

Operation Is Simple:

Squeeze the bulb

to create a vacuum to draw up fluid

Push the control lever forward (down)

to dispense the liquid

Push the control lever backward (up)

to aspirate liquid through a pipette

Press the central bulb on the unit

to blow out any remaining fluid if necessary

Inline Ergonomic Design

reduces stress with ease of handling and operation

1
Year
Warranty

Filter not included
in warranty

Squeeze

to create a vacuum

Stay Organized

with three additional nose cones (green, white, blue) for color identification for both the Pipette Controller & RF1000



Squeeze
to create a vacuum

1
Year
Warranty

Filter not included
in warranty

RF1000

Lightweight and easy to use

The RF1000 manual pipette controller is uniquely designed for one-handed operation

Simply squeeze the main bulb and push the toggle forward to aspirate or back to dispense. Push the blowout bulb to expel any liquid remaining in the pipette

- Tapered silicone rubber pipette adapter allows a large range of pipette sizes ranging 1–100 mL and is replaceable
- Includes a replaceable .45 µm hydrophobic filter to protect your samples and internal components



Nose cone assembly only

Item No.	Color	Volume Range mL	L x W x H in	L x W x H cm
HS1000RFU	Blue	1–100	7.1 x 1.8 x 4.7	18 x 4.5 x 12

Replacement Parts

Item No.	Description	Color
HS5000RFFIL	Replacement Filters, Pack of 5	
120604	Three Noses Cone with Adapter	Green, White, Blue



Ergonomic

Slight downward slant to tip ejector button reduces hand fatigue, incidence of repetitive strain injury (RSI)

Volume Adjustment
and low-force dispensing
all in one button for
convenience

**Finger
hook**

**3-Position
volume display**

Volume Locking Ring
prevents accidental
volume changes

Pearl® Adjustable Volume Pipettes

Precision liquid handling

A complete range of single- and multi-channel pipette volumes to cover all your possible needs

The sleek, ergonomic design and reduced plunger pressure help reduce the chance of Repetitive Stress Injury (RSI), as the pipettes require 2/3 the operation pressure of most standard pipettes

Each pipette is calibrated and a calibration certificate is supplied. Manufactured in accordance with ISO8655 and GMP regulations

- Internal metal parts for increased accuracy and calibration stability
- Separate tip ejector with the exception of the 500 µL (120105)
- Large digit volume display



Color Coded & Labeled
with volume range (Both top & side
colored components for quick visual
confirmation)

Precision Metal Components
for both shaft and barrel ensure reliable
performance and long life



No. of Channels	Item No.	Color	Volume Range µL
Single	120097	Green	0.1–2 µL
	120098	Lt. Blue	0.5–10 µL
	120099	Lt. Purple	2–20 µL
	120100	Orange	5–50 µL
	120101	Blue	10–100 µL
	120102	Yellow	20–200 µL
	120103	Purple	30–300 µL
	120104	Blue	100–1,000 µL
	120105	Pink	500–5,000 µL
	120106	Red	1,000–10,000 µL
Eight	120215	Orange	5–50 µL
	120216	Purple	30–300 µL
Twelve	120217	Orange	5–50 µL
	120218	Purple	30–300 µL

Pipette Starter Kit

Universal Carousel Pipette Stand

Hold and store pipettes in style with the Universal Carousel Pipette Stand by Heathrow Scientific that holds up to 8 single-channel, 4 multi-channel, or a combination of both pipettes

A sleek and modern design enables this stand to securely hold virtually every known brand of pipette in the marketplace making it truly 'universal' while adding a pop of color to a bench top

Reagent Reservoirs

Multi-channel reservoirs by Heathrow Scientific are superior in every way with thicker walls, crush-resistant "turrets" at the corners, easy pour-spouts, 25 mL graduation marks, and tabs to prevent stacked reservoirs from sticking together

Pearl® Pipettes

Pearl Pipettes by Heathrow Scientific were specifically engineered for ease of use and comfort

The sleek, ergonomic design and reduced plunger pressure help reduce the chance of Repetitive Stress Injury (RSI), as the pipettes require 2/3 the operation pressure of most standard pipettes



Each pipette calibrated and a calibration certificate is supplied

Manufactured in accordance with ISO8655 and GMP regulations






Pipette Starter Kit Item No. 120554

Universal Carousel Pipette Stand For more information see page 45


Upper Holder Color	Lower Holder Color	Base Width in	Top W x H in	Base Width cm	Top W x H cm
Assorted/Clear 	Clear 	6.2	9 x 16.1	15.8	22.9 x 40.8

*Assorted colors (Purple, Blue, Green, Grey, Red Orange, Yellow and Clear)

Pearl® Pipettes - Single Channel For more information see page 18

Color	Volume Range μL
Lt. Purple 	2–20 μL
Yellow 	20–200 μL
Blue 	100–1,000 μL

Reagent Reservoirs For more information see page 51

Color	Volume mL	Material	Sterility	L x W x H in	L x W x H cm
White 	50	Polystyrene	Sterile	4 x 1.1 x 1.2	10.2 x 2.8 x 3

*All parts available separately



eCount

Rapid, accurate colony counter

Combines a multi-function electronic counter with a Sharpie® Pen for marking to prevent missing or double-counting colonies

- Light, ergonomic body is balanced for easy manipulation
- Select between an audible beep or an LED flash to acknowledge count
- Store up to 32 separate counts and display the total of all counts
- Count up/down verification
- Quick and easy to change between different colored markers to check different colony counts in the same view

Package Includes:

- eCount Colony Counter
- Black Sharpie® Pen (Page 26)
- Long-life button battery included
- Hand tally Adapter
- Weighted Stand



US Patent No. D564,385
EP Registered Design No. 705751-0001

Item No.	L x W in	L x W cm
120000	7.4 x 1.2*	17 x 3*
120228	CR1225 Battery	

*Measurements shown are with Sharpie® Pen inserted



Battery and Pen not included in warranty



Lock Counter Off

to allow use of the marker without registering more counts

Hand Tally Adapter

allows you to use the eCount as a digital tally counter



Weighted Stand
for bench top use





Mini-Light Box and Mini-Magnifier

Perfect accessories for the eCount or use as stand-alone items



Mini-Light Box

- LED Light source. Color temperature 5,300K +/- 270K
- Viewing area (125 x 100 mm) is perfect for quick screening of 96-well plates, petri dishes, electrophoresis gels, slides, and transparencies
- Powered by six 1.5V AAA batteries (included), not included in warranty

Mini-Magnifier

- Swiveling molded lens 1.75x magnification
- Enlarges objects for easier viewing

Item No.	Description	L x W x H in	L x W x H cm
120661	Mini-Light Box*	6.2 x 5.9 x 0.6	15.8 x 15 x 1.4
HS15092D	Mini-Magnifier	6.3 x 5.3 x 8.1	16.1 x 13.5 x 20.7

*Batteries included



Hand Tally Counter

Use for counting blood cells, bacterial colonies, drops of liquid or any repetitive event

- Tallies up to 9999
- Features a quick-reset knob
- Lens for easy reading, and finger ring for a secure grip

Item No.	Dia. x W in	Dia. x W cm
HS6594	1.8 x 1.6	4.6 x 4.1



Benchtop Counters

Mechanical, benchtop counters are perfect for blood cells (label included), bacterial colonies or repetitive counting

- Each key records up to 999
- Blood cell key labels
- Total window keeps track of total number of strokes made on all other keys; bell rings when total reaches 100

Item No.	No. of Keys	L x W x H in	L x W x H cm
HS5360A	6-Key	9.8 x 3 x 2.2	25 x 7.6 x 5.6
HS5360B	9-Key	12.5 x 3 x 2.2	32 x 7.6 x 5.6

Workstation Essentials

Products necessary for workflow efficiency and keeping your work space clean and organized



Lab Mats

Notebooks & Timers

Drying Racks

Organization

Curiosity drives us beyond yesterday's way of thinking:

We are an inquisitive bunch focused on understanding the lab science world and how to make a difference.

We change, we adapt and we expect creative thinking and solutions from everyone we align with.

We don't just make lab products, we make the work lab professionals do easier, more efficient and more fun.



Reusable Lab Mats

Sustainable, Safe Solutions

Reusable lab mats are an essential lab solution to help keep benchtops organized and clean

They help facilitate workflow processes and are ideal for every type of laboratory

Save Money & Eliminate Waste

- Since the lab mats are reusable, they reduce cost and waste making it a more environmentally friendly option over the alternative paper mats

Adds Stability & Safety

- These mats offer a large working surface that is anti-skid which keeps delicate glassware and expensive equipment in place
- The mat's surface is also noise dampening, reducing sound from the vibration of common benchtop equipment
- Made from durable food grade silicone, they are chemically resistance and have a raised edge to contain chemical spills

Easy to Clean & Store

- When the mats are done being used, they can be easily cleaned using any standard disinfectant by simply wiping it, or putting it in a dishwasher and it's ready for repeated use



Workstation Lab Mat

Safely facilitate workflow processes

- Protect counter tops from hot items up to 200°C (392°F)
- Lip around the edge can hold approx. 250 mL of liquid
- Small footprint, 12 x 14.9 x 0.1 in (30.5 x 37.9 x 0.2 cm)



Gridded Surface
for tube staging, preventing movement of objects

Reversible
Side One: combination gridded and smooth surface
Side Two: smooth surface only



Item No.	Color	L x W x H in	L x W x H cm
120747	Blue	12 x 14.9 x 0.1	30.5 x 37.9 x 0.2
120748	Purple		

Silicone Lab Mat

Large working space

The lab mat can be used with items such as centrifuges, vortex mixers, hot plates, stirrers, pipettes and so much more

- Mat features designated magnetic area to keep stir bars in place
- Protect countertops from hot items up to 200°C (392°F)



Raised Edges
enable liquid retention of approximately 300 mL of fluid

Large Working Surface
accommodates benchtop equipment, weigh boats, stir bars, scales, etc.



Item No.	Side #1 Color	Side #2 Color	L x W x H in	L x W x H cm
120506	Yellow	Grey/Blue	14 x 23.5 x 0.1	35 x 60 x 0.2
120507	Purple	Purple/Grey		
120599	Blue	White/Blue		

Laboratory Notebooks

Keep information safe and well organized with Laboratory Notebooks

Record research data, experimental procedures, design and development ideas and more

- Quality Smyth Sewn binding holds pages securely and lies flat when opened
- Made of acid-free paper, pages are numbered and include space for title, date, project number, book number, date and witness signature
- User data, notebook instructions, table of contents, periodic table, and conversion tables are included in every book

Keep Information Safe
with a durable hardbound, water-resistant cover



Marker Ribbon
to keep your place

Easy Identification
with a foil stamped title

Item No.	Color	Pages	Page Format.	L x W <i>in</i>	L x W <i>cm</i>
HS8610A	Black	100	Grid	8.6 x 11.1	21.8 x 28.1
HS8610B	Blue	200			
HS8610D	Red	200	Lined		
HS8610E	Black	100			



2 Page Formats

- Grid: 6 x 6 mm grid size
- Lined: college-ruled pages for easy report writing



Twin Tip & Fine Tip Sharpie® Permanent Pens

Ink dries quickly and resists both fading and water

- Pen has a fine and an ultra fine tip and comes in a bright selection of colors
- Ink dries quickly and resists both fading and water
- AP certified

Item No.	Description	Color
HS15094	Twin Tip Marker Set, Pack of 8	Assorted
HS15092ZP	Fine Point Marker, Pack of 12	Black

*Assorted colors (Black, Turquoise, Navy, Lime, Orange, Purple, Sky Blue, Berry)
Sharpie® is a registered trademark of Sanford, A Newell Rubbermaid Company

Lab Alert® Clocks & Timers

a



Pocket Timer/Stopwatch/Clock

- Four independent timing channels with visual and audible alarms
- Use as a stopwatch, a count up/down timer, or a clock
- Includes spring clip and magnetic back

b



Triple Timer/Clock

- Three channels to count up/down with
- Simultaneous display, time-out, and 60 second audible alarm
- Channels may be started independently or simultaneously
- Timer automatically counts up after counting down to zero
- Memory recall

c




4-way Timer with Whiteboard/Pen

- Clock feature count up or down to 24 hours 59 minutes
- 4 timer channels
- Stands up or attaches to a magnetic surface (magnets on back of unit)
- Whiteboard and marking pen included
- Sounds when time expires

	Item No.	No. of Channels	Timing Capacity	Count Direction	Display Digits, H mm	Battery Included	L x W x H in	L x W x H cm
a	HS24490D	Four	100 hr.	Up/Down	6	One LR44 1.5V	2.4 x 2.8 x 1	6.1 x 7.1 x 2.5
b	HS24490M	Three	20 hr.	Up/Down	10.5	One AAA	3.4 x 2.6 x 0.9	8.6 x 6.6 x 2.2
c	120365	Four	24 hr.	Up/Down	10	Two AAA	4.8 x 3.4 x 0.9	12.3 x 8.6 x 2.1
d	HS24645	One	24 hr.	Up	8.5	One LR44 1.5V	3.4 x 3.1 x 0.8	8.7 x 7.8 x 2.1
e	HS24670	Four	99 hr.	Down	14	Two LR44 1.5V	2.9 x 2.3 x 0.5	7.3 x 5.9 x 1.3

d



Multi-function Stopwatch/Clock

- Quick start, quick stop timer features
- time, split, hourly chime, calendar, alarm
- 100 cm neck cord

e



Four Channel Timer/Clock

- Independently or simultaneously set time on four channels
- Audible and visual alarms signify elapsed time
- Easily toggle between channels and clock with top buttons
- Includes spring clip with magnetic back

Lab Drying Rack

Bringing color to a mundane task

Most Versatile Rack in the market


- Channels designed for more useful configurations
- Large and Small pegs for all glassware varieties
- One of a kind Drying Basket to collect smaller items

Designed to protect against breakage

- Peg angle allows for better hold - especially large, heavy pieces
- Patent Pending SureGrip Hook securely grips glassware into place

Included:

- 72 Pegs ~ 36 Small (Green) ~ 36 Large (Blue)
- 2 Drying Baskets
- Drainage Tubing
- Mounting Hardware



Drying Basket

to keep small items and slotted
to hold pipets safe while drying



SureGrip Hook

Patent-pending design
securely grips glassware
into place




Use on counter or wall

Stands alone on integrated
hinged stand. Mounts easily on
wall with included guide and
hardware

Can hold up to 77 lbs

(34.9 kg)

Rack made of durable
polystyrene



Lab Drying Rack

Item No.	L x W x H in	L x W x H cm
120791	17.4 x 26 x 5.5	44.1 x 66 x 13.9

Replacement Parts

Item No.	Description
120793	Peg 10 Piece Kit, (5) Small, (5) Large
120792	Drying Basket Kit, 2 ea
120813	Replacement Draining Set, Tube and Elbow

Wire Drying Rack

Saves glassware from chips and cracks

Drain and dry test tubes, beakers, flasks, and other labware

- Place self-standing rack on a table or near a sink, or use hooks to hang rack on the wall
- Holds up to 90 pieces of labware on 18 x 127 mm pins and 72 x 76 mm pins
- Constructed of spot-welded, HDPE-coated steel wire for durability with polystyrene tray and back
- Easily accessible and removable bottom drainage pan



Removable Drip Pan
helps retain and dispose of excess water

Item No.	L x W x H in	L x W x H cm
HS23243A	18.2 x 7.4 x 20.7	46.2 x 18.2 x 52.5



Workstation Lab Mat

Keep glassware in place while drying



Workstation Lab Mat, Page 25

Peg Racks

Perfect for drying small items

Fiberglass-reinforced polypropylene racks are ideal for bench work or tube storage

- Use in sink or on countertop
- Sturdy pegs hold items in place



Open tear drop wells allow for water to drain

Item No.	Color		Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS24311A	White	<input type="checkbox"/>	10	50	14–17	7.4 x 4.1 x 2.7	18.8 x 10.4 x 6.9
HS24311B	Blue	<input checked="" type="checkbox"/>					
HS24312A	White	<input type="checkbox"/>	7	96	10–13	3.2 x 5 x 2.7	8.2 x 12.6 x 6.9
HS24312B	Blue	<input checked="" type="checkbox"/>					

Pack of 2



MagLab

Magnetic lab storage series

Easily organize common lab products off your bench top with the MagLab storage system

Vibrant & Functional

- Made of robust ABS plastic in bright translucent colors
- Secure direct to metal with built-in powerful magnets or non-metal surfaces with included self-adhesive mounting plate
- Available individually or as a pack of four

MagRack holds four 50 mL tubes and nine 15 mL tubes

MagWipe holds a standard (110 x 210 mm) Kimwipes® box





MagPette holds two pipettors (fits most popular styles)

MagBox holds miscellaneous accessories



Magnetic

Includes self-adhesive mounting plates that adhere to any surface

Item No.	Description	Color	L x W x H in	L x W x H cm
HS23500	MagLab Series, Pack of 4	Assorted 	—	—
HS23502	MagWipe	Green 	5.1 x 3.5 x 3.8	13 x 9.1 x 9.8
HS23503	MagPette	Purple 	5.1 x 4.4 x 2.9	13.1 x 11.2 x 7.6
HS23504	MagBox	Orange 	5.1 x 3 x 4	13 x 7.8 x 10.2

*Assorted colors: Teal, Green, Purple and Orange

Workspace Organization

Studies* show that by clearing clutter from your work environment, you are better able to focus and process information, resulting in increased productivity.

*Princeton University Neuroscience Institute



Workstation Organizer

Keeps supplies organized and easily accessible

For easy access to items on counter, shelf or inside drawer

- Protect fragile items and contain small items

Item No.	Color		L x W x H in	L x W x H cm
HS2538A	White	<input type="checkbox"/>	17.9 x 12.1 x 2.5	45.5 x 30.7 x 6.3

Compartments	L x W in	L x W cm
2 Small Compartments	17 x 3.1	43 x 8
1 Large Compartment	17 x 4.3	43 x 11



Workstation Dispensing Bins

Clear acrylic bins let you see your small items inventory

- See inventory levels of pipette tips, tubes, fittings, finger cots, and other bulk small items at a glance
- Flip-top covers on generous openings make filling and dispensing easy and keep products dust-free
- 4 sizes available to cater to a wide variety of storage needs
- Includes non-skid rubber feet



Item No.	Color	Compartments	L x W x H in	L x W x H cm
HS23401	Clear <input type="checkbox"/>	1	7.2 x 7.6 x 9.1	18.4 x 19.4 x 23.2
HS23402		2	5.9 x 6.9 x 14.4	15 x 17.7 x 36.8
HS23403		3	8.5 x 6.6 x 8.5	21.7 x 17 x 21.7
HS23404		4	11.3 x 6.6 x 8.5	28.7 x 17 x 21.7

Tubby® Storage Container

Sturdy storage tub with lid and dividers

Our handy Tubby® stackable storage container is designed to help organize a laboratory

- The removable dividers can create as big or as small a compartment as needed
- Perfect for keeping tubes, tips, pipettes, gloves, and other small items organized and out of the way
- The lid snaps in place to secure the contents and for safe transport

Pack of 5

Item No.	L x W x H in	L x W x H cm
HS23453	13 x 7.9 x 4.5	33 x 20 x 11.5



Constructed of high-impact polypropylene

More Dispensers for organization and to keep safety PPE products at the ready



Workstation Storage Bin, Page 116



Pipette Workstation, Page 48



Modular Glove Box Holders, Page 117

Test Tube Basket

Ideal for washing and storing test tubes

Strong basket for use with mechanical washers and autoclaves

- Basket for use in washers, autoclaves or for storage of lab consumables
- Made from aluminum



Item No.	L x W x H in	L x W x H cm
HS20341B	10 x 6 x 6	25.6 x 15.2 x 15.2



Sample Collection & Transport

Make everyday easier and safer with products for sample and blood collection



Portable Trays
Totes & Transports
Clinical Tube Racks
Sample Trays & Racks



Ratcheting Handle
folds out of the way for easy access,
even over tall sharps containers

Tray Inserts
come in two configurations and
include an extra-deep compartment
for a sharps container



Disposable Insert A
for 13 mm tube rack



Disposable Insert B
for 16 mm tube rack

Droplet Blood Collection Tray

Designed with phlebotomists in mind

Keeping a clean and organized tray has never been easier. Wash the ABS tray as needed, and replace the insert periodically for a fresh start

- Multiple compartments keep a variety of supplies neatly at hand, and organized
- Built in side pockets offer additional storage space
- Handle includes soft rubber grip to reduce hand and finger strain
- Each tray includes either a 13 or 16 mm Heathrow OneRack® tube rack (page 78) and two disposable polystyrene inserts

Droplet Configurations

US Patent No. D544703

Item No.	Tray Insert	Tube Rack mm	L x W x H in	L x W x H cm	Handle Height in	Handle Height cm
HS2200A	HS2201A	13	16.8 x 10.4 x 6.7	42.9 x 26.5 x 17.2	11	28
HS2200B	HS2201B	16				

Droplet Inserts

Pack of 4

Item No.	Description	L x W in	L x W cm
HS2201A	Tray insert for 13 mm tube rack	14.2 x 8.2	36 x 21
HS2201B	Tray insert for 16 mm tube rack		



Keeps contents clean and secure

Droplet Lid

Cover for the Droplet blood collection tray

Protect the contents of the Droplet Tray with this Droplet lid. Custom designed sturdy lid is a snug fit, protects contents and provides access to the tray handle in its vertical position

- Made of polystyrene
- Can be secured with cable ties for extra security
- Stackable



Item No.	Color	L x W x H in	L x W x H cm
HS2202	White	16.5 x 9.4 x 8	42 x 24 x 20.5


X-tourn® Tourniquets

Hi-viz orange alert for safer practice

High-quality disposable tourniquets eliminate roping or curling when used

- Made from supple and smooth polyisoprene
- Latex free
- Less abrasive than other brands
- No rubber proteins, resins or fatty acids that can produce hypersensitive reactions
- Available in bulk or on a perforated roll

Item No.	Color	Packaging	L x W in	L x W cm
120108	Green	100/Bulk	18 x 1	45.7 x 2.5
120110	Orange	250/Bulk	18 x 1	45.7 x 2.5
120111		25 Rolls		



Uniform Thickness
prevents shredding and
breaking



Blood Collection Tray

Light weight, heavy on organization

Organize your phlebotomy, I.V. start or catheter insertion supplies with this spacious, robust and easy to clean carry tray

- Lightweight easy-to-clean ABS construction
- Customize the compartments to suit your requirements with the dividers provided
- Comfort grip on durable acetal plastic handle tested to hold up to 10 kg

Item No.	L x W x H in	L x W x H cm	Handle Height in	Handle Height cm
120261	16.4 x 12.2 x 3.3	41.9 x 31.1 x 8.5	9.4	24.1



For safe disposal of syringes, needles and other small sharps see our Sharps Containers

1 L Capacity
with a clear lid to view fill level



Compatible with
the Droplet blood collection tray &
Totemax® blood collection tray

Sharps Containers, Page 114

Totemax® Blood Collection Tray

Keep off the floor, work top or bed

Unique phlebotomy tray design and patented hook feature allows the tray to be on or over most bed rails, foot-boards and handrails found in a typical patient room

- Made from food grade HDPE resins for a durable non-breakable design
- Deep tray design helps reduce spills




Item No.	Color	L x W x H in	L x W x H cm
120021	Blue	17 x 11 x 11	43.2 x 27.9 x 27.9

Hooks Lock

in place when not in use

Hangs on bed rails,
footboards and handrails



Carry handle folds

neatly into body space allowing the units to be stacked

Molded Slots

to hold tubes while working

**Silicone Water-tight Seal**

and three lid clasps provide secure closure

**Lid Opens a Full 180°**

to expedite loading and emptying

**DuraPorter®****Specimen or Sample Transport Box**

Ideal for transporting biological and clinical samples, as well as instruments and products that you want to keep clean and dry under testing environmental conditions

The DuraPorter® transparent design is OSHA 29 CFR 1910.1031 and UN3373 compliant

- Polycarbonate construction and available in three colors for user/owner identification
- Accommodates 13 mm, 16 mm and 20 mm tubes with our Half and Full Size OneRack®
- Use the separators provided to create your own compartment sizes



Item No.	Color		L x W x H in	L x W x H cm
120052	Clear with Blue Handle	■	15 x 7.7 x 6.3	38 x 19.6 x 16
120077	Red with Red Handle	■		
120078	Yellow with Yellow Handle	■		

Load



Lock



Transport



Store

can be placed in a freezer or refrigerator to keep samples chilled






Extra Security Feature
Accommodates lockout tags
(not included)



A Perfect Fit
OneRack® Tube Rack

Accommodates 13 mm, 16 mm and 20 mm tubes with our Half and Full Size OneRack®
See page 78



For field sample processing, learn more about our Portable Centrifuge Kit

Includes Mini-Centrifuge and 12V Car Adapter for processing outside the laboratory




Clinical 4-Way Tube Rack

The only tube rack specifically designed to hold round, flat & conical blood collection and culture tubes

- Accommodates all brands of 13 and 17 mm round, flat and conical bottom blood collection & culture tubes
- Customize the rack to fit the number of wells required, infinite units can be linked together or use one standalone unit
- Molded alphanumeric identification for quick location and retrieval of tubes
- Robust, made from polypropylene that is chemical resistant and easy to clean
- Racks can be used in water baths and placed in freezers
- Patent Pending

Customizable

Racks can slide together to create many combinations



Holds Multiple Tube Sizes without multiple racks!

US Patent No. D900337


Pack of 5

Item No.	Color	Item No.	Color
120684	Blue	120687	Natural
120685	Purple	120683	Assorted*
120686	Green		

*Assorted colors (Blue, Purple, Green, Red, Natural)

Individual Rack Specifications

Side	Wells	Well Dia. mm	L x W x H in	L x W x H cm
A	12	13	9.1 x 2.8 x 3.8	23.2 x 7.0 x 9.5
B	18	13		
C	12	17		
D	18	17		




Mega Rack®

High-capacity storage for your tubes

Fully assembled, polypropylene mega racks are great for high throughput labs requiring large sample storage

- Sturdy, material for easy and secure transport of blood collection and other tubes
- Alphanumeric identification for quick location and retrieval of tubes
- Place different test groupings into one rack eliminating the need for multiple racks



Mega Racks® 10-13 mm - Holds 5-7 mL sample tubes

Item No.	Color	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS24322A	Blue	216	13	13.3 x 9.3 x 3.6	33.7 x 23.5 x 9.1
HS24322E	White				
120112	White	432	13	26.5 x 9.3 x 3.5	67.3 x 23.5 x 9.0
120113	Blue				

Mega Racks® 13-16 mm - Holds 10 mL sample tubes

Item No.	Color	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS24323A	Blue	120	16	13.3 x 10.5 x 3.6	33.7 x 26.7 x 9.1
HS24323E	White				
120114	White	240	16	26.5 x 10.5 x 3.5	67.3 x 26.7 x 9.0
120115	Blue				

Mega Racks® 18-20 mm - Holds 15 mL conical/swab tubes and 20 mm round bottom tubes

Item No.	Color	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120834	White	120	20	14.9 x 11.9 x 3.6	37.7 x 30.1 x 9.1
120835	Blue				



Clinical Tube Racks

Learn about more of our Tube Racks that hold 10-17 mm blood collection tubes

Rota-Rack®

Rotate and customize to suit your application

Page 76



OneRack® Multi

Ultimate in versatility handles multiple tube diameters and heights


Page 78



OneRack®

Tough and highly resistant, available in multiple materials and sizes


Page 83



HDPE Coated Wire Racks

Durable racks for sample collection, handling and transport


Page 84



Fold and Snap Tube Rack

Economical, easy to assemble


Page 82



Mega Rack®

High-capacity storage for 5-7 mL tubes

Page 86



Adapt-A-Rack®

Flexible rack holds variety of tube sizes

Page 79



Peg Rack

Sturdy design for quick location and retrieval of tubes

Page 86






Droplet Corrugated Sample Storage Tray

Affordable sample organization and storage

Corrugated storage trays can be used with the included dividers for specimen containers or remove them for general lab product storage or transport

- Made from corrugated polypropylene
- Stackable with other trays
- Moisture resistant, significantly reducing the chance of mold
- Pack includes: 10 trays and 10 dividers

Item No.	Color	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120375	White	24	58	16 x 10.2 x 3.5	40.6 x 26 x 8.9



Trays are supplied "flat-packed"
Eliminating the hassle of stacking and reduces the amount of space needed for storage



Removable Inserts

allow outer shell to be used for storage or transport of laboratory products



Lightweight and Sturdy

to hold up to 24 full specimen containers



Droplet Sample Storage Tray

Sample organization for the busy lab

Hard wearing storage trays that can be used for specimen containers with the insert provided, or for general lab products without the insert

- Bin and tray are made from polystyrene
- Lightweight and sturdy to hold up to 24 full specimen containers (Specimen containers not included)
- Stackable with other trays to support as many as four high
- Can be easily cleaned with common disinfectants and germicides




Removable Insert

allows outer shell to be used for storage or transportation of laboratory products

Integrated
Label Holder

Item No.	Color	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120172	Blue/White	24	58	18.9 x 13 x 3.3	48 x 33.1 x 8.4



Coated Wire Rack - Urine Containers

Durable racks for transport or sample handling

- Chemically-resistant, HDPE-coated steel
- Mesh bottom support wire rack
- Fits bottles with 55-58 mm diameters



Item No.	Color	Tube Array	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120091	Blue	6 x 4	24	58	15.5 x 10.3 x 2	39.3 x 26.2 x 5

Liquid Sample Handling

Whether it be dispensing or aspirating microliters or milliliters, we offer you precision at your fingertips



Pipettors, Fillers & Bulbs
Reagent Reservoirs & Chillers
Funnels & Lead Rings
Pipettor Racks & Accessories

High-quality Pipettors offering precision at your fingertips

RF3000®

The powerful pump offers rapid, yet sensitive aspiration and dispensing, has adjustable valves and is microprocessor controlled for constant reliability and accuracy

Page 16



RF1000

Manual pipette controller is lightweight, easy to use and uniquely designed for one-handed operation

Page 17



High-quality Pipettors offering precision at your fingertips

Pipette Controller

Manual inline pipette controller with ergonomic grip and bulb for simple and precise operation

Page 17



Pearl® Adjustable Volume Pipettes

Complete range of single- and multi-channel pipette volumes to cover all your possible needs

Page 18



Fillers and Bulbs

1 Latex Dropper Bulbs

- Designed for use with Pasteur pipettes, small pipettes, and medicine droppers

2 Dropper Bulbs

- Designed for use with Pasteur pipettes, small pipettes, and medicine droppers

3 Rubber Bulbs

- Contoured shape and heavy-wall neck increase grip for easier use with pipettes and syringes

4 Vacuum/Pressure Bulb

- Ideal for gas sampling or for pumping gas through apparatus
- Two rubber valves with PVC connections (tubing not supplied)



5 Safety Bulb Pipette Filler


- Features a white ABS attachment containing a silicone conical adapter that seals around pipettes without twisting
- Pipette does not contact hard plastic, which reduces the risk of cracking and chipping
- Elastic band for hanging on reagent bottle or wrist
- Filler fits standard size pipettes

6 Pipette Filler

- Efficient one-handed operation
- Control air-evacuation, liquid uptake, and liquid dispensing with three valves
- Pinch release-valve for smooth liquid flow
- Red pipette filler is safe for pipetting of caustics and does not contain carbon black and will not smudge


Latex Dropper Bulbs *(Natural Rubber)*

Pack of 72


Item No.	Color	Volume mL	Length in	Length mm	Outer Diameter in	Outer Diameter mm	Hole Diameter in	Hole Diameter mm
HS20622A	Yellow 	1	1.7	43	0.51	13	0.20	5
HS20622B		2	2.0	52	0.51	13	0.24	6

Dropper Vial Bulbs (*Silicone Rubber*)

Pack of 24

Item No.	Color	Volume mL	Length in	Length mm	Outer Diameter in	Outer Diameter mm	Hole Diameter in	Hole Diameter mm
HS20629A	Blue 	1	1.1	28.7	0.47	12	0.24	6.0
HS20629B		3	1.6	41	0.70	17.8	0.26	6.5

Rubber Bulbs (*Natural Rubber*)

Item No.	Color	Volume mL	Length in	Length mm	Outer Diameter in	Outer Diameter mm	Hole Diameter in	Hole Diameter mm
HS20630B	Red 	15	1.9	50	1.3	34	0.16	4
HS20630C		30	2.6	66	1.8	46	0.24	6
HS20630D		60	3.1	80	2.1	54	0.24	6
HS20630E		150	4.1	103	2.9	73	0.24	6




Vacuum/Pressure Bulb (Natural Rubber)

Item No.	Color	Volume mL	Length in	Length mm	Outer Diameter in	Outer Diameter mm	Hole Diameter in	Hole Diameter mm
HS20631B	Black 	85	4.2	106	2.1	52	0.29	7.4

Safety Bulb Pipette Filler (*Natural Rubber*)

Item No.	Color	Volume mL	Length in	Length mm	Outer Diameter in	Outer Diameter mm	Hole Diameter in	Hole Diameter mm
HS20623	Red 	15	4.2	106	1.6	41	0.39	10

Pipette Filler (*Natural Rubber*)

Item No.	Color	Volume mL	Length in	Length mm	Outer Diameter in	Outer Diameter mm	Hole Diameter in	Hole Diameter mm
HS20632A	Black 	35	5.1	130	2.1	54	0.24	6.0
HS20632B	Red 							
HS20632C	Blue 							

Transfer Pipettes

One-piece, unbreakable polyethylene plastic pipettes

Well-suited for blood banking, urinalysis, hematology, wet chemistry, microbiology, and serology

- Bulk packed versions come in handy dispenser box
- Sterile versions available single wrapped or 20 wrapped at a time
- Lot-to-lot traceability on sterile versions

	Item No.	Vol. mL	Working Vol. mL	Length mm	Graduations mL	Packaging 500 per box	Sterile
a	HS206371A	5	1	150	0.25	Bulk	
	HS206371B					20/bag	✓
	HS206371C					1/single wrap	✓
b	HS206372A	5	2	155	0.5	Bulk	
	HS206372B					20/bag	✓
	HS206372C					1/single wrap	✓
c	HS206373A	7	3	155	0.5	Bulk	
	HS206373B					20/bag	✓
	HS206373C					1/single wrap	✓

STERILE R



Universal Carousel Pipette Stand

Pipettes never looked so good, on your bench top


Stylish pipette stand holds and stores up to 8-single channel, 4-multi channel, or a combination of both pipettes

Versatile

- Patent pending design enables this stand to securely hold virtually every known brand of pipette in the marketplace
- Modular upper holder inserts accommodate the fit of different pipette widths
- Assorted color inserts allow for color coding configurations

Durable Utility

- Stable weighted base allows for easy one handed removal and replacement of pipettes
- Rotating head makes pipette selection fast and easy
- Moving the stand while keeping pipettes in place is simple with a quick grip handle allowing a user to grab and go



Inserts

on the upper holders can be added or removed to accommodate the fit of different pipette widths

Item No.	Upper Holder Color	Quantity	Lower Holder Color	Quantity
120480	Blue/Green	8 ea/color	Blue	8 ea
120481	Purple/Clear	8 ea/color	Purple	8 ea
120499	Assorted/Clear*	1 ea/color	Clear	8 ea
		8 clear		

*Assorted colors (Purple, Blue, Green, Grey, Red, Orange, Yellow and Clear)

Base Width in	Top W x H in	Base Width cm	Top W x H cm
6.2	9 x 16.1	15.8	22.9 x 40.8

Universal Pipette Holder

Wall mounted storage solution


The pipette holder comes with two wall mounting brackets that each hold two pipettes

Secure

- Upper and lower segment inserts easily slide in place to securely hold virtually every known brand of pipette
- Segments unique design includes traction grips that center pipettes allowing for easy one handed removal and replacement

Adaptable

- Inserts come in two separate colors allowing for color coding, quick identification or just to add a pop of color to your workspace




Inserts

on the upper holders can be added or removed to accommodate the fit of different pipette widths

Item No.	Lower Holder		Upper Holder		L x H in	L x H cm
120575	2 Blue	■	2 Blue	■	6.9 x 1.8	17.5 x 4.5
	2 Green	■	2 Green	■		
	2 Clear	□	1 Purple	■		
			1 Yellow	■		

Easy To Clean
Made of ABS

Patent Pending




Clamp-ette® Pipette Holder

Create an instant pipette holding point within your lab

Secures to tabletops, benches and shelves with a twist

- Chemical-resistant polypropylene and fits most major pipette brands



Secures Clamping
for tabletops, benches
and shelves

Tri Clamp-ette®

Item No.	Color	L x W x H in	L x W x H cm
HS20162A	Blue ■	6.2 x 4.7 x 3.7	15.7 x 12 x 9.5



Flip & Grip® Pipette Holder

Giving labs a choice

Unique pipette holder that gives you the option of using on the bench top or hanging under shelf to help maximize your working space

- Will hold up to 12 single-channel pipettes, or a combo of 3 single- and 3 multichannel pipettes
- Holds most common brands
- Rubber grips prevent scratching pipette
- Weighted base
- Installation hardware included

US Design Patent No. D754,872 & D763,464

Item No.	Color	Bench Top Style W x H in	Bench Top Style W x H cm	Under Shelf Style W x H in	Under Shelf Style W x H cm
120026	Blue				
120056	Lime	13.1 x 7.6	33.2 x 19.3	7.4 x 7.6	18.7 x 19.3
120057	Magenta				





Pipette Workstation

Organize your pipettes and tips

Keep pipette and disposable tips together in one convenient stand

- Acrylic workstation is opaque on sides with clear front and bin lids for easy content identification
- Holds up to four pipettes on top shelf
- Dual bins hold multiple size tips
- Self-closing lid provides dust-free storage
- Oversized lip allows for easy access when wearing gloves

Item No.	Color		L x W x H in	L x W x H cm
HS234631	Clear	<input type="checkbox"/>	9.6 x 6.5 x 9.7	24.4 x 16.5 x 24.6



Pipette Tip Disposal Box

Quick flip of port cover makes tip ejection easy

Convenient for one-handed disposal of used pipette tips, eliminates any need to handle pipette tips used with hazardous fluids

- The 0.375 in (10 mm) thick acrylic provides protection from Beta radiation
- The large, hinged lid allows easy disposal of contents
- Use in conjunction with small biohazard disposal bags, if hazardous fluids are not acrylic friendly

Item No.	Lid Opening mm	L x W x H in	L x W x H cm
HS234634	30	6.1 x 6.7 x 8.1	15.4 x 17.1 x 20.6



3-Place Stand Acrylic Clear

4-Place Stand Acrylic Clear



6-Place Stand Acrylic Blue



6-Place Stand ABS Blue

Pipette Stands

Convenient and safe storage

3-, 4-, and 6-Place Pipette Stands are designed to hold most major pipette brands

- Acrylic stands are transparent/clear; ABS stands are solid opaque colors

Item No.	Color		Material	Pipette Capacity	Slot Upper in	Slot Upper mm	Slot Lower in	Slot Lower mm
HS206203	Clear	<input type="checkbox"/>	Acrylic	Three	1	25	0.8	19
HS206204				Four				
HS20620				Six				
HS20613C	Blue	<input checked="" type="checkbox"/>	Acrylic	Six	1	25	0.8	20
HS20613E	Blue	<input checked="" type="checkbox"/>	ABS	Six				

Cargo Pipette Rack

Meets the need of any lab

Racks securely lock together horizontally or vertically to desired configuration or can be used individually

Can be used on top of benchtop or conveniently hung on any metallic surface with strong magnets that are designed into the racks

- Angled shelves store 0.1 to 50 mL pipettes
- Protective lid encloses the unit to prevent dust build-up and provides additional shelf space

Pack of 3

Item No.	Color	L x W x H in	L x W x H cm
HS20615M	Blue	11.9 x 3.9 x 5	30.1 x 10 x 12.8




Magnets

connect securely to most metallic surfaces

Units Lock Securely

together horizontally or vertically




Pipette Racks

Protects and organizes glass and plastic pipettes

Four angled-shelved compartments keep your 0.1 to 50 mL pipettes or cryo vial canes within easy reach while preventing damage and clutter

- Heavy-duty acrylic version provides crystal-clear viewing from any angle
- ABS versions have cut-out viewing window and come with two removable attachments that hold pipette bulbs, fillers, and pumps
- HS20612M magnetized ABS version has strong magnets on both sides and mounts securely to steel surfaces

Item No.	Color	Material	L x W x H in	L x W x H cm
HS20611	Clear	Acrylic	11.4 x 3.5 x 15.8	29.1 x 8.8 x 40.1
HS20612	Blue	ABS		
HS20613	Red	ABS		
HS20612M	Blue	ABS (mounting magnets)		



Pipette Attachments

ABS versions includes 2 removable pipette attachments

Reservoir Chiller

Ice-free, non-toxic way to extend working period and preserve integrity for temperature sensitive reagents

- Reagents are kept at less than 5°C (41°F) for approximately 1 hour
- Engineered custom fit design for Heathrow Scientific 50 or 100 mL reservoirs to ensure a maximum contact for longer uniform cooling
- Eliminates the risk of contamination from messy ice bath
- Not for use with liquid nitrogen



Stackable
Saves storage and work space

Symmetrical cooler design

Allows reagent reservoirs to be mounted in either direction



Pack of 2

Item No.	Color	Volume mL	L x W x H in	L x W x H cm
120730	Purple	50	6.1 x 2.9 x 1.7	15.5 x 7.4 x 4.2
120731	Blue			
120732	Assorted*			
120733	Purple	100	6.2 x 9 x 1.7	15.7 x 9 x 4.2
120734	Blue			
120735	Assorted*			

*Assorted colors (Purple and Blue)



50 mL
Reservoir Chiller



100 mL
Reservoir Chiller



**2 Basins
in One!**



Dual Solution Reservoir

Innovative shape guides and protects pipette tips

Multi-functional reagent reservoir has a single compartment 50 mL basin on one side and 12 slotted 5 mL basins on the other side

- Angled wells guide pipette tips toward a center position to help avoid cross contamination and waste
- Individual 5 mL basin wells numerical imprinted to identify samples
- Cover included to prevent contamination

Item No.	Volume mL	L x W x H in	L x W x H cm
HS20821A	5/50	4.9 x 2 x 1.5	12.5 x 5 x 3.7

Pack of 25



Reagent Reservoirs

Reagent Reservoirs by Heathrow Scientific offer a fresh and modern design with enhanced features to streamline pipetting procedures

- Superior durability with enhanced walls and crush-resistant "turrets" at the corners
- Easy handling with tabs to prevent stacked reservoirs from sticking together

Sterile Versions

- Manufactured and gamma irradiated in an ISO certified facility
- Certified DNase, RNase and pyrogen free
- Packaged individually in an easy open HDPE package
- Traceability Certificate available upon request

STERILE R

Divided Reservoir Design

allows pipetting of two different reagents or maximizes recovery if using only single channel and small amount of reagent

Graduation Markings

on the top edge and inner walls allows measurements of liquid to be easily identified and reduces wastage of expensive reagents

Pour-off Spouts

conveniently located on all four corners enables easy spill-free pouring and reduces chance of dripping

Item No.	Volume mL	Graduation Marks mL	Material	Sterile	Packaging	Interior Dimensions L x W x H in	L x W x H cm
120637	5	5	Polystyrene	Yes	50/Single Bag	2.8 x 0.8 x 0.7	7.2 x 2.1 x 1.8
120638	10	5 and 10	Polystyrene	No	100/Bulk	2.8 x 0.8 x 0.7	7.2 x 2.1 x 1.8
120639	10	5 and 10	Polystyrene	Yes	100/Single Bag		
120640	25	15 and 25	Polystyrene	Yes	200/40 x 5pk	4.4 x 0.9 x 1.1	11.3 x 2.4 x 2.7
120641	25	15 and 25	Polystyrene	Yes	100/Single Bag		
120642	25	15 and 25	Polystyrene	No	100/Bulk		
120884	25	15 and 25	Polystyrene	Yes	100/Single Bag		
120643	50	25 and 50	Polystyrene	Yes	200/40 x 5pk	4.5 x 1.1 x 1.4	11.5 x 2.9 x 3.5
120644	50	25 and 50	Polystyrene	Yes	100/Single Bag		
120645	50	25 and 50	Polystyrene	No	100/Bulk		
120646	50	25 and 50	PVC	Yes	100/20 x 5pk		
120690	50	25 and 50	PVC	No	100/Bulk		
120647	100	50 and 100	Polystyrene	Yes	200/40 x 5pk	4.6 x 1.9 x 1.0	11.7 x 4.8 x 2.5
120648	100	50 and 100	Polystyrene	Yes	100/Single Bag		
120649	100	50 and 100	Polystyrene	No	100/Bulk		
120650	100	50 and 100	PVC	No	100/Bulk		
120651	25 (10 and 15 mL divided)	10/15	Polystyrene	Yes	200/40 x 5pk	4 x 1 x 0.8	10.2 x 2.5 x 2.1
120652	25 (10 and 15 mL divided)	10/15	Polystyrene	Yes	100/Single Bag		

Patent Pending



* Not recommended
for hot liquids

Eco smartFunnel™

Paper Funnel

A clean transfer every time

- Supplied flat-pack to minimize storage
- With a minimal amount of finger pressure the funnel takes its working shape
- Maximum flexibility – cut the stem to any size or shape to meet with your needs
- Fits a wide variety of openings and can be used with virtually any liquid* or powder

Pack of 25

Item No.	Capacity mL	Top Dia. in ²	Top Dia. cm ²	Bottom Dia. in ²	Bottom Dia. cm ²	Height in	Height cm
120219	186	2.4	6	0.3	0.6	5	12.8

Eco smartFunnel™ is a trademark of LevGo Inc.



Lead Rings

Stabilize vessels and labware on the bench or in water baths

- Vinyl-coated lead rings will not mark or scratch benchtops and labware
- C-shape format allows some manual size variation for a customized fit

Item No.	Shape	Color	Flask Sizes mL	ID in	ID cm	OD in	OD cm	Weight g
HS8881A	C	Blue ■	125–500	1.9	4.2	3	7.5	215
HS8882A	O			1.9	4.8	3	7.7	210
HS8881B	C	Blue ■	250–1000	2	5.1	3.5	9	465
HS8882B	O			2	5.1	3.5	9	490
HS8881C	C	Blue ■	500–2000	2.1	5.4	3.9	9.8	595
HS8882C	O			2.2	5.7	3.9	9.8	610
HS8881D	C	Blue ■	1000–4000	2.6	6.6	4.4	11.3	890
HS8882D	O			2.8	7	4.5	11.5	935

*Weights may vary by +/- 3%

Sample Handling

Complete range of products and tools for gathering, weighing, and working with various types of samples



Pill Counting Trays

Weigh Boats, Dishes, Papers

Sampling Scoops

Utensils: Spatulas, Forceps, Cell Spreaders, Loops & Picks

Pill Counting Trays

Counting trays that offer safety and convenience when counting pills and capsules

Intended for single use to prevent cross-contamination for batch sampling or distribution of pills and capsules

Safety

- Disposable to ensure no chemical cross-contamination for patient safety
- Made from FDA food grade polystyrene in an ISO Class 8 clean room
- Meets quality control and GMP requirements for manufacturing batch sampling
- Durable and shatterproof trays constructed of inert, BSE/TSE free material
- Lot traceable


Convenience

- Pour spout in corner for dispensing control
- Ergonomic spatula (length - 6.9 in/17.6 cm)
- Available in both, sterile* and non-sterile for intended usage

* Electron beam sterilized and individually wrapped with Certificate of Sterility available upon request



STERILE R



For round pills



For capsules from size 2 to 00

Item No.	Description	L X W X H in	L x W x H cm	UOM
120845	Pill Counting Tray Rectangular Sterile	10.3 x 5.2 x 0.5	26.2 x 13.1 x 1.2	Case 40; 4 x 10/Single Bagl
120846	Pill Counting Tray Rectangular Non-Sterile	10.3 x 5.2 x 0.5	26.2 x 13.1 x 1.2	Case 40; 4 x 10/Bulk
120847	Pill Counting Tray Triangular Sterile	9.1 x 9.1 x 0.5	23.0 x 23.0 x 1.3	Case 30; 3 x 10/Single Bag
120848	Pill Counting Tray Triangular Non-Sterile	9.1 x 9.1 x 0.5	23.0 x 23.0 x 1.3	Case 30; 3 x 10/Bulk

Aluminum Weighing Boats & Dishes

Dishes for general-purpose weighing, dispensing, storage, food processing, evaporating or drying

Aluminum Weighing Boats

- Sturdy yet lightweight boats feature a smooth, flat bottom for sample handling and a flanged, crimped wall for rigidity
- Finger tab handle can be etched for identification prior to weighing



Aluminum Weighing Boat



Aluminum Weighing Dish

Aluminum Weighing Dish

- Designed for moisture determination balances with slightly tapered walls for stacking

Item No.	Size	Dia. x H in	Dia. x H mm	Volume ml	Packaging
HS14521A	Small	1.7 x 0.5	43 x 12	20	100
HS14521B	Medium	2.2 x 0.6	57 x 14	39	100
HS14521C	Large	2.8 x 0.6	70 x 16	50	100
HS14522	Large	4 x 0.3	101 x 8	30	50

smartPans™

Ideal for microbalances, organic micro elementary and thermal analysis, and pharmaceutical use

- Made of Ultra High Vacuum (UVH) Foil - "A" wettable and completely free from oils normally found on other grades of aluminum foil
- Foil meets ASTM B-479, specification 3.1.4 and 10.3.1
- "Cahn style" accessory pans are compatible with most commercially available analyzers and microbalances



Pack of 20

Item No.	Dia. x H in	Dia. x H mm
120627	1 x 0.1	25 x 3.6

smartPans™ is a trademark of LevGo Inc.

Weighing Papers

Smooth surface for good sample recovery

- Nitrogen-free
- Strong enough to support dense material




Pack of 500

Item No.	L x H in	L x H cm
120116	3 x 3	7.6 x 7.6
120117	4 x 4	10.2 x 10.2
120118	6 x 6	15.2 x 15.2

Weighing Boats

Flexible polystyrene boats feature smooth surfaces to provide accurate pour-outs with minimal sample loss

- Select white, black, or blue for best color contrast to ensure sample is completely removed
- Biologically inert and resistant to dilute and weak acids, aqueous solutions, alcohols, and bases
- Temperature range -10 to 70°C (14 to 158°F)
- Antistatic versions available making the boats nonreactive with most substances and suitable for weighing static-affected samples
- Wide flat bottoms ensure stability and prevent tipping
- Rounded corners allow simple flexing for controlled pouring of liquid, powdered or granular samples
- Can be used with microwave moisture analyzers
- Made from FDA food grade polystyrene



Standard Weighing Boats

Manufactured
and packaged
in a clean
room



Sterile Weighing Boats

Ideal for applications where there is concern of cross contamination and 'clean' is a necessity in life science research, food, biopharmaceutical, environmental and molecular research environments

- Electron beam sterilized
- Lot number traceability for quality assurance
- Certificate of Sterility available upon request


STERILE R



Square Weighing Boats



Hexagonal Weighing Boats



Diamond Weighing Boats



Pour-boat Weighing Dishes

Weigh Boats

Pack of 500

Shape	Regular Material Item No.	Antistatic Material Item No.	Color	Size	Volume mL	L x W x H in	L x W x H cm
Standard	HS1420A	HS1420AA	White <input type="checkbox"/>	Small	10	1.8 x 1.8 x .3	4.6 x 4.6 x .8
	HS1421A	HS120222	Blue <input checked="" type="checkbox"/>				
	HS1423A	HS1423AA	Black <input checked="" type="checkbox"/>				
	HS1420B	HS1420BB	White <input type="checkbox"/>	Medium	85	3.4 x 3.4 x .9	8.5 x 8.5 x 2.4
	HS1421B	HS120223	Blue <input checked="" type="checkbox"/>				
	HS1423B	HS1423BB	Black <input checked="" type="checkbox"/>				
	HS1420C	HS1420CC	White <input type="checkbox"/>	Large	280	5.5 x 5.5 x .9	14 x 14 x 2.2
	HS1421C	HS120224	Blue <input checked="" type="checkbox"/>				
	HS1423C	HS1423CC	Black <input checked="" type="checkbox"/>				
Hexagonal	—	HS14251A	White <input type="checkbox"/>	Small	6	1.5 x 1 x .4	3.8 x 2.5 x .9
	—	HS14251B		Medium	46	2.5 x 2 x .7	6.4 x 5.1 x 1.9
	—	HS14251C		Large	194	4 x 3 x .9	10.2 x 7.6 x 2.2
	—	HS14251D		X-Large	364	5 x 3.5 x 1.2	12.7 x 8.9 x 3.1
Pour-boat	HS1419A	HS120225	White <input type="checkbox"/>	Small	14	1.7 x 2.3 x .5	4.3 x 5.8 x 1.3
	HS1419B	HS120226		Medium	137	3.3 x 5.2 x 1	8.3 x 13.2 x 2.6
	HS1419C	HS120227		Large	259	4.3 x 7.2 x 1	10.8 x 18.3 x 2.6
Square	HS1422	—	White <input type="checkbox"/>	Medium	51	3.5 x 3.5 x .3	9 x 9 x .7
Diamond	HS1424AA	HS1425AA	White <input type="checkbox"/>	Small	5	2.1 x 1.3 x .2	5.3 x 3.3 x .6
	HS1427AA	HS1426AA	Black <input checked="" type="checkbox"/>				
	HS1424A	HS1425A	White <input type="checkbox"/>	Medium	25	2.6 x 1.9 x .6	6.7 x 4.7 x 1.4
	HS1427A	HS1426A	Black <input checked="" type="checkbox"/>				
	HS1424B	HS1425B	White <input type="checkbox"/>	Large	30	3.1 x 2.2 x .6	7.8 x 5.6 x 1.4
	HS1427B	HS1426B	Black <input checked="" type="checkbox"/>				
	HS1424C	HS1425C	White <input type="checkbox"/>	X-Large	100	4.5 x 3.3 x 0.7	115 x 85 x 17
	HS1427C	HS1426C	Black <input checked="" type="checkbox"/>				

Sterile Weigh Boats

Pack of 200/40 x 5pk

Shape	Regular Material Item No.	Antistatic Material Item No.	Color	Size	Volume mL	L x W x H in	L x W x H cm
Standard	—	120880	White <input type="checkbox"/>	Small	10	1.8 x 1.8 x .3	4.6 x 4.6 x .8
Standard	—	120710	White <input type="checkbox"/>	Medium	85	3.4 x 3.4 x .9	8.5 x 8.5 x 2.4
Hexagonal	—	120881	White <input type="checkbox"/>	Small	6	1.5 x 1 x .4	3.8 x 2.5 x .9
Hexagonal	—	120882	White <input type="checkbox"/>	Medium	46	2.5 x 2 x .7	6.4 x 5.1 x 1.9

Micro Weighing Boats

Ideal for microbalance weighing and/or distribution of very small amounts of sampling material

Easy to Handle

- Large tabs make boats easy to use, position and hold, especially with forceps

Superb Sample Recovery

- Inert material enables the boat and contents to be added to solutions to aid 100% sample recovery, where protocols allow
- Anti-static properties prevent sample retention from static electricity
- Small surface area reduces liquid sample evaporation


Micro Weigh Boats

Item No.	Color	Dimensions in	Dimensions mm	UOM
120836	White	.98 x 0.4 x 0.3	25 x 8.9 x 7.5	250/pk



Easy to Handle Tabs

Made of Polystyrene
Food Grade Compliant



Actual size

Ideal for very small volume pipette calibration and checks and micro sampling





Functional in the Lab

Available in a small and medium in both brown for light colored samples and white for dark samples



The mark of responsible forestry

smartBoats™

A truly green solution to sample weighing

Made from FSC® Certified (FSC® Mix Credit) paper (Eco-IN) and when disposed of in the standard waste stream they will naturally biodegrade (Eco-OUT)

- A cost-effective solution to reduce plastics in the lab
- Made of grease-proof paper that meets FDA 21 CFR 176.170 (components of paper in contact with aqueous and fatty foods) and with 21 CFR 176.180 (contact with dry foods)
- Control contents like a polystyrene weigh boat because they are easier to pour
- An alternative to regular weighing papers because they are easier to handle

Item No.	Color	Size	Packaging	L x W x H in	L x W x H cm
120312	Brown	Small	Pack of 500	1.4 x 1.4 x 0.6	3.5 x 3.5 x 16
120314	White	Small	Pack of 500	1.4 x 1.4 x 0.6	3.5 x 3.5 x 16
120313	Brown	Medium	Pack of 500	1.9 x 1.9 x 0.6	4.8 x 4.8 x 14
120315	White	Medium	Pack of 500	1.9 x 1.9 x 0.6	4.8 x 4.8 x 14
120626	Brown	Large	Pack of 250	3.8 x 3.8 x 1	9.5 x 9.5 x 2.5
120195	White	Large	Pack of 250	3.8 x 3.8 x 1	9.5 x 9.5 x 2.5

smartBoats™ is a trademark of LevGo Inc.

Forceps

Corrugated handles and serrated tips

- ABS plastic

* Forceps shown at actual size

Pack of 6

Item No.	Color	Length in	Length cm	Tip
HS15998C	Blue	4.4	11.1	Straight



HS15906



HS15907



HS15909

Spatulas

Transfer, lift, and handle powders and granular materials

- Made from nickel-plated, stainless steel wire
- Glazed spatula HS15906 has a Hayman-type end



Item No.	End Type	End L x W in	End L x W mm	Overall Length in	Overall Length cm	Finish
HS15906	Flat	1.2	30	6	15	Glazed
	Spoon	0.5 x 0.2	12 x 4			
HS15907	Round Flat	2 x 0.3	50 x 8	7.7	19.5	Mirror
	Tapered Flat	1.8 x 0.3	45 x 6.5			
HS15909	Flat	2	50	8.9	22.5	Mirror
	Spoon	1.5 x 0.6	38 x 15			

smartSpatula®


Designed to fit perfectly with standard
0.2, 0.5, 1.5/2.0, 15 and 50 mL tubes

- Made of virgin, food-grade polypropylene
- Ideal for dispensing or sample collection of powders or soft solids
- Will not shatter when used with liquid nitrogen frozen samples
- Hollow core minimizes heat transfer—good for hot and cold liquids
- Antistatic, sterile and Eco-Friendly versions available
- Complies with FDA standards for food contact

AUTOClavABLE
121° C

STERILE EO

Always there
and always
clean!



Lab tools that scoop, stir and scrape.



Item No.	Description	Color		Suggested Tube mL	Packaging	L x W in	L x W cm	Spoon (Flare) in	Spoon (Flare) mm	Scoop/Pick L x W in	Scoop/Pick L x W mm	Taper
a 120005	Micro Antistatic	Opaque	<input type="checkbox"/>	0.2 and 0.5	300/Bulk	5.5 x .2	14 x .4	—	—	4.3	11 (V-Tip)	15°
b 120008	Eco-Friendly	Opaque	<input type="checkbox"/>	15 and 1.5	300/Bulk	5.5 x .3	14 x .7	9.8 x 5.5	25 x 14	11.8 x 1.2	30 x 3	
c 120003	Standard	Blue	<input checked="" type="checkbox"/>		300/Bulk							
d 120004	Standard	Opaque	<input type="checkbox"/>	15 and 1.5	300/Bulk							
d 120007	Std. EtO Sterile	Opaque	<input type="checkbox"/>		100/10 x 10pk	8.3 x .3	21 x .7	9.8 x 5.5	25 x 14	11.8 x 1.2	30 x 3	
d 120625	Std. EtO Sterile*	Opaque	<input type="checkbox"/>		100/10 x 10pk							
e 120006	Macro	Green	<input checked="" type="checkbox"/>	50	150/Bulk	12.2 x .4	31 x 1	16.5 x 9.3	42 x 23.5	11.8	30	35°

*Non Toxic/Non Pyrogenic smartSpatula™ is a trademark of LevGo Inc.

Dual Purpose Scoops

Just scoop to retrieve, place on scale to weigh

Unique economical solution that combines the function of a scoop and a weigh boat into one, streamlining the process

- Made from anti-static polystyrene that is waterproof and reduces the chance of materials sticking to the scoop
- Ships and stores flat, easy to assemble for use
- Locking tab feature enables a scoop shape to help retain materials
- Coated to stand up to powders, gels, and semi-solids



Case of 100 10 x 10pk

Item No.	Color		Size	L x W in	L x W cm
120656	White	<input type="checkbox"/>	Large	5 x 3.5	12.7 x 9
120657	White	<input type="checkbox"/>	Small	5 x 2.5	12.7 x 6.3



Sampling Scoops

Our line of durable Sampling Scoops designed for easy sample retrieval of granules, powders and pastes

Available in multiple sizes and intended for single use to prevent contamination of bulk sample material

Durable Design

- Large, wide openings for easy sample retrieval
- Flat bottom to stabilize sample on scale or flat surface
- Lot traceable
- Shatterproof constructed of inert, BSE/TSE free material

Safe Sample Recovery

- Disposable to ensure no chemical cross-contamination
- Made from FDA food grade polystyrene in an ISO 8 Class clean room
- Meets quality control and GMP requirements for sampling
- Available in both, sterile* and non-sterile for intended usage

*Electron beam sterilized and individually wrapped
Certificate of Sterility available upon request



Manufactured
and packaged
in a clean
room

STERILE R

60 mL

125 mL

250 mL

500 mL





Small Sampling Scoop

Item No.	Description	Volume	L x W x H in	L x W x H mm	Sterile	UOM
120858	Small Sampling Scoop Non-Sterile	60 mL/2 oz.	8.1 x 1.6 x 1.1	205 x 40 x 27		Case 200; 20 x 10/Bulk
120857	Small Sampling Scoop Sterile	60 mL/2 oz.	8.1 x 1.6 x 1.1	205 x 40 x 27	✓	Case 200; 20 x 10/Single Bag



Medium Sampling Scoop

Item No.	Description	Volume	L x W x H in	L x W x H mm	Sterile	UOM
120856	Medium Sampling Scoop Non-Sterile	125 mL/4 oz.	8.6 x 1.9 x 1.3	218 x 49 x 38		Case 200; 20 x 10/Bulk
120855	Medium Sampling Scoop Sterile	125 mL/4 oz.	8.6 x 1.9 x 1.3	218 x 49 x 38	✓	Case 200; 20 x 10/Single Bag



Large Sampling Scoop

Item No.	Description	Volume	L x W x H in	L x W x H mm	Sterile	UOM
120854	Large Sampling Scoop Non-Sterile	250 mL/8 oz.	7.7 x 2.3 x 1.7	195 x 58 x 44		Case 200; 20 x 10/Bulk
120853	Large Sampling Scoop Sterile	250 mL/8 oz.	7.7 x 2.3 x 1.7	195 x 58 x 44	✓	Case 200; 20 x 10/Single Bag



XLarge Sampling Scoop


Item No.	Description	Volume	L x W x H in	L x W x H mm	Sterile	UOM
120852	XLarge Sampling Scoop Non-Sterile	500 mL/16 oz.	10.0 x 3.6 x 2.1	255 x 92 x 53		Case 200; 20 x 10/Bulk
120851	XLarge Sampling Scoop Sterile	500 mL/16 oz.	10.0 x 3.6 x 2.1	255 x 92 x 53	✓	Case 200; 20 x 10/Single Bag

smartPicks® Inoculation Picks

The 'smart' way to isolate individual colonies

Available in three lengths and two needle sizes to suit a wide range of application needs

- Naturally biodegradable, made from 100% U.S.A. grown white birch wood
- Manufactured in an FDA-registered factory and complies with FDA standards for food contact
- Available in semi or needle point styles for delicate operations



Mini Sized SmartPicks

fit neatly into a 1.5 mL microfuge tube allowing the lid to be closed for full sample transfer/recovery

Drop it in & Close it up!



Item No.	Size	Point	Tip	Suggested Tube mL	Packaging	Length in	Length cm
120010	Mini	Needle Point	Single Tip	Fits 1.5 tube (closed)	2,000	1.3	3.3
120011	Standard	Semi Point			1,600	2.5	6.4
120012	Standard	Needle Point					

smartPicks® is a trademark of LevGo Inc.



Inoculating Needle/Loops

Combination provides choice to spread or pick colonies and cultures

Transfer bacterial cultures without risk of contamination

- Use loop end to uniformly streak a gel surface
- Use needle end for removal of colonies
- Choose between flexible polypropylene or the more rigid ABS versions
- Lot-to-lot traceability

STERILE EO



Item No.	Color	Material	Loop μL	L x D in	L x D cm	Packaging 1,000 per Box	Sterile
HS81121A	Blue	Flexible PP	1	7.9 x .2	20 x .4	50 x 20pk	✓
HS81121B		Rigid ABS					
HS81121C	Blue	Flexible PP	10	7.9 x .2	20 x .4	50 x 20pk	✓
HS81121D		Rigid ABS					

Cell Spreaders

Delta™

- Made from ABS
- Flexible rippled handle and domed spreading bar to help reduce agar tears

STERILE R

L-Shape

- Made from ABS
- Lot-to-lot traceability
- 10 Pack comes in a zip bag

STERILE R

Stainless Steel

- Flame-sterilizable
- Reusable

**AUTOCCLAVABLE
121° C**

Domed Headpiece
promotes smooth spreading

Item No.	Shape	L x W x H in	L x W x H cm	Packaging	Sterile
HS8151	Delta™	6.8 x 1.5 x .2	17.2 x 3.7 x .4	500/Single Wrap	✓
HS8171A	L-Shape	5.7 x 1.5 x .2	14.6 x 3.8 x .4	500/Single Wrap	✓
HS8171B				500, 50 x 10 Pack	✓
HS86655	Stainless Steel	8 x 1.3 x .2	20.3 x 3.3 x .4	1/Single Wrap	

Need help counting cells?

For accurate cell counts, check out our Digital and Manual cell counting solutions



Hand Tally Counter, Page 21



eCount, Page 20



Benchtop Counters, Page 21

Petri Dish Dispenser

- Sturdy acrylic holds up to 30 x 100 mm dishes in two easy-access partitions
- Set on counter or wall mount with included hardware

Item No.	L x W x H <i>in</i>	L x W x H <i>cm</i>
HS23452	8.7 x 4.2 x 9.3	22.2 x 10.7 x 23.6



Tubes & Vials

Reliable sample storage



Micro & Conical Centrifuge Tubes

Screw-top Tubes

Cryogenic Vials

Micro & Conical Centrifuge Tubes

Available in natural and pure black for light sensitive samples


Microcentrifuge Tubes:

0.2, 0.5 and 1.5 mL (a-d)

- Attached flat caps for easy one-handed use
- 0.2 mL is suitable for thermal cyclers
- 0.5 and 1.5 mL versions have a frosted writing area and a syringe pierceable cap
- 100% polypropylene



Certified DNase,
RNase & Pyrogen Free
(a, b, c only)



Microcentrifuge Tubes

	Item No.	Color		Vol. mL	Outer Dimension mm	Closed Height mm	Graduations mL	Max RCF Rating	Packaging
a	HS4423	Natural	<input type="checkbox"/>	0.2	6.9	22.2	None	15,000 x g	1,000/Bulk
b	HS4422	Natural	<input type="checkbox"/>	0.5	10.1	31.7	0.1, 0.3, 0.5	15,000 x g	1,000/Bulk
c	HS4323	Natural	<input type="checkbox"/>	1.5	10.8	39.0	0.1, 0.5, 1	16,400 x g	500/Bulk
c	HS4325	Natural	<input type="checkbox"/>	1.5	10.8	39.0	0.1, 0.5, 1	16,400 x g	1,000/Bulk
d	HS4323K	Black	<input checked="" type="checkbox"/>	1.5	10.2	40.6	0.1, 0.5, 1, 1.5	18,000 x g	500/Bulk

Conical Centrifuge Tubes:

15 and 50 mL (e-i)

- Lot-to-lot traceability with easy tear packaging
- Writing area and flat polyethylene lids for easy labelling
- Imprinted graduations and a self-standing version (only in clear bulk pack)
- Recyclable chipboard used for the racked versions, making them more environmentally friendly



Skirted

Certified DNase, RNase
& Pyrogen Free (f & h)

Conical Centrifuge Tubes

	Item No.	Color		Vol. mL	Outer Dimension mm	Closed Height mm	Graduations mL	Max RCF Rating	Packaging 500 per Case
e	HS4426B	Natural	<input type="checkbox"/>	15	17	120	0.5	11,000 x g	Bulk/10 x 50pk
	HS4426R	Natural	<input type="checkbox"/>	15	17	120	0.5	11,000 x g	Rack/10 x 50pk
f	HS4428	Black	<input checked="" type="checkbox"/>	15	17	120	0.5	9,400 x g	Bulk/10 x 50pk
g	HS4427B	Natural	<input type="checkbox"/>	50	29	114	5.0	9,500 x g	Bulk/20 x 25pk
	HS4427R	Natural	<input type="checkbox"/>	50	29	114	5.0	9,500 x g	Rack/20 x 25pk
h	HS4429	Black	<input checked="" type="checkbox"/>	50	29	114	5.0	9,400 x g	Bulk/20 x 25pk
i	HS4427S	Natural	<input type="checkbox"/>	50	29 - skirted	114	5.0	9,500 x g	Bulk/20 x 25pk

Screw-top Tubes with O-ring Cap

Safer long-term storage

- Screw-on polypropylene lids with TPE O-ring (white)
- Clear polypropylene tubes provide excellent view of samples
- Skirted, self-standing base makes filling and retrieval easy
- Non-sterile



Frosted
writing area

TPE O-ring
to create an excellent seal
for secure storage

Item No.	Color	Vol. mL	Outer Dimensions mm	Closed Height mm	Graduations mL	Max RCF Rating	Packaging
HS10060	Natural <input type="checkbox"/>	2.0	13.3	47.1	0.5, 1, 1.5	18,000 x g	500/Bulk

Cryogenic Vials

Innovative design saves space in your freezer

- Puncture/rupture-resistant polypropylene
- Printed writing area and graduations
- Silicone lid seal & patented lid with ergonomic 1.5 twist to open and close
- Recommended for mechanical freezers and vapor phase LN₂ (-196°)
- Flush cap saves space during storage



STERILE R



Patented
easy-grip faceted lid

Internal Thread

Writing Area
& Graduations



Traceable Sterility Certificate

Sterile, non-pyrogenic, non-cytotoxic, DNase/RNase-free vials

Item No.	Color	Vol. mL	Packaging	Outer Dim. mm	Grad. mL	Closed Height mm
HS23202A	Assorted <input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input checked="" type="checkbox"/> *	2.0	500/50 x 10pk	12	12	49
HS23202N	Natural <input type="checkbox"/>					

US Patent No. D595,424

*Assorted colors: Blue, Red, Green, Purple, Grey

Explore our Tube Racks & Storage Products

Largest and most innovative portfolio in the Industry



Tube Racks & Storage

Classic, adaptable, interlocking racks and everything in between



PCR & Microtube
Adaptable & Multi-Tube
Classic Single Racks

We like to push boundaries

We make new and better ways of lab work possible by supplying everyday lab-ware and equipment that solve problems, improves productivity, and saves space.



We strive to create versatile, adaptable and multi-functional products that deliver value.



With a spectrum of colors and modern styles, our products stimulate the mind.



PCR Tube Racks 96-Well

Perfect for pre- or post-PCR sampling

- Made of polypropylene
- Accommodates individual tubes or strips
- Removable hinged lid lies flat
- Alphanumeric identification for quick location and retrieval of tubes



Pack of 5

Item No.	Color	Item No.	Color
HS23461A	Assorted*	HS23461C	Blue
HS23461B	Natural	HS23461D	Green

*Assorted colors (Blue, Green, Pink, Yellow, Orange)

Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
0.2	96	5.2	5 x 3.9 x 1.2	12.7 x 9.9 x 3.0



PCR Tube Racks Reversible

Removable hinged lid fits both sides

- Made of polypropylene
- PCR side of the rack holds 8 or 12 x 0.2 mL tube strips or individual tubes
- Reverse side holds 0.5 mL tubes
- Both sides of the rack have 12 wells for 1.5 mL tubes

Pack of 5

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in (cm)
HS2344A		0.2	168	5.5	8.1 x 5.2 x 2.1
		0.5	40	7.6	(20.6 x 13.1 x 5.3)
		1.5/2.0	12	11.4	

*Assorted colors (Blue, Green, Pink, Yellow, Orange)



Engineering, OEM & Industrial Design

- We have in-house product engineering and R&D teams that are dedicated to delivering quality, innovative, efficient, and compact solutions that customers can rely on.
- We have a range of capabilities from private labeling to product development, design, packaging and manufacturing.



PCR Cooler, 96-Well

Ice-Free Cooling System for PCR Preparation and Storage

Cools industry standard format PCR tubes, strips and single format as well as most standard and low profile 96-Well plates

Consistent & Reliable

- Modern design cools faster in freezers and stays cold longer on the bench top
- Maintains the entire 96-Well plate at a temperature of 0 °C for over an hour
- Maintains sample temperature of less than 5°C for up to three hours (with lid off)
- When the temperature ceiling of 6°C has been exceeded, users are signaled by a contrasting color change of the cooler or a single well

Alphanumeric Well Identification

for quick location and retrieval of tubes




Lid doubles as a stand
enables the rack to stay cold longer as there is no conductive heat transfer from surface beneath



Two Color Contrasting Options Available


purple to pink or dark blue to light blue allowing easy monitoring of samples



US Patent Pending, Community Design No: 004101756-0001-0002 Pack of 2

Item No.	Warm Color	Cold Color	Wells	Tube Sizes mL	L x W x H in	L x W x H cm
120727	Pink	Purple	96	0.2 Assay Plates	4.3 x 5.8 x 1.9	11 x 14.8 x 4.8
120728	Light Blue	Dark Blue				
120729	Assorted*					

*Assorted colors (Light Blue/Dark Blue & Pink/Purple)



Low-Temp PCR Rack: 96-Well

Slide, set, work & store PCR samples

Multi-purpose rack converts from compact workstation to lidded storage solution accommodating individual and strip 0.2 mL tubes

Versatile & Robust Solution

- Polycarbonate lid serves as cover and elevated rack base - just slide open and pivot lid under rack for compact workstation
- Lid can be used in right or left-handed orientation and is detachable
- Stackable saving valuable bench-top and freezer space
- Multiple colors offered for coding options



US Patent Pending, Community Design No: 004101756-0001-0002 Pack of 5

Item No.	Color	Item No.	Color
120538	Assorted*	120540	Green
120539	Blue	120541	Purple

*Assorted colors (Purple, Blue, Green, Orange, Yellow)





Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
0.2	96	5.2	4.9 x 3.6 x 1.5	12.5 x 9.2 x 3.7



Tabs on each end
of rack snap lid into place

Raised Rings
around the 1.5 to 2.0
mL wells make tubes
easier to remove



Molded Grid Lines

0.5 mL wells separated by molded-in
grid lines, making them easier to see

Reversible Tube Racks 96-Well

Easy-to-grasp tabs for sample handling and storage

96 wells for 0.5 mL tubes on one side and 96 wells for 1.5/2.0 mL tubes on the reverse side

Optimum Function & Identification

- Wells arranged in standard 8 x 12 array
- Molded in-grid lines separate 0.5 mL wells making them easier to see
- Raised rings around the 1.5 mL/2.0 mL wells make tubes easier to remove
- Alphanumeric identification for quick location and retrieval of tubes

Secure lid

- Translucent lid snaps onto rack tabs for secure, visible storage
- Spacious enough to fit most screw-top micro-centrifuge tubes
- Stackable, saving valuable space



Pack of 5


Item No.	Color	Item No.	Color
HS2345A	Assorted*	HS2345D	Green
HS2345B	Natural	HS2345E	Pink
HS2345C	Blue	HS2345F	Yellow

*Assorted colors (Blue, Green, Yellow, Orange, Pink)

Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
1.5/2.0	96	11	9.6 x 4.7 x 2	24.4 x 12 x 5.1
0.5		7.6		

Easy Identification

Alphanumeric well identification



Rota-Rack® Duo Tube Rack

Dilution and concentration studies made easy

Transfer solutions between 0.2 mL PCR tubes and 0.5 and 1.5 mL microcentrifuge tubes

- Purple section holds 48 x 0.5 mL and 48 x 1.5/2.0 mL microcentrifuge tubes on the reverse side
- Blue section holds 96 x 0.2 mL tubes and 48 x 1.5/2.0 mL tubes on the reverse side



Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS29041A	Blue/	0.2	96	5	11.6 x 4.6 x 1.3	29.6 x 11.8 x 3.4
	Purple	0.5	48	8		
		1.5/2.0	96	11		

Two Tier Microtube Rack 20-Well

Interlock multiple racks to increase tube storage


Improved Experience

- Tiered rack allows for improved staging and visibility for easier, more accurate sample filling and retrieval
- Made of easy to clean polypropylene
- Easily remove tubes using finger slot on front of the rack
- Embossed Index Numbers for quick sample ID and retrieval



Pack of 6

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS4284A	Blue	1.0	20	10	8.3 x 3.5 x 2	21.1 x 9 x 5
		1.5				
		2.0				



Microtube Racks 80-Well

Popular fraction-collector-plate style rack

Compact 5x16 array makes filling and transfer easy

- Frosted labeling area on bottom of rack
- Embossed index numbers for quick sample ID and retrieval
- Made of polypropylene



Pack of 5

Item No.	Color	Part No.	Color
HS29025A	Natural	HS29025E	Yellow
HS29025B	Blue	HS29025F	Orange
HS29025C	Green	HS29025G	Assorted*
HS29025D	Pink	HS29025J	Red

*Assorted colors (Blue, Green, Yellow, Orange, Pink)

Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
1.5/2.0	80	11	8.9 x 2.6 x 1	22.5 x 6.7 x 2.8



Storage Racks for 80-well Racks

Strong, clear acrylic construction

These racks help you keep large numbers of microtubes organized on your benchtop.

- Holds up to eight 80-well microtube racks filled with either 1.5 or 2.0 mL microtubes
- Plenty of room to slide the racks in and out

Item No.	Description	Color	L x W x H in	L x W x H cm
HS23263A	Vertical	Clear	6.2 x 9.3 x 9.2	15.7 x 23.6 x 23.4
HS23263B	Horizontal	Clear	12.2 x 9.3 x 4.8	31 x 23.6 x 12.2

Microtube Racks

Durable, Simple to Use

- Made from reinforced nylon offering superior chemical and heat-resistance over most other plastic based racks
- Handles provide easy and safe transportation
- Stable design can be used in water baths
- Alphanumeric identification for quick location and retrieval of tubes



Item No.	Color	Tubes mL	Wells	Well Dia. mm	Packaging	L x W x H in	L x W x H cm
120241	Blue	1.5/2.0	24	11.1	8	5.6 x 3.5 x 2	14.3 x 9 x 5.2
120243			96		4	8.3 x 5.2 x 2.1	21 x 13.1 x 5.4



Learn more about our convenient multi-functional working rack solutions:


Work 2 Store® Expanding Storage Box

Single rack, double function

Streamline and reduce sample handling time with this convenient working and storage rack, eliminating the need several racks when working with sample for cold temperature storage

Holds 64 x 1.5/2.0 mL cryo and/or microtubes

For more details and ordering information see page 92



TracRack® Tube Rack


Take what you need and store the rest

Comprised of a storage box and eight removable tube segments allowing you to easily take what you need - an entire segment or individual sample - and store the rest!

Convenient Rack System

- Eight polypropylene segments hold 1.5/2.0 mL cryo and microcentrifuge tubes (64 tubes)

For more details and ordering information see page 93





Alphanumeric Identification
for quick location and retrieval of tubes

Microtube Storage Boxes

- Hold 1.5 and 2.0 mL microcentrifuge tubes and cryogenic vials




Stylish & Identifiable

- Made of translucent, vibrantly colored polypropylene allows for visibility inside box and color coding identification
- Tube capacity is indicated on the outside the box within a distinctive wave design
- Quick visual identification with uniform wells and alphanumeric location identification along with a frosted labeling area

Robust

- Easy open and secure lid fitted with either friction fit or snap clasp designs
- Stackable, saving valuable space in the lab



50-Well Microtube Storage Boxes

Pack of 5

Item No.	Color
120032	Natural
120033	Assorted*

*Assorted colors (Blue, Green, Purple, Yellow, Orange)

Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
1.5/2.0	50	12.3	5.6 x 3.6 x 2.2	14.2 x 9.2 x 5.6



81-Well Microtube Storage Boxes

- Black, light blocking version available for light sensitive samples

Pack of 5

Item No.	Color	Item No.	Color
120034	Natural	120038	Yellow
120035	Blue	120039	Orange
120036	Green	120040	Assorted*
120037	Pink	120041	Black

*Assorted colors (Blue, Green, Yellow, Orange, Pink)

Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
1.5/2.0	81	12.3	5.1 x 5.1 x 1.9	13 x 13 x 4.7



100-Well Microtube Storage Boxes

- Black, light blocking version available for light sensitive samples

Pack of 5

Item No.	Color	Item No.	Color
120042	Natural	120202	Green
120043	Black	120203	Yellow
120044	Assorted*	120204	Orange
120201	Blue		

*Assorted colors (Purple, Blue, Green, Yellow, Orange, Pink)

Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
1.5/2.0	100	12.3	5.6 x 6 x 2.2	14.1 x 15.1 x 5.7

Rota-Rack® Tube Racks

Get organized with this interlocked rack

Each module is interlocked and adjustable with just a twist allowing you to work with multiple sized tubes in the same rack

- Suitable for variety of culture, centrifuge, micro-centrifuge, round, flat or conical bottomed tubes
- Two versions to accommodate small and large tube sizes





Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS29040A	Yellow	50	2	29	8.8 x 4.3 x 4.3	22.4 x 11 x 11
	Green	10-15	8	20		
	Blue	15-20	6	16		
	Purple	10-13	10	13.5		
HS29040B	Yellow	15	6	16.4	8.2 x 3.9 x 3.9	21 x 10 x 10
	Green	1.5/2.0	9	11.2		
	Blue	0.5	12	8		
	Purple	0.2	32	6.5		



Rotate and Customize
to suit your application needs

Version available for PCR and Microtubes:

Rota-Rack® Duo Tube Rack



Rota-Rack, Page 72

Cube Tube Rack

Interlocking multi-tube racks

- Made of polypropylene
- Allows you to work with different sized tubes all in one unit



Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS29050A	Green	50	4	30	4.3 x 4.3 x 4.5	10.8 x 10.8 x 11.3
	Yellow	15	10	17		
	Pink	3-7	12	12.5		
	Orange*	1.5/2.0	16	11.5		

*Assorted colors (Blue, Green, Pink, Yellow, Orange)

Pack of 5



4-Way Tube Rack

Link together to meet your application needs

Versatile

- Use as standalone unit or infinite units can be linked to create configurations to fit number of wells and tube sizes needed

Robust Functional Design

- Molded alphanumeric identification for quick location and retrieval of tubes
- Made from polypropylene that is chemical resistant and easy to clean
- Racks can be used in water baths and placed in freezer



Pack of 5

Item No.	Color	Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS29022A	Natural	HS29022E	Yellow	50	4	30	6.9 x 3.7 x 2.1	17.4 x 9.5 x 5.2
HS29022B	Blue	HS29022F	Orange	15	12	17.5		
HS29022C	Green	HS29022G	Assorted*	0.5	32	8.2		
HS29022D	Pink	*Assorted Colors (Blue, Green, Pink, Yellow, Orange)			1.5/2.0	32	11.6	

Check out our Clinical Version

Accommodates all brands of 13 and 17 mm round, flat and conical bottom blood collection & culture tubes



Transformer Cube Tube Rack

Low temp storage of large samples,
15 mL & 50 mL tubes

- Quick change over with 2 snap-in inserts that hold 15 and 50 mL tubes
- Rack keyed and numeric wells for tube identification
- Clear Lid included for viewable, stackable storage
- Open space available under insert for extra storage



Pack of 5

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H Dia.	L x W x H cm
HS23271A	Purple	15	16	18	4.9 x 4.9 x 5	12.5 x 12.5 x 12.9
		50	9	31		

OneRack® Multi

Exclusive design, optimal functionality

The most flexible tube rack designed for use with a
multitude of tube diameters and heights

Holds virtually any brand tube 10-30 mm of
varying heights

Multiple Configurations for Workflow Efficiency


- Interchangeable tube rack inserts hold multiple tube sizes firmly in place when alternating between sample vessels
- Inserts can be pivoted in either an upper or lower position to securely accommodate the size and height of tube being used

Multi-Functional Design

- Open design allows for better filling visibility of tube graduations when filling
- Rack wells have a unique teardrop design that enables better retention of tubes over traditional designs
- Rotate tube rack insert to an angle position to facilitate more comfortable pipetting

Unique Well Design

Tear drop shape allows
for better retention of
tubes



Lower position is ideal height
when working with shorter
microcentrifuge tubes

Alphanumeric grid

for quick sample identification
and retrieval


Ergonomic handles

for ease of transporting
the rack




Angled position for
more comfortable pipetting

Included:
1 Base + 6 Tube
Rack Inserts



Two 10-13 mm Tube Racks Fits
5 mL, 7.5 mL, 10 mL Tubes



Two 13-17 mm Tube Racks
Fits 5 mL, 15 mL Tubes



Two 25-30 mm Tube Racks
Fits 25 mL, 50 mL Tubes

OneRack® Multi

US Patent No. 9776188

Item No.	Color	L x W x H in	L x W x H cm	# Inserts Included	Well Dia./Insert	Wells/Insert
120817	Blue/Natural	12.2 x 4.2 x 3	30.9 x 10.7 x 7.7	6 (2 ea/ Diameter)	13 mm	20
					17 mm	20
					30 mm	6
120818	Purple/Natural	12.2 x 4.2 x 3	30.9 x 10.7 x 7.7	6 (2 ea/ Diameter)	13 mm	20
					17 mm	20
					30 mm	6

























Adapt-a-Rack® Adaptable Multi-tube Rack

Each well has fingers that flex that adapt to accommodate a variety of tubes sizes while holding them in a firm upright position

- Made of POM
- Racks can be linked together for additional capacity
- Open sides allow for label viewing
- Flat base to hold tubes upright
- Available in solid or a combination of colors for easy coding identification



Pack of 2

Item No.	Rack #1 Color		Rack #2 Color	
120184	Blue/Green	 	Green/Blue	 
120185	Blue/Yellow	 	Yellow/Blue	 
120186	Pink/White	 	White/Pink	 
120187	Blue/Blue	 	Blue/Blue	 
120188	White/White	 	White/White	 
120189	Pink/Pink	 	Pink/Pink	 
120190	Yellow/Yellow	 	Yellow/Yellow	 
120191	Green/Green	 	Green/Green	 

Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
7-50	4	12-30	7.1 x 2.2 x 3	18.1 x 5.6 x 7.6

Looking for floating racks for use in water or ice baths:




One-Well Connecting Tube Racks

Form a line or use as a single, stand-alone polypropylene rack

- Open Design enables easy identification of samples
- Molded circular base ring keeps conical or round-bottom tubes in place

Pack of 5

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS23052		50	1	29	3.2 x 2.8 x 2.6	8.3 x 7 x 6.7

*Assorted colors (Blue, Green, Yellow, Orange, Pink)



One-Well Friction-Fit Tube Rack

- Made from polystyrene
- Open Design enables easy identification of samples
- Holds conical and round bottom tubes
- Slides smoothly across the benchtop

One-Well Wire Tube Rack

- HDPE-coated steel wire
- Holds conical and round bottom tubes
- Open Design enables easy identification of samples
- Base ring keeps tubes in place



One-Well Friction-Fit Tube Rack

Pack of 5

Item No.	Color	Tubes mL	Wells	Well Dia. mm	W x H in	W x H cm
HS23053A	Blue	50	1	29	3.6 X 1.2	9.2 x 3

One-Well Wire Tube Rack

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS23051	Yellow	50	1	29	3.5 x 2.6 x 2.6	9 x 6.7 x 6.7




Snap-Together Conical Tube Rack

Space Saver - shipped flat and easy-to-assemble

Conveniently holds 15 and 50 mL tubes simultaneously

- ABS racks snap together securely and are designed for stable use in water baths
- Alphanumeric identification for quick location and retrieval

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in (cm)
HS24306	Orange	15	30	17	8.1 x 6.7 x 2.4
		50	20	28	(17.1 x 17 x 6.1)



Universal Tube Rack

Ideal for holding a wide range of tubes

- Made of polypropylene
- Holds up to 12 similar size tubes
- Alphanumeric identification for quick location and retrieval of tubes
- Side handles for easy transport

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS24330A	White/Purple	15-50	12	16-30	7.6 x 5.5 x 2	19.4 x 14 x 5.1

Pop-Up Tube Racks

Collapse for compact storage and “pop-up” for quick use

Designed to save space!

- Easy to open, holds 15 and 50 mL tubes
- Easy to close, stackable for compact storage
- Alphanumeric well identification for quick location and retrieval of samples

Pack of 2

Item No.	Color	Tubes mL	Wells	Well Dia. mm
HS24320B	Blue	15, 50	21, 12	17, 30
HS24320P	Pink			
HS24318	Purple	15	45	17
HS24319	Green	50	18	30
L x W x H (Expanded) in	L x W x H (Collapsed) in	L x W x H (Expanded) cm	L x W x H (Collapsed) cm	
10 x 5.4 x 2.8	10 x 5.4 x 0.8	25.5 x 13.7 x 7.2	25.5 x 13.7 x 2.1	



Save Space

Rack easily “Pops” open for easy use



Fold and Snap Tube Rack

Economical durable rack

- Polypropylene racks ship flat in one piece
- Easy-to-assemble, just fold and snap together
- Stable design can be used in water baths
- Open design allows samples to be incubated at the same temperature
- Color code for easy identification




Easy to assemble


Just fold and snap together!




90-Well



60-Well



40-Well



Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS243073B	Blue	5–10	90	13	9.6 x 4.1 x 2.4	24.5 x 10.5 x 6.3
HS243073W	White					
HS243073R	Red					
HS243073Y	Yellow	5–15	60	16	9.6 x 4.2 x 2.8	24.5 x 10.8 x 7.2
HS243077B	Blue					
HS243077W	White					
HS243077R	Red	5–30	40	20	9.8 x 4.2 x 2.8	24.9 x 10.8 x 7.2
HS243077Y	Yellow					
HS243071B	Blue					
HS243071W	White					
HS243071R	Red					
HS243071Y	Yellow					

OneRack® Tube Rack

Tough and highly resistant to acids, bases, solvents and heat

Durable racks are molded in a single, continuous piece so no assembly required and no detachable pieces

Intuitive Design

- Unique tear drop shaped openings offer better retention of tubes
- Flat labeling area for easier identification
- Stable design for use with water baths
- Molded in coordinates for tube orientation

Available in two materials:

Fiberglass-Reinforced Polyoxymethylene (POM)

- Available in half and full rack sizes for 5 different tube diameters
- Offers the highest chemical and heat resistance (-40 to +121°C)

Polypropylene

- Available in full rack sizes for 2 different tube diameters
- Offers good chemical and heat resistance (0 to 110°C)

OneRack® Full-Size POM



Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in (cm)
HS27511A	White	5-10	72	13	7.9 x 4 x 2.3 (20.2 x 10.4 x 5.9)
HS27511B	Blue				
HS27511C	Green				
HS27511D	Magenta				
120205	Yellow				
HS27512A	White	5-10	72	16	9.8 x 5 x 2.7 (25 x 12.7 x 7)
HS27512B	Blue				
HS27512C	Green				
HS27512D	Magenta				
120206	Yellow				
HS27513A	White	10-18	40	20	9.9 x 3.9 x 3.2 (25.2 x 10 x 8.3)
HS27513B	Blue				
HS27513C	Green				
HS27513D	Magenta				
HS27514A	White	10-18	40	25	11.8 x 4.7 x 3.6 (30 x 12 x 9.2)
HS27514B	Blue				
HS27514C	Green				
HS27514D	Magenta				
HS27515A	White	15-50	24	30	11.1 x 4.3 x 3.3 (28.2 x 11 x 8.5)
HS27515B	Blue				
HS27515C	Green				
HS27515D	Magenta				
120207	Yellow				

OneRack® Half-Size POM



Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in (cm)
HS27501A	White	5-10	36	13	4.1 x 4.1 x 2.3 (10.4 x 10.4 x 5.9)
HS27501B	Blue				
HS27501C	Green				
HS27501D	Magenta				
120398	Yellow				
HS27502A	White	5-10	36	16	5 x 5 x 2.7 (12.7 x 12.7 x 7)
HS27502B	Blue				
HS27502C	Green				
HS27502D	Magenta				
120399	Yellow				
HS27503A	White	10-18	20	20	4 x 5 x 3.2 (10 x 12.7 x 8.3)
HS27503B	Blue				
HS27503C	Green				
HS27503D	Magenta				
HS27504A	White	10-18	16	25	4.7 x 4.8 x 3.6 (12 x 12.2 x 9.2)
HS27504B	Blue				
HS27504C	Green				
HS27504D	Magenta				
HS27505A	White	15-50	9	30	4.3 x 4.3 x 3.3 (11 x 11 x 8.5)
HS27505B	Blue				
HS27505C	Green				
HS27505D	Magenta				


OneRack® Full-Size Polypropylene

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in (cm)
HS27551A	White	5-10	72	13	7.9 x 4 x 5.9 (20.2 x 10.4 x 5.9)
HS27551B	Blue				
HS27551C	Green				
HS27551D	Magenta				
HS27552A	White	5-10	72	16	9.8 x 5 x 2.7 (25 x 12.7 x 7)
HS27552B	Blue				
HS27552C	Green				
HS27552D	Magenta				

Coated Wire Racks



Chemically resistant HDPE Coated Wire Racks


- Holds all brands and types of conical, round or flat bottom, glass or plastic, tubes and/or bottles
- Square wells and three deck levels help keep tubes secure



***120089**
6x12, 72 Wells







Fits 10-13 mm Tubes

Item No.	Color		Tube Array	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120089*	Blue		6x12	72	13	8.1 x 4.1 x 2.6	20.5 x 10.5 x 6.5
120759	Blue		9x12	108	13	8 x 6.4 x 2.6	20.4 x 16.1 x 6.5



***120761**
4x12, 48 Wells



Fits 13-16 mm Tubes


Item No.	Color		Tube Array	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120760	Blue		4x10	40	16	7.9 x 3.3 x 2.4	20 x 8.3 x 5.9
120761*	Blue		4x12	48	16	9.4 x 3.3 x 2.4	23.7 x 12.7 x 5.9
HS23072	Blue		6x12	72	16	9.5 x 5 x 2.6	24.1 x 18.8 x 6.5
120763	Blue		9x12	108	16	9.3 x 5 x 2.6	23.7 x 12.7 x 6.5
HS23108	Red		9x12	108	16	9.3 x 5 x 2.6	23.7 x 12.7 x 6.5

Fits 16-18 mm Tubes – 1 5mL


Item No.	Color	Tube Array	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS23215*	Blue	3x5	15	20	5 x 3 x 3	12.7 x 7.6 x 7.6

Fits 16-20 mm Tubes – 15 mL

Item No.	Color	Tube Array	Wells	Well Dia. mm	L x W x H in	L x W x H cm	
120764	Blue		4x10	40	20	9.5 x 3.9 x 3.1	24 x 9 x 8
120765	Blue		8x10	80	20	9.5 x 7.6 x 3.1	24 x 19 x 8




***HS23215**
3x5,
15 Wells







***120768**
8x10, 80 Wells

Fits 20-24 mm Tubes

Item No.	Color		Tube Array	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120766	Blue		4x10	40	24	10.2 x 4.2 x 3.2	26 x 10.7 x 8

Fits 20-25 mm Tubes

Item No.	Color		Tube Array	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120769	Blue		3x8	24	25	10.7 x 4.1 x 3.2	27.3 x 10.5 x 8
120767	Blue		4x10	40	25	11.4 x 4.7 x 3.1	29 x 11.8 x 8
120768*	Blue		8x10	80	25	11.4 x 9.1 x 3.1	29 x 23 x 8

Coated Wire Racks

Chemically resistant HDPE Coated Wire Racks



***HS232461**
2x4, 8 Wells

Fits 25-30 mm Tubes - 50 mL


Item No.	Color	Tube Array	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS232461*	Blue	2x4	8	30	5.8 x 3.2 x 3.5	14.7 x 8 x 8.9
120770	Blue	6x8	48	30	10.9 x 8.3 x 3.9	27.7 x 20.9 x 9.8

Fits 30-40 mm Tubes

Item No.	Color	Tube Array	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120771*	Blue	3x8	24	40	14.1 x 5.5 x 3.9	35.7 x 13.9 x 9.8
120772	Blue	6x8	48	40	14 x 10.6 x 3.8	35 x 26 x 9.8



***120771**
3x8, 24 Wells




***120086**
5x10, 50 Wells

Fits 23-25 mm Bottles - Bijou Bottles


Item No.	Color	Tube Array	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120085	Blue	2x6	12	25	6.4 x 2.5 x 1.3	16.4 x 6.4 x 3.2
120086*	Blue	5x10	50	25	10.7 x 5.7 x 1.3	27.1 x 14.5 x 3.2

Fits 30-33 mm Bottles - Universal & McCartney Bottles

Item No.	Color	Tube Array	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120087*	Blue	2x6	12	33	8.3 x 3.1 x 2.4	21.2 x 8 x 6
120088	Blue	5x10	50	33	8.7 x 3.2 x 2.4	21.8 x 8.2 x 6



***120087**
2x6, 12 Wells



***120091**
4x6, 24 Wells

Fits 55-58 mm Bottles - Urine Containers

Item No.	Color	Tube Array	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120091*	Blue	4x6	24	58	15.5 x 10.3 x 2	39.3 x 26.2 x 5

Fits 55-60 mm Bottles - Centrifuge Bottles

Item No.	Color	Tube Array	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS23224*	Blue	2x3	6	60	8.1 x 5.3 x 4.1	20.5 x 13.5 x 10.4




***HS23224**
2x3, 6 Wells

Mega Rack®

High-capacity storage for your tubes

Fully assembled, polypropylene mega racks are great for high throughput labs requiring large sample storage

- Sturdy, material for easy and secure transport of blood collection and other tubes
- Alphanumeric identification for quick location and retrieval of tubes
- Place different test groupings into one rack eliminating the need for multiple racks



Mega Racks® 10-13 mm - Holds 5-7 mL sample tubes

Item No.	Color	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS24322A	Blue	216	13	13.3 x 9.3 x 3.6	33.7 x 23.5 x 9.1
HS24322E	White				
120112	White	432	13	26.5 x 9.3 x 3.5	67.3 x 23.5 x 9.0
120113	Blue				

Mega Racks® 13-16 mm - Holds 10 mL sample tubes

Item No.	Color	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS24323A	Blue	120	16	13.3 x 10.5 x 3.6	33.7 x 26.7 x 9.1
HS24323E	White				
120114	White	240	16	26.5 x 10.5 x 3.5	67.3 x 26.7 x 9.0
120115	Blue				

Mega Racks® 18-20 mm - Holds 15 mL conical/swab tubes and 20 mm round bottom tubes

Item No.	Color	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120834	White	120	20	14.9 x 11.9 x 3.6	37.7 x 30.1 x 9.1
120835	Blue				



Peg Tube Racks


Strong, tough, versatile racks for bench use or tube storage

Fiberglass-reinforced polypropylene racks are ideal for bench work or tube storage

- Spacious peg design allows easy viewing of tube contents
- Pegs remain straight even after repeated autoclaving
- Embossed alphanumeric identification for quick location and retrieval of tubes





Pack of 2



Unique tear drop well design
allows for better retention of tubes


Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS24311A	White	10	50	14-17	7.4 x 4.1 x 2.7	18.8 x 10.4 x 6.9
HS24311B	Blue					
HS24312A	White	7	96	10-13	3.2 x 5 x 2.7	8.2 x 12.6 x 6.9
HS24312B	Blue					

View our complete line of Freezer Storage Boxes



Cold Products & Storage

From cool-looking ice containers to freezer necessities, our innovative range of cold products continues to expand



Cooling Containers

Low Temp Racks & Storage

Freezer Boxes, Racks & Accessories



Cool Containers

Durable containers for chilling temperature-sensitive samples

Suitable for low temperature use. Ideal for ice, ice slurries and CO₂ pellets


- Made from polyurethane with metal inserts for stability
- Close-fit lid keeps evaporation to a minimum and features a comfortable recessed knob
- Tapered design for easy stacking
- Precise pour spout and large handles



Lid nests under bucket
for convenience and to save space

Item No.	Color		Style	Capacity liters	L x W x H in	L x W x H cm
HS28714K	Black	■	Bucket	4	7.1 (H) x 9.2 (D)	18.1 (H) x 23.4 (D)
HS28721G	Green	■	Pan	1	10.5 x 7.7 x 4.5	26.7 x 19.6 x 11.4
HS28724B	Blue	■	Pan	4	15 x 10 x 5.5	38.1 x 25.4 x 14
HS28729R	Red	■	Pan	9	20.5 x 11.5 x 6.4	52.1 x 29.2 x 16.3

Learn more about our open rack systems that work perfectly with our Cooling Containers



Floating Microtube Rack

For use on bench and in ice and water baths

- Made of polypropylene
- Attach legs and rack stands upright for use on bench top
- Detach legs and samples can be quickly removed by pushing down onto any surface



Pack of 4

Item No.	Color	Tubes mL	Wells	Well Dia. mm	Shape/Dia. in
HS2132A	White	1.5/2.0	16	10.9	Square 4
HS2134A			20		Round 3.8

Floating Tube Racks

Reusable and flexible

- Ideal for floating tubes in water baths, tubs or beakers
- Brightly colored foam racks made of high-density polyethylene that resists water absorption
- Racks include detachable carrying handles to easily lift samples out of baths
- Temperature Range 0° to 70° C



Pack of 5

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS2165A	Green	50	4	29	6.4 x 5.3 x 2.5	16.5 x 13.7 x 6.5
HS2165B	Blue	15	8	16	6.7 x 2.2 x 2.5	17.2 x 5.7 x 6.5
HS2165C	Yellow	1.5/2.0	24	9.3	6 x 3.2 x 2.7	15 x 8.3 x 7
HS2166	Blue	0.2/0.5/1.5/2.0	18	9.3	3.2	8.2 dia.



Alphanumeric Well Identification

for quick location and retrieval of tubes



Lid doubles as a stand
enables the rack to stay cold longer as there is no conductive heat transfer from surface beneath



Vented Lid

enables faster cooling of samples while in freezer

PCR Cooler, 96-Well

Ice-Free Cooling System for PCR Preparation and Storage

Cools industry standard format PCR tubes, strips and single format as well as most standard and low profile 96-Well plates

Consistent & Reliable

- Modern design cools faster in freezers and stays cold longer on the bench top
- Maintains the entire 96-Well plate at a temperature of 0 °C for over an hour
- Maintains sample temperature of less than 5°C for up to three hours (with lid off)
- When the temperature ceiling of 6°C has been exceeded, users are signaled by a contrasting color change of the cooler or a single well




Ergonomic Textured Finger Grips

provide a secure hold, reducing the chance of dropping the PCR rack, and makes handling and transportation easy, especially when using gloves

Two Color Contrasting Options Available

purple to pink or dark blue to light blue
allowing easy monitoring of samples



				US Patent Pending, Community Design No: 004101756-0001-0002			Pack of 2
Item No.	Warm Color	Cold Color	Wells	Tube Sizes mL	L x W x H in	L x W x H cm	
120727	Pink	Purple	96	0.2 Assay Plates	4.3 x 5.8 x 1.9	11 x 14.8 x 4.8	
120728	Light Blue	Dark Blue					
120729	Assorted*						

*Assorted colors (Light Blue/Dark Blue & Pink/Purple)

Need to keep your reagents cool while working?

Check out our ice-free Reagent Reservoir Chillers with temperature color changing feature



Reservoir Chiller, Page 50

Gridded and Keyed
one-position lids for easy
sample identification




0.5 mL Tube Adapters
supplied as standard with the
0.5–2.0 mL tube size versions

Mini Coolers


Protect temperature-sensitive samples at the bench or in the freezer

Ideal for keeping samples or reagents cool during testing


- Polycarbonate construction with rubber non-slip feet
- Coolers 120072 and 120075 are filled with a self-insulating gel to improve cooling capacity of unit



Tube Array
Options



8 x 12



Item No.	Color	Tubes mL	Lid	Tube Array	Wells	L x W x H in	L x W x H cm	Temp.	Temp. Retention
120070	Green	0.5–2.0	Clear	3 x 4	12	5.7 x 4.1 x 3.9	14.5 x 10.5 x 10	0°C	1°C for 3 hr
120071			Clear	4 x 8	32	9.6 x 4.1 x 3.9	24.5 x 10.5 x 10	0°C	1°C for 2 hr
120072			Gel Filled	4 x 8	32	9.6 x 4.1 x 3.9	24.5 x 10.5 x 10	0°C	1°C for 5 hr
120073	Blue	0.5–2.0	Clear	3 x 4	12	5.7 x 4.1 x 3.9	14.5 x 10.5 x 10	-20°C	-15°C for 1 hr
120074			Clear	4 x 8	32	9.6 x 4.1 x 3.9	24.5 x 10.5 x 10	-20°C	-15°C for 1 hr
120075			Gel Filled	4 x 8	32	9.6 x 4.1 x 3.9	24.5 x 10.5 x 10	-20°C	-15°C for 1 hr
120076	Blue	0.2	Clear	8 x 12	96	8.9 x 6.5 x 4.3	22.5 x 16.5 x 11	-20°C	-15°C for 2 hr



Work 2 Store® Expanding Storage Box

Single rack, double function

Streamline and reduce sample handling time with this convenient working and storage rack, eliminating the need several racks when working with samples requiring cold temperature storage

Efficient Handling

- Expand rack to opening position instantly creating more working space.
- Enough room for sample manipulation with no interference from open microfuge lids
- Collapse rack into closed position when it's ready for freezer storage
- Compact design will fit into standard 53 mm freezer racking system
- Alphanumeric identification for quick location and retrieval of tubes

Safe and Stackable Storage

- Clear lids locks the rack together for safe handling and transportation
- Low profile, easily stackable

Holds 64 x 1.5/2.0 mL cryo and/or microtubes

- Expanding wells have a center depression to stabilize conical tubes and a cryogenic vial lock for one handed cap removal



Work 2 Store® Expanding Tube Storage Box

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H (Open) in	L x W x H (Closed w/lid) in	L x W x H (Open) cm	L x W x H (Closed w/lid) cm
120018	White	1.5/2.0	64	12.3	12.2 x 5 x 1.3	5.2 x 5.2 x 2	31 x 12.8 x 3.4	13.3 x 13.3 x 5.3


TracRack® Tube Rack

Take what you need and store the rest

Comprised of a storage box and eight removable tube segments allowing you to easily take what you need - an entire segment or individual sample - and store the rest!

Convenient Rack System

- Eight polypropylene segments hold 1.5–2.0 mL cryo and microcentrifuge tubes (64 tubes)
- Wells feature a center depression stabilizing conical tubes and a cryogenic vial lock for one-handed cap removal
- Printed index letters on each segment make for quick sample ID and retrieval




TracRack® Outer Box

- Made from thin film polypropylene material
- Moisture and mold resistant
- 2 in (5.3 cm) boxes fit in standard freezer racks
- Not recommended for autoclaving

TracRack® Tube Rack

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120356	Purple	1.5/2.0	64	12.3	5.1 x 5.1 x 1.4	12.9 x 12.9 x 3.5



Low-Temp PCR Rack: 96-Well

Slide, set, work & store PCR samples

Multi-purpose rack converts from compact workstation to lidded storage solution accommodating individual and strip 0.2 mL tubes

Versatile & Robust Solution

- Polycarbonate lid serves as cover and elevated rack base - just slide open and pivot lid under rack for compact workstation
- Lid can be used in right or left-handed orientation and is detachable
- Stackable saving valuable bench-top and freezer space
- Multiple colors offered for coding options



US Patent Pending, Community Design No: 004101756-0001-0002

Pack of 5

Item No.	Color	Item No.	Color
120538	Assorted*	120540	Green
120539	Blue	120541	Purple

*Assorted colors (Purple, Blue, Green, Orange, Yellow)



Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
0.2	96	5.2	4.9 x 3.6 x 1.5	12.5 x 9.2 x 3.7



Alphanumeric Identification
for quick location and retrieval of tubes

Microtube Storage Boxes

- Hold 1.5 and 2.0 mL microcentrifuge tubes and cryogenic vials



Stylish & Identifiable

- Made of translucent, vibrantly colored polypropylene allows for visibility inside box and color coding identification
- Tube capacity is indicated on the outside the box within a distinctive wave design
- Quick visual identification with uniform wells and alphanumeric location identification along with a frosted labeling area

Robust

- Easy open and secure lid fitted with either friction fit or snap clasp designs
- Stackable, saving valuable space in the lab



50-Well Microtube Storage Boxes

Pack of 5

Item No.	Color
120032	Natural
120033	Assorted*

*Assorted colors (Blue, Green, Purple, Yellow, Orange)

Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
1.5/2.0	50	12.3	5.6 x 3.6 x 2.2	14.2 x 9.2 x 5.6



81-Well Microtube Storage Boxes

- Black, light blocking version available for light sensitive samples

Pack of 5

Item No.	Color	Item No.	Color
120034	Natural	120038	Yellow
120035	Blue	120039	Orange
120036	Green	120040	Assorted*
120037	Pink	120041	Black

*Assorted colors (Blue, Green, Yellow, Orange, Pink)

Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
1.5/2.0	81	12.3	5.1 x 5.1 x 1.9	13 x 13 x 4.7



100-Well Microtube Storage Boxes

- Black, light blocking version available for light sensitive samples

Pack of 5

Item No.	Color	Item No.	Color
120042	Natural	120202	Green
120043	Black	120203	Yellow
120044	Assorted*	120204	Orange
120201	Blue		

*Assorted colors (Purple, Blue, Green, Yellow, Orange, Pink)

Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
1.5/2.0	100	12.3	5.6 x 6 x 2.2	14.1 x 15.1 x 5.7

Transformer Cube Tube Rack

Low temp storage of large samples,
15 mL & 50 mL tubes

- Quick change over with 2 snap-in inserts that hold 15 and 50 mL tubes
- Rack keyed and numeric wells for tube identification
- Polypropylene construction capable of handling extreme storage temperatures
- Clear Lid included for viewable, stackable storage
- Open space available under insert for extra storage



Fill with ice
and use as a temporary low
temperature workstation

Pack of 5

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS23271A	Purple ■	15	16	18	4.9 x 4.9 x 5	12.5 x 12.5 x 12.9
		50	9	31		



Foam Tube Racks

High storage, minimal footprint

- Two-piece racks are constructed of expanded polystyrene (EPS) foam
- Lightweight microtube racks are ideal for long-term studies or tube storage
- Extra spacing between wells allows for easy gripping with fingertips
- Racks stack securely with nesting features on lid and base
- Must use tubes capable of being stored at the required temperature of the user



Pack of 100

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS4282	White □	1.5/2.0	100	11.5	13.2 x 3.7 x 2.9	33.6 x 9.5 x 7.3
HS4283P	White □	1.5/2.0	50	11.5	8.3 x 4.3 x 2.7	21 x 11 x 7.1

Store microslides in low temperatures with our Low Temperature Microscope Slide Boxes



Cryogenic Vial Rack

Safely store and transport cryogenic vials

- Made from polycarbonate
- Well tabs allow one-handed opening and closing of vials that have a slot feature
- Embossed index numbers and squared, upper left corner provide easy orientation
- Drainage holes in each well aid in cleaning
- Alphanumeric identification for quick location and retrieval of tubes



Embossed Index Numbers and squared, upper left corner provide easy orientation



Ribbed Grips make handling with gloves easier

Pack of 4

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
HS23229B	Blue	2-5	50	12.9	7.8 x 4.1 x 0.9	19.9 x 10.5 x 2.8

Cryogenic Vial Canes

Organize cryogenic samples in dewars


Aluminum canes hold vials for freezer storage or immersion in LN₂ dewars

- Use in liquid and vapor phase; use approved vials only
- Canes bend to fit all major cryogenic vial brands



Pack of 12

Item No.	Tubes mL	No. of Vials	Length in	Length cm
HS23205	1.5/2.0	5	11.4	29
HS23206	1.5/2.0	6	11.8	30



PVC Cryogenic Cane Sleeves

For extra security and sample safety

Clear plastic sleeves protect vials in aluminum canes during handling and storage

- Safe for cryogenic freezing
- More durable than cardboard sleeves
- Transparency allows sample identification without removing sleeve



Pack of 100

Item No.	Diameter in	Length in	Diameter cm	Length cm
HS23210	0.6	10.8	1.5	27.5



Don't miss our line of Cryogenic Vials


Our innovative flush design saves space in your freezer



Cryogenic Vial Boxes & Partitions

For storage of small samples and tubes

- Made from cardboard
- Colors allow for easy identification and sorting purposes
- All box footprints are 13.3 x 13.3 cm, select from two different heights



Vent Holes

allow faster cooling and removal of LN₂

Cardboard Cryogenic Vial Boxes

Pack of 12

Item No.	Color		Tubes mL	L x W x H in	L x W x H cm
HS2860A	White	<input type="checkbox"/>	1.5/2.0	5.2 x 5.2 x 2	13.3 x 13.3 x 5
HS2860B	White	<input type="checkbox"/>	1.5/2.0	5.2 x 5.2 x 3	13.3 x 13.3 x 7.5
HS2860C	Assorted	<input checked="" type="checkbox"/> <input checked="" type="checkbox"/> <input checked="" type="checkbox"/>			
HS2860CB	Blue	<input checked="" type="checkbox"/>	1.5/2.0	5.2 x 5.2 x 2	13.3 x 13.3 x 5
HS2860CR	Red	<input checked="" type="checkbox"/>			
HS2860CY	Yellow	<input checked="" type="checkbox"/>			

Partitions

Partitions must be ordered separately

Pack of 12

Item No.	Array	Tubes mL	Wells	Well Dia. mm
HS2860E	10 x 10	1.5/2.0	100	11.5
HS2860H	9 x 9	1.5/2.0	81	12
HS2860F	8 x 8	1.5/2.0	64	15
HS2860G	7 x 7	1.5/2.0	49	17.5

Microtube Storage Boxes

Safe microtube storage

Low profile boxes are available in two sizes to accommodate standard 0.5 mL and 0.2 mL microtubes and vessels of a similar size

- Partitions included in the boxes
- Will fit into most standard freezer racking




Pack of 10

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120220	White	<input type="checkbox"/> 0.2	169	9.5	5.2 x 5.2 x 1.1	13.3 x 13.3 x 2.9
120221	White	<input type="checkbox"/> 0.5	100	12	5.2 x 5.2 x 1.6	13.3 x 13.3 x 4.1

15/50 mL Tube Storage Boxes & Partitions

- Write-on surface for descriptions or self-numbering systems
- Partitions available in two sizes



Cardboard Tube Storage Box & Lid

Pack of 10

Item No.	Color		L x W x H in	L x W x H cm
HS2860M	White	<input type="checkbox"/>	5.7 x 5.7 x 4.8	14.5 x 14.5 x 12.2

Partitions

Partitions must be ordered separately

Pack of 10

Item No.	Array	Tubes mL	Wells	Well Dia. mm
HS2860T	6 x 6	15	36	23
HS2860P	4 x 4	50	16	35

Cryogenic Vial Tube Storage Box

- Two-piece set-up boxes with finger holds
- Boxes made from paper coated cardboard
- 5 x 5 format for tubes up to 14 mm in diameter
- Partitions are included in each box



Case of 200

Item No.	Color	Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H mm
120519	Blue	<input checked="" type="checkbox"/> 2.0	25	14	3.4 x 3.4 x 1.5	87x87x38

Freezer Boxes

Reusable and Affordable

Easy to assemble boxes made from polypropylene material

- Moisture resistant, more durable than cardboard boxes
- Dividers included with the boxes
- Corrugated Polypropylene: Store down to temperatures of -80°C with appropriate temperature viable tubes
- Thin Film: Can be used in mechanical freezers & vessels down to -196°C with appropriate temperature viable tubes
- Vapor phase only

Freezer boxes are supplied "flat-packed"

Eliminating the hassle of stacking and reduces the amount of space needed for storage

Corrugated Polypropylene 0.2/0.5/1.5–2.0/5 mL/Vial



Pack of 10

Item No.	Color		Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120379	White	<input type="checkbox"/>	0.2	144	7.9	5.2 x 5.1 x 1.2	13.3 x 13 x 3
120383	White	<input type="checkbox"/>	0.5	81	12	5.2 x 5.1 x 1.6	13.3 x 13 x 4.2
120262	White	<input type="checkbox"/>	1.5/2.0	81	12	5.2 x 5.1 x 2	13.3 x 13 x 5.1
120376	Blue	<input checked="" type="checkbox"/>					
120387	White	<input type="checkbox"/>	5	25	22.6	5.2 x 5.1 x 2.9	13.3 x 13 x 7.5
120388	Blue	<input checked="" type="checkbox"/>					
120428	White	<input type="checkbox"/>	Vial	100	17	7.4 x 7.8 x 1.6	18.8 x 19.8 x 4.1

Corrugated Polypropylene 15 mL and 50 mL



Pack of 10

Item No.	Color		Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120391	White	<input type="checkbox"/>	15	36	19.8	5.8 x 5.7 x 4.7	14.7 x 14.5 x 12
120392	Blue	<input checked="" type="checkbox"/>					
120366	White	<input type="checkbox"/>	50	16	31.6	5.8 x 5.7 x 4.7	14.7 x 14.5 x 12
120367	Blue	<input checked="" type="checkbox"/>					

Thin Polypropylene film 1.5–2.0 mL



Pack of 10

Item No.	Color		Tubes mL	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120370	White	<input type="checkbox"/>					
120371	Blue	<input checked="" type="checkbox"/>	1.5/2.0	81	12.1	5 x 4.9 x 1.9	12.5 x 12.4 x 4.9
120372	Lt. Green	<input checked="" type="checkbox"/>					



Arctic Squares® Polycarbonate Cryostorage Boxes

Safely store cryogenic vials

Our polycarbonate boxes are designed with unique features to help orient, identify and access vials with ease

- Alphanumeric labeling on lid or base (depending on the unit) for quick location and retrieval of tubes
- Liquid nitrogen vent and drainage holes
- Magnet on top of forceps enables easy storing on metal surfaces
- Stackable, saves valuable space in the lab



US Patent No.: D529,622

EP Registered Design No.: 00551429-0001, 00551429-0002, 00551429-0003

Item No.	Color	Tube Array	Tubes mL	Wells	Well Dia. mm	Lid	Packaging	L x W x H in	L x W x H cm	Forceps
HS23325A	Ruby Red	5 x 5	2.0	25	12.5	Slide On	8	3 x 3 x 2	7.6 x 7.6 x 5.3	
HS23381B	Assorted*	9 x 9	2.0	81	12.5	Hinged	4	5.2 x 5.2 x 2	13.3 x 13.3 x 5.3	✓
HS23300D	Arctic Blue	10 x 10	2.0	100	12.5	Hinged	4	5.2 x 5.2 x 2	13.3 x 13.3 x 5.3	✓

*Assorted colors (Ruby Red, Amethyst Purple, Arctic Blue, Emerald Green)



Arctic Squares® Upright Freezer Rack

Innovative horizontal drawer system for upright freezers

- Constructed of corrosion-resistant stainless steel with four drawers to hold 16 boxes
- Drawers feature a unique stop, that when fully extended pivot down allowing for easy box access
- Steel drawer handles provide a user-friendly grip
- Freezer rack also includes a pull-out wire handle for easy removal off the freezer shelf
- Compatible with most standard 2 in (5.3 cm) freezer storage boxes

Unique drawer stop that pivots down allowing for easy box access

Steel drawer handles

have a space to insert a label, and provide a user-friendly grip

Item No.	Config.	Box Height in	Box Height cm	L x W x H in	L x W x H cm
HS2864	4 x 4	2	5.3	5.4 x 22.4 x 9.3	13.9 x 57.1 x 23.8

Freezer Racks

A variety of racks to safely store samples in your freezer

Choose from a selection of corrosion-resistant stainless steel racks that are compatible with most brands of upright and chest freezers

- Racks store boxes up to a 135 x 135 mm footprint and accommodate 2" tall boxes (up to 5 cm) as well as 3" tall boxes (up to 7.6 cm) depending on the rack model



Cube Upright Freezer Rack For 15 and 50 mL tube box storage

- Constructed of corrosion-resistant stainless steel with 6 rack spaces
- Provides safe, convenient storage of your samples
- Handles at each end accommodate freezers with right or left-handed doors
- Compatible with 6 in (15.2 cm) or smaller boxes

Handles for easy transfer

Robust mounting brackets assure handles stay securely in place for easy use



Shown with our

15/50 mL Tube Storage Box, Page 97
Transformer Cube Tube Rack, Page 78

Item No.	Config.	Box Size in	Box Size cm	L x W x H in	L x W x H cm
HS28622	2 x 3	6	15.2	18 x 6.2 x 10.6	45.9 x 15.9 x 27.1

Chest Freezer Racks

Item No.	Config.	Box Size in	Box Size cm	L x W x H in	L x W x H cm
HS2862D	8 x 1	3	7.6	5.5 x 5.5 x 25.5	14.1 x 14.1 x 64.9
HS2862DB	9 x 1			5.5 x 5.5 x 29	14.1 x 14.1 x 73

Handles for easy transfer

Robust mounting brackets assure handles stay securely in place for easy use

Locking Rod

Keeps tube storage boxes safely in place



HS2862AB

Upright Freezer Racks: 2 inch Boxes

Item No.	Config.	Box Size in	Box Size cm	L x W x H in	L x W x H cm
HS2862HA	4 x 3	2	5.3	16.4 x 5.5 x 8.9	41.8 x 14.1 x 22.5
HS2862H	4 x 4			21.8 x 5.5 x 8.9	55.6 x 14.1 x 22.5
HS2862HB	5 x 4			22 x 5.5 x 11	55.6 x 14.1 x 28

Upright Freezer Racks: 3 inch Boxes






Item No.	Config.	Box Size in	Box Size cm	L x W x H in	L x W x H cm
HS2862KB	4 x 4	3	7.6	22 x 5.5 x 12.8	55.6 x 14.1 x 32.5
HS2862KC	5 x 4			22 x 5.5 x 16	55.6 x 14.1 x 40.6



HS2862H

Shown with our 81-Well Microtube Storage Boxes, Page 94

Upright Freezer Rack Configurations



Tips for Keeping your Laboratory Freezer Running Efficiently

Monitor Temperatures

Installing an alarmed thermometer monitoring system, like Heathrow Scientific's Freezer Thermometer, can alert you early to any fluctuations in temperatures that could indicate the need for cleaning or maintenance.

Organize to Ensure Proper Air Flow

Our Freezer Rack systems not only use all optimal storage space available; they also allow for air flow through the unit. With materials stacked and organized, it makes it easier to locate samples quickly, greatly reducing the amount of time the door remains open.

Remove Ice and Frost Routinely

An accumulation of ice and frost on racks and doors can interrupt air flow and overwork the freezer. Defrost on a routine basis to keep your unit in proper working order.

Clear Debris from Compressor Coils

Vacuum or use a brush to clean the compressor coils on the outside of the unit to remove any buildup of dust and grime on a regular basis to keep your unit running smoothly.

Also, be sure not to stack any items on or in front of your units venting system. This could cause the compressor to work overtime to maintain temperatures.

Clean Door Gaskets

To ensure that cold air is not escaping, clean any debris or ice build-up from the door gasket to keep an air-tight seal.

Freezer Thermometer

Protect frozen and chilled samples from potential spoilage

Monitor temperatures inside your freezer or refrigerator with the added convenience of sound and flashing LED lights to alert out-of-range temperatures

- Suction cup attaches the probe inside freezer/refrigerator while the thermometer secures to the outside with hook-and loop tape
- Large digital display shows freezer/refrigerator temperature plus ambient temperature
- Switches toggle for freezer/refrigerator monitoring and °F/°C
- Set freezer range: 14.9° to -9.4°F (-9.5° to -23°C);
Set refrigerator range: 35.6° to 46.4°F (2° to 8°C)
- Includes one AAA battery



Item No.	L x W x H in	L x W x H cm
HS28630	4.1 x .79 x 0.9	1.4 x 1.9 x 2.8

Probe Cable measures 0.9 m, Probe measures 5 mm OD x 17 mm L

Slide Products & Accessories

Products meeting the demanding standards of clinical, research and educational labs



Boxes & Storage Racks

Folders & Mailers

Holders & Staining



Sliding Lock


provides secure storage with easily identified lock and unlock symbols.

Slide Rail Design

ensures a secure fit allowing easy placement and removal of slides

Centrally Positioned Numbering

enables simple slide organization with easy viewing of location numbers



100-Place Premium Plus Slide Box

Modern design for maximum functionality

Optimal Organization

- Organize by "project" and "date" molded in on lid
- Inventory card index corresponds to slide numbers and is removable for trouble-free writing on flat surface
- Vibrant translucent colors are easy to identify and enable coding options

Functional Design

- Flexible storage, store vertically or horizontally to efficiently use space
- Molded feet allow secure horizontal stacking of multiple units
- Unique side storage can accommodate a pen, desiccant or other small items needed for working
- Hinged lid smoothly opens 180° sitting flat on surfaces and prevents misplacement

Safe & Secure

- EVA foam lines the bottom of the box helping to hold slides securely in place and protect against damage
- Stylish wave surface texture is easy to hold for safe sample transport

100-Place Premium Plus Slide Box

Pack of 5

Item No.	Color	Item No.	Color
120576	Blue	120579	Red
120577	Orange	120592	Purple
120578	Green	120593	Assorted*

*Assorted colors (Purple, Blue, Green, Red, Orange)

Lining	Lock	L x W x H in	L x W x H cm
EVA Foam	Sliding Latch	8.2 x 7.7 x 1.3	20.7 x 19.5 x 3.2

Replacement Inventory Card

Pack of 10

Item No.	L x H in	L x H cm
120618	8.5 x 11	21.6 x 28

Premium Microscope Slide Boxes

Slide boxes are perfect for storage or transportation of standard size slides (76 x 26 mm and 3 x 1 in)

Organize & Identify

- Lining and slots are numbered to correspond to the slide inventory sheet
- Organize by "Project" and "Date" with pre-molded labeling area
- Available in multiple colors for project coding

Durable & Safe

- ABS plastic boxes available with either cork or foam lining
- Designed to securely store slides in grooved slots
- Molded feet secure stacking for easy transport and storage



25-Place Slide Box

Easily access slides with friction-fit cover



Item No.	Color	Lining	Lock	L x W x H in	L x W x H cm
HS15989A	Blue	Cork	Friction Fit	5.6 x 3.5 x 1.4	14.1 x 8.8 x 3.5
HS15989B	Green				
HS15989C	Red				
HS15989E	White				
HS15989M	Blue	Foam	Friction Fit	5.6 x 3.5 x 1.4	14.1 x 8.8 x 3.5
HS15989P	Red				
HS15989R	White				

50-Place Slide Box

Thumb-latch lock mechanism



Item No.	Color	Lining	Lock	L x W x H in	L x W x H cm
HS15996A	Blue	Cork	Thumb-Latch	8.2 x 3.4 x 1.4	20.8 x 8.6 x 3.5
HS15996B	Red				
HS15996C	White				

100-Place Slide Box

Closes securely with a rust-resistant nickel-plated clasp and pin hinges



Item No.	Color	Lining	Lock	L x W x H in	L x W x H cm
HS15994A	Blue	Cork	Nickel-plated Clasp	8.2 x 6.9 x 1.3	20.8 x 17.5 x 3.4
HS15994B	Green				
HS15994C	Red				
HS15994D	Yellow				
HS15994E	White				
HS15994F	Grey				
HS15994G	Black	Foam	Nickel-plated Clasp	8.2 x 6.9 x 1.3	20.8 x 17.5 x 3.4
HS15994M	Blue				
HS15994P	Red				
HS15994R	White				

100-Place Tall Slide Box

Additional height for 3 x 2 in slides
Stackable with 100-Place Premium Slide Boxes



Item No.	Color	Lining	Lock	L x W x H in	L x W x H cm
HS15979A	Blue	Cork	Nickel-plated Clasp	8.2 x 6.9 x 2.4	20.8 x 17.5 x 6



Storage Rack for Microscope Slide Boxes

Create a file system for 50-Place and 100-Place slide boxes

Storage rack makes organization of samples easy

- Clear acrylic construction allows the labels to be seen on the slide boxes
- Place boxes in rows or columns by color to create an easy coding system

Item No.	Color	Capacity	L x W x H in	L x W x H cm
HS159890B	Clear	15 Bins	28 x 6.7 x 9	71.1 x 17.1 x 22.8

Standard Microscope Slide Boxes

Constructed of durable polypropylene



Slide boxes hold either 25 or 100 standard slides (76 x 26 mm and 3 x 1 in)




Improved Identification

- Molded slot numbers correspond directly to the slide inventory card
- Organize with "Project" and "Date" pre-molded labeling area




Secure & Stackable

- Hinged lids are held securely closed with a latch lock
- Molded feet secure stacking for easy transport and storage

25-Place Slide Box

Item No.	Color	Lining	Lock	L x W x H in (cm)
HS15990A	Blue 	None	Latch Lock	5.6 x 3.6 x 1.2 (14.1 x 9.2 x 3)
HS15990B	Red 			
HS15990C	White 			
120740	25-Place Replacement Inventory Card, Pack of 10			5.3 x 3.1 (13.4 x 8)

100-Place Slide Box

Item No.	Color	Lining	Lock	L x W x H in (cm)
HS15991A	Blue 	None	Latch Lock	8.3 x 6.7 x 1.5 (21 x 16.9 x 3.7)
HS15991B	Red 			
HS15991C	White 			
120741	100-Place Replacement Inventory Card, Pack of 10			3.9 x 6.1 (10 x 15.5)

Low Temperature Microscope Slide Boxes

Store in ultra-low temperatures with confidence



Holds standard slides (76 x 26 mm and 3 x 1 in)

Durable Construction

- Made of polycarbonate to store in ultra-low temperatures
- Uniform grooved slide slots separate slides preventing slide contact
- Molded feet enable secure stacking for easy transport and storage

Durable Construction

- Low temperature resistant numbered foam lining corresponds with identification sheet
- Organize with "Project" and "Date" pre-molded labeling area



Low Temperature Microscope Slide Boxes

Item No.	Color	Capacity	Lock	L x W x H in (cm)
120316	Blue	25-Place	Friction Fit	5.6 x 3.5 x 1.4 (14.1 x 8.8 x 3.5)
HS15988H	Purple	100-Place	Stainless Steel Clasp	8.2 x 6.9 x 1.3 (20.8 x 17.5 x 3.4)



Cardboard Slide Tray

Cardboard slide trays hold standard slides (76 x 26 mm and 3 x 1 in)

Cardboard slide folders offer a convenient solution for safely working, transporting and short-term storing of slides

Organize & Protect

- Index on the front cover corresponds with the number grid inside for identification
- Recessed slide compartments and inside covers for specimen protection
- Slides cannot be accessed when cover is closed reducing contamination

Sturdy Construction

- Compartment holder is made from a double layer of densely pressed cardboard
- Cover is single layer cardboard
- Front protection flaps are attached with durable tape



Item No.	Color	Capacity	L x W x H in	L x W x H cm
HS9902	Black	20-Place	8.2 x 13.4 x 0.4	20.8 x 34.1 x 0.9
HS9913	Blue			
HS9914	Green			
120196	Red			
120197	Yellow	30-Place	8.2 x 12.9 x 0.4	28.1 x 32.8 x 1
HS9930	Black			

Slide Tray Clasps

Effortlessly locate, retrieve and keep slide folders organized

- Assorted color clasps are ideal for color coding identification and to keep slide folder flaps securely closed
- Clasps can be labeled with a permanent marker or a label can be attached

Slide Tray Clasps

Pack of 100

Item No.	Color	L x W x H in	L x W x H cm
120653	Assorted*	2.2 x 1.7 x 0.6	5.6 x 4.4 x 1.6
120654	Clear		

*Assorted colors (Clear, Grey, Blue, Green, Red)



Have a Question?

Our experienced and responsive CUSTOMER SERVICE and TECHNICAL SUPPORT teams are ready and available to help.



PHONE:
+1.847.816.5070
+44 (0) 845.450.7316



EMAIL:
info@heatsci.com

Slide Mailers

Safe and secure transport of standard microscope slides
(76 x 26 mm and 3 x 1 in)

Cardboard Slide Mailers

Item No.	Color	Capacity	UOM	L x W x H in	L x W x H mm
HS9903	Natural	1-Place	Pack of 25	3.7 x 1.7 x 0.2	95 x 42 x 5
HS9904	Natural	2-Place	Pack of 36	4.1 x 3.1 x 0.2	103 x 79 x 5
120557	White	4-Place	Pack of 20	4.7 x 5.7 x 0.3	120 x 140 x 6

- Heavy cardboard mailers for safe transport
- Cover that opens for viewing and thumb grooves for easy slide removal
- Securely holds slides with or without cover glass
- Flat surface area is ideal for labeling or writing necessary information



Plastic Slide Mailers

Item No.	Color	Capacity	UOM	L x W x H in	L x W x H mm
120556	Natural	1-Place	Pack of 50	3.3 x 1.6 x 0.2	83 x 41 x 5
HS15983B	Natural				
HS15983A	Assorted*	2-Place	Pack of 25	3.7 x 2.7 x 0.2	93 x 68 x 5

*Assorted colors (Blue, Purple, Red, Orange, Yellow)

- Polypropylene protectors ideal for secure shipping and storing of slides
- Large, easy-to-secure snaps, provide safe storage for valuable slides
- Variety of colors brings some fun to the lab and enables color coding



Plastic 5-Place Slide Mailers

Item No.	Color	Capacity	UOM	L x W x H in	L x W x H mm
HS15982	Natural	5-Place	Pack of 25	3.2 x 0.7 x 1.2	81 x 19 x 31
HS15983G	Assorted*	5-Place	Pack of 25	1.7 x 0.9 x 3.5	43 x 24 x 88
HS15986	Natural				

*Assorted colors (Blue, Purple, Red, Orange, Yellow)


- Polypropylene containers ideal for secure shipping, storing or staining slides
- Polypropylene containers hold five standard slides (76 x 26 mm and 3 x 1 in)
- Inside wall features slotted channels to keep slides separated and to make removal easier
- Lid closes securely to keep slides safe



Slide Holders

Hold slides during evaluation, cleaning and drying

- Made from chemical-resistant polyethylene making the unit resist stains and easy to clean
- Solid, strong construction with rubber feet for bench-top stability
- Angled slide holder slots allow easy access to individual slides and viewing of slide labels



Item No.	Color	Capacity	L x W x H in	L x W x H cm
HS15981A	Natural <input type="checkbox"/>	40-Place	11.7 x 4.1 x 0.9	29.8 x 10.5 x 2.4
HS15981B	Natural <input type="checkbox"/>	78-Place	11.8 x 7 x 0.8	29.9 x 17.8 x 2.1
HS15981C	Natural <input type="checkbox"/>	15-Place	6 x 4 x 0.8	15.3 x 10.2 x 2.1

Slide Staining Tray Set

Reduce possible cross-contamination and messy work areas

This rugged staining tray is an economical alternative to metal versions

- Made from a polypropylene and polyethylene blend for stability
- Stain, rinse, and dry slides on a single working tray
- A clean and fresh working surface area every time
- Up to eight slides fit comfortably onto the base and the deep well holds up to 38 mL
- Compact size and recessed handles for easy transport from work area to sink
- Pack includes 4 trays and 1 lid



Dark Lid
protects slides for light-sensitive applications

Deep Well
holds up to 38 mL of liquid

Item No.	Color	Capacity	L x W x H in	L x W x H cm
HS15951A	Black <input checked="" type="checkbox"/>	8-Place	12 x 5 x 1.2	30.5 x 12.7 x 3

2 Convenient Pour Spouts
for quick disposal of unwanted liquid waste

Safety Products & Dispensers

Protect and comply with our range of safety products and PPE dispensing solutions



Wash Bottles

Biohazard Containment

PPE Dispensers & Holders

Wash Bottles Color-Coded






Color-coded for easy identification

Pack of 6

The closure and stem are molded in one piece to make them leak-proof

The angled neck and tip provide a uniform stream of liquid and the tip can be cut back if required to increase the flow


- Wide bottle mouth makes filling easier
- Bottle is made from LDPE, with a PP closure
- Color-coding for user, application or location identification

Item No.	Cap/Lid Color	Capacity
120245	White 	500 mL
120246	Red 	
120247	Yellow 	
120248	Green 	
120249	Blue 	



Create a safer working surface with our reusable Silicone Lab Mats

- Anti-skid to keep glassware in place
- Protect counter tops from stain, spills and wear
- Available in multiple colors and sizes



Wash Bottles Safety Labeled and Self Venting

Clear and informative product safety information

The information provided on the front of each bottle is specific to the product contained and includes the following — the Chemical Abstract Service (CAS) No., National Fire Protection Association (NFPA) codes and symbols, Protective Equipment Symbol (PES), Target Organ Information and Hazard Codes

There is also the option of using other liquids within the bottles and creating the appropriate safety information on the outside of the bottle by using the self-labeling versions

The closure and stem are molded in one piece to make them leak-proof, the angled neck and tip provide a uniform stream of liquid and the tip can be cut back to increase the flow

- Wide bottle mouth makes filling easier
- Sodium Hypochlorite bottle made from natural LDPE for light protection and has a HDPE closure
- All other bottles are made from LDPE, with a PP closure
- Non-metallic vent that prevents the risk of solvent leakage
- Six named chemicals, a self-labeling** and assorted pack of bottles to choose from

Pack of 6

Item No.	Cap/Lid Color	Description	Capacity
120250	Red	Acetone	500 mL
120251	Clear	Distilled Water	
120252	Yellow	Isopropanol	
120253	White	Ethanol	
120254	Green	Methanol	
120255	Blue	Sodium Hypochlorite	
120256	White	Self-labeling**	
120257	Assorted*	Assorted*	


* Assorted Pack includes one of each: acetone, distilled water, isopropanol, ethanol, sodium hypochlorite & a self-labeling bottle

**Self-labeling bottles must only be used with chemicals that are safe and compatible with LDPE. Labeling stickers included



smartSlicer™**Laboratory Safety Razor Blade****Say no to dangerous sharps**

- At just 0.27 mm the US Patent Pending product is as thin as a single edge industrial razor blade, but much safer
- Cleanly slice agarose gels right on the light box without scratching the surface
- Ideal for cutting soft tissue and plant matter
- Adjust the size required by simply cutting with scissors
- Made from food grade PET



Pack of 100

Item No.	L x W x H in	L x W x H mm
120259	0.8 x 1.5 x 0.01	19 x 38 x 0.3
120260	1 x 2 x 0.01	25 x 50 x 0.3

smartSlicer™ is a trademark of LevGo Inc.

Safely seal PP and PE bags, sheets or flat tubing with our Heat Sealer

- Heating element measures 31 cm long and is only hot when sealing
- Remove excess material safe with integrated sliding blade

*Heat Sealer with Cutter, Page 120***Lead Rings****Stabilize vessels and labware on the bench or in water baths**

- Vinyl-coated lead rings will not mark or scratch benchtops and labware
- C-shape format allows some manual size variation for a customized fit

Item No.	Shape	Color	Flask Sizes mL	ID in	ID cm	OD in	OD cm	Weight g
HS8881A	C	Blue	125–500	1.9	4.2	3	7.5	215
HS8882A	O			1.9	4.8	3	7.7	210
HS8881B	C	Blue	250–1000	2	5.1	3.5	9	465
HS8882B	O			2	5.1	3.5	9	490
HS8881C	C	Blue	500–2000	2.1	5.4	3.9	9.8	595
HS8882C	O			2.2	5.7	3.9	9.8	610
HS8881D	C	Blue	1000–4000	2.6	6.6	4.4	11.3	890
HS8882D	O			2.8	7	4.5	11.5	935

**Weights may vary by +/- 3%*

Biohazard Disposal Bags

For personal or general lab use, sizes to suit all needs

- Polypropylene
- Includes an indicator patch that darkens when autoclaved
- A range of biohazard collection bags printed with biohazard warning symbol and precautions in English, Spanish and French



Biohazard Disposal Bags

Item No.	Color	Thickness mil	Capacity	Packaging	L x W in	L x W cm
HS10320	Red ■	1.6	2 L	500	8 x 12	20.3 x 30.5
HS10321			9 L	200	14 x 19	35.6 x 48.3
HS10322			19 L	200	19 x 23	48.3 x 58.4
HS10323			40 L	200	25 x 33	63.5 x 83.8

HDPE Coated Steel Wire Bag Holder

Item No.	L x W x H in	L x W x H cm
HS10044H	3.6 x 5 x 8.4	9.1 x 12.8 x 21.3

**Holds Biohazard Disposal Bag HS10320*



Epoxy Coated
HDPE Wire Holder



Autoclavable Disposal Bags

Clear bags imprinted with blue biohazard labeling designed for safe disposal

- Constructed of 2 mil thick polypropylene, allowing for repeated autoclaving
- Front write-on panel
- 4 inner packs of 50 bags
- Strong side seams for securely holding items



Case of 200

Item No.	Color	Capacity	L x W in	L x W cm
HS1002A	Clear □	14 L	12.2 x 26	31 x 66
HS1002B		26 L	16.3 x 23.6	41.5 x 60
HS1002C		69 L	24 x 31.9	62 x 81

Sharps Containers

Easy to use and transport

- Clearly labeled with the Biohazard warning on one side
- Provides safe disposal of syringes, needles and other small sharps
- Stackable to reduce storage space before use
- 1 L Capacity (quart) with a clear lid to view fill level
- Slide access lid provides temporary and final closures for safety and security
- Built-in needle removal port, prevents direct contact with sharps
- Compatible with the Droplet blood collection tray, Totemax® and blood collection tray, Pages 34 and 35

Item No.	Color	L x W x H in	L x W x H cm
120177	Red ■	4.5 x 4.5 x 7.5	11.5 x 11.5 x 19.2
120178	Yellow ■		



PPE Dispensers & Holders



Soft Cover Dispenser

Store and dispense with ease

Translucent acrylic dispensers provide easy access to safety accessories

- Hinged lid protects contents from dust and debris
- Set on counter or wall mount (screws included)
- Great for holding other personal items



Item No.	Color		L x W x D in	L x W x D cm
HS1040C	Clear	□	8 x 8 x 16	20.3 x 20.3 x 40.6

Earplug Dispenser

Translucent acrylic dispensers provide easy access to safety accessories

- Holds approximately 200 pairs of foam earplugs
- Set on counter or wall mount (screws included)



Item No.	Color		L x W x D in	L x W x D cm
HS1040B	Neon Green	■	8 x 8 x 16	20.3 x 20.3 x 40.6

3-in-1 Dispenser

Three separate lift-up lids make loading and dispensing supplies easy

Three-compartment dispenser can be used to organize safety glasses, earplugs, soft covers (including bouffant caps, beard covers, shoe covers) or any small items that need to be kept within reach

The unit includes screws for wall mounting or can be set on a counter

- Middle compartment has reach-in hole for dispensing
- Durable acrylic is see-through so product levels are always visible

3
Compartment



Item No.	Color		L x W x D in	L x W x D cm
HS1042	Grey	■	16 x 16 x 8	40.6 x 40.6 x 20.3

Safety Glasses Dispensers

Ergonomically angled for easy dispensing

Translucent acrylic dispensers provide easy access to safety accessories

- Hinged lid protects glasses from dust and debris
- Set on counter or wall mount (screws included)
- Large dispenser holds approx. 20 glasses, small approx. 10 glasses
- Great for holding other personal items

Item No.	Color		L x W x D in	L x W x D cm
HS1041	Neon Blue	■	8 x 6 x 6	20.3 x 15.2 x 15.2
HS1040A	Neon Blue	■	8 x 8 x 16	20.3 x 20.3 x 40.6



PPE Dispensers & Holders

Workstation Storage Bin

The essential combination accessory holder

- Organize your most commonly used laboratory disposables in one central location
- Acrylic unit has room to store gloves, Kimwipes® and Parafilm® M
- Use bins to store tubes, tips, etc.

Item No.	Color		L x W x H in	L x W x H cm
HS234523	Clear	<input type="checkbox"/>	10.4 x 5.4 x 17.3	26.5 x 13.7 x 44

Parafilm® is a registered trademark of Bemis Inc.

Kimwipes® is a registered trademark of Kimberly-Clark



Lab Wipe Holders

Holds Kimwipes® and other major brands

Clear acrylic holders fit boxes of small or large lab wipes

- The smaller version attaches with double-sided adhesive tape and the larger version with screws
- Mount on a shelf or wall
- Mounting materials included

Item No.	Description	Color		L x W x H in (cm)
HS234521	Small Holder	Clear	<input type="checkbox"/>	4.6 x 4.8 x 3.3 (11.8 x 12.2 x 8.4)
HS234522	Large Holder	Clear	<input type="checkbox"/>	12.6 x 5 x 3.8 (32 x 12.7 x 9.7)

Kimwipes® is a registered trademark of Kimberly-Clark


Benchtop Dispensing Bins

Work great with gloved hands, bright colors

Ergonomically designed with large, angled openings and oversized lip

- Lid remains open when tilted back
- Nonskid rubber feet and pre-drilled holes for wall mounting
- Green bin includes a removable divider
- Mounting screws included

Item No.	Color		Compartments	L x W x H in	L x W x H cm
HS23411	Neon Green	<input checked="" type="checkbox"/>	2	7.1 x 6.1 x 6.5	18 x 15.5 x 16.5
HS23412	Neon Blue	<input checked="" type="checkbox"/>	1	6.1 x 6.1 x 6.7	15.5 x 15.5 x 17



Clear Acrylic Glove Box Holders

Keep gloves organized & ready for immediate use

- Free up bench-top space by mounting holders on wall for easy dispensing of gloves
- Universal fit for most major brands of glove boxes
- Mounting screws included
- Includes a polyester foam insert to securely hold a range of glove box sizes—snug fit ensures easy removal of gloves from the box

Item No.	Capacity	L x W x H in	L x W x H cm
HS23471A	1	10.3 x 5.9 x 4.5	26.3 x 14.6 x 11.4
HS23471B	2	11.4 x 10.3 x 4.5	29.2 x 26 x 11.4
HS23471C	3	11.6 x 15.2 x 4.3	29.4 x 38.8 x 11.6



HDPE Coated Wire Glove Box Holders

Keep gloves organized & ready for immediate use

- HDPE-coated steel wire racks are an economical solution
- Free up bench-top space by mounting holders on wall for easy dispensing of gloves
- Universal fit for most major brands of glove boxes
- Mounting screws included

Item No.	Capacity	L x W x H in	L x W x H cm
HS23457	1	5.5 x 8.2 x 4.3	14 x 20.9 x 10.8
HS23458	3	9.8 x 17.9 x 4.3	25 x 45.4 x 10.8

Modular Glove Box Holders

Available infused with Microban® for anti-microbial protection

The anti-microbial version (blue) is ideal for use in clean rooms, hospitals, and microbiological work areas. The Microban additive prevents growth of molds and bacteria

- High-impact, durable resin is resistant to harsh cleaning solutions
- Universal fit for most major brands of glove boxes
- Modular design allows single boxes or horizontal stacking for multiple box holders
- Mounting screws included

US Patent No. 6669157 Pack of 3

Item No.	Description	Color	L x W x H in	L x W x H cm
120258	Standard	White <input type="checkbox"/>	10 x 6.1 x 3.8	25.4 x 15.6 x 9.7
HS23491A	Anti-microbial	Blue <input checked="" type="checkbox"/>		

Microban® antimicrobial protection has been built-in to our anti-microbial glove boxes during manufacturing. When microbes, such as bacteria, mold, and mildew that can cause stains, odors and product deterioration come in contact with the glove box surface, Microban protection penetrates the cell wall of the microbe and disrupts key cell functions so that the microbe cannot function, grow, or reproduce. Microban protection is an intrinsic part inside and at the surface and will not wash off or wear away for the useful lifetime of the glove box.





Flexible Glove Box Holder

Glove box holder ideal for gloves, face masks, shoe covers and more

The adjustable holding points provide true flexibility to hold up to four different boxed protective products and other lab essentials (size dependent)

- Ideal for holding different boxed sizes of gloves, mop caps, ear plugs, face masks, shoe covers
- Adjustable universal fit
- Easy to clean, with a low profile
- Wall- or door-mountable
- Mounting hardware included

Adaptable

Store different size boxes
all in one holder

Item No.	Color	L x W x D in	L x W x D cm
120107	White <input type="checkbox"/>	21.9 x 5.2 x 2.7	55.8 x 13.2 x 6.8

Delivery

We stock our warehouse to ensure product availability and on-time delivery.

By offering low minimum order quantities, we help dealers efficiently maintain lower stocking levels.



Sealing & Labeling

Sample protection and identification for a wide range of tubes and vials



Heat Sealer

Parafilm & Dispensers

Films, Tapes & Foils

Labels

Heat Sealer with Cutter

Quickly and efficiently seal polyethylene or polypropylene bags, sheets or flat tubing from 0.01 to 0.8 mm thick

Sliding blade removes excess material 19 mm from the seal

- Can seal a bag up to 12 in (30.5 cm) long
- Heating element measures 31 cm long and is only hot when sealing
- Includes one replacement part kit
- Additional replacement part kits can be ordered separately



Replacement Part Kit:

- 2 extra blades,
- 2 extra heating elements,
- 2 extra PTFE cloths



Timing knob ensures safe sealing



Item No.	Power Supply		L x W x H in	L x W x H cm
HS28270	220-230 VAC, UK	⏏	17.6 x 3.4 x 6.5	44.6 x 8.6 x 16.6
HS28271	220-230 VAC, Cont. Europe/Israel	⏏		
HS28272	110-120 VAC, US/Taiwan/Japan	⏏		
HS28280	Replacement Part Kit			

Parafilm® M Sealing Film

Protect the integrity of your laboratory and research applications

Parafilm® M is a semi-transparent, flexible, thermoplastic film with two-way stretch

- Protect the integrity of research applications with a tight moisture-barrier that holds a strong seal and doesn't release from glassware, despite off-gassing
- Reduce the risk of tampered data with a covering that is odorless, colorless and semi-transparent
- Prevents moisture and volume loss due to material composition and tight seal
- Flexible and self-sealing material allows it to cling around irregular shapes and surfaces
- Stretches more than 200% of original length
- Resistant up to 48 hours against many polar substances such as saline solutions, inorganic acids and alkaline solutions
- Two color options available in 2" rolls (natural and purple) for color coding purposes



Item No.	Color		Film Width		Roll Length	
			in	cm	ft	m
HS234526A	Natural	□	2	5	250	76.2
120762	Purple	■	2	5	250	76.2
HS234526B	Natural	□	4	10	125	38.1
HS234526C	Natural	□	4	10	250	76.2
120265	Natural	□	20	50.8	50	15.2

Parafilm® M is a registered trademark of Amcor Limited

Box Top Dispenser for Parafilm® M Sealing Film

Easy application

Simply place the dispenser on top of the Parafilm® M box and pull the sealing film through, there is no need to remove rolls from the packaging

- Accommodates two 50 mm or one 100 mm Parafilm® M rolls
- Built in safety blade slides and cuts with no danger to fingertips and eliminates the need for scissors
- Supplied with two safety blades, additional replacement blades are available for purchase
- Made with ABS plastic making the unit easy to clean



Item No.		L x W x H in	L x W x H cm
120667	Blue ■		
120668	Purple ■	4.9 x 6.8 x 3.9	12.3 x 17.2 x 9.9
120758	Natural □		

Item No.	
120722	Replacement Blades for Box Top Dispenser, Pack of 3 <i>Parafilm® M is a registered trademark of Amcor Limited</i>



ABS Dispenser for Parafilm® M Sealing Film

Store, dispense, and cut Parafilm® M, tape, labels, and Tough-Tags®

- ABS plastic dispenser designed with integrated safety razor that slides smoothly creating a clean cut
- Supplied with two safety blades, additional replacement blades are available for purchase
- Dispenser holds one or two 50 mm (2-inch) rolls or one 100 mm (4-inch) roll

Item No.		L x W x H in	L x W x H cm
HS234525B	Blue ■		
HS234525C	Green ■	4.7 x 6.1 x 6.7	12 x 15.6 x 17.1
HS234525D	Red ■		

Item No.	
HS234525R	Replacement Blades for ABS Dispenser, Pack of 5 <i>Parafilm® M is a registered trademark of Amcor Limited</i> <i>Tough-Tags® is a registered trademark of Diversified Biotech</i>

Acrylic Dispenser for Parafilm® M Sealing Film

Store, dispense, and cut Parafilm® M, tape, labels, and Tough-Tags®

- Strong, crystal-clear acrylic dispenser makes your film supply visible and features a serrated cutting blade
- Rubber feet secures in place while cutting
- Dispenser holds one or two 50 mm (2-inch) rolls or one 100 mm (4-inch) roll

Item No.	L x W x H in	L x W x H cm
HS234524	4.7 x 6.6 x 6.9	11.9 x 16.8 x 17.5

Parafilm® M is a registered trademark of Amcor Limited
Tough-Tags® is a registered trademark of Diversified Biotech



PetriSEAL™ and ContainerSEAL™

Sealing tape for petri dishes, microplates, lids and caps

- Tape stretches, conforms and seals to irregular perimeters
- Seal remains flexible even with intermittent exposure to 100°C
- Waterproof and chemically resistant to most solvent and caustic agents

PetriSEAL™

Item No.	Color	Film Width		Roll Length	
		in	cm	ft	m
HEAPSEAL-108B	Blue				
HEAPSEAL-108C	Clear				
HEAPSEAL-108R	Red	0.5	1.3	108	33
HEAPSEAL-108Y	Yellow				
HEAPSEAL-108W	White				



ContainerSEAL™

Item No.	Color/Size	Film Width		Roll Length	
		in	cm	ft	m
HEACSEAL-75R	Red				
HEACSEAL-75W	White	0.8	1.9	108	33
HEACSEAL-75Y	Yellow				
HEACSEAL-XLR	Red XL				
HEACSEAL-XLW	White XL	1	2.5	108	33
HEACSEAL-XLY	Yellow XL				

PetriSeal™ and ContainerSeal™ are trademarks of Diversified Biotech, Inc



Foil Rolls

Designed for easy use and minimal waste

- Two sizes available for covering small surface such as beakers, flasks, tubes and Petri dishes
- Withstands autoclaving, incubating and freezing temperatures

Item No.	Thickness	Roll Width		Roll Length	
	mm	in	cm	ft	m
HS23534A	0.015	4	10.2	500	152.4
HS23536A	0.015	6	15.3	500	152.4




Cryo-Babies® and Cryo-Tags® Low Temperature Labels

For secure uncompromised labeling

- Chemically inert
- Adheres to most plastics (100% virgin PP tubes recommended), glass and metals without cracking, peeling or degrading
- Flexible polyolefin material and acrylic adhesive expand and contract with temperature fluctuations



Cryo-Babies® on a Roll

Ideal Use	0.5 mL Tubes	1.5 mL Tubes	General Purpose
			
Label Dimensions	0.94 x 0.50 in (24 x 13 mm)	1.28 x 0.50 in (33 x 13 mm)	1.50 x 0.75 in (38 x 19 mm)
UOM	1,000/roll	1,000/roll	1,000/roll
Item No.			
White	HEACRY-0500	HEACRY-1000	HEACRY-2000
Blue	-	HEACRY-1000-B	HEACRY-2000-B
Green	-	HEACRY-1000-G	HEACRY-2000-G
Orange	-	HEACRY-1000-O	HEACRY-2000-O
Red	-	HEACRY-1000-R	HEACRY-2000-R
Yellow	-	HEACRY-1000-Y	HEACRY-2000-Y

Cryo-Babies® and Cryo-Tags® are registered trademarks of Diversified Biotech, Inc

For write-on labels, see our selection of Sharpie® Pens



Cryo-Babies® and Cryo-Tags® Low Temperature Labels

For secure uncompromised labeling

- Chemically inert
- Adheres to most plastics (100% virgin PP tubes recommended), glass and metals without cracking, peeling or degrading
- Flexible polyolefin material and acrylic adhesive expand and contract with temperature fluctuations



Laser Cryo-Tags® and Cryo-Babies® on Sheets formatted for laser printers

Ideal Use	0.5 mL Tubes	1.5 mL Tubes	1.5 mL Tubes	Cryogenic Vials	Bottles & Beakers	Microplates
Label Dimensions	0.94 x 0.50 in (24 x 13 mm)	1.28 x 0.50 in (33 x 13 mm)	1.28 x 0.50 in (33 x 13 mm)	1.69 x 0.75 in (43 x 19 mm)	3.0 x 2.0 in (76 x 51 mm)	1.50 x 0.25 in (38 x 6 mm)
Labels/Sheet	119	51	85	52	12	156
Labels/pk	2,380	1,020	1,700	1,040	240	3,120
UOM	20/pk	20/pk	20/pk	20/pk	20/pk	20/pk
Item No.						
White	HEALCRY-2380	HEAMINI-1700	HEALCRY-1700	HEALCRY-1100	HEALCRY-3020	HEASIDE-1000
Blue	HEALCRY-2380-B	-	HEALCRY-1700-B	HEALCRY-1100-B	-	HEASIDE-1100
Green	HEALCRY-2380-G	-	HEALCRY-1700-G	HEALCRY-1100-G	-	HEASIDE-1200
Orange	HEALCRY-2380-O	-	HEALCRY-1700-O	HEALCRY-1100-O	-	HEASIDE-1300
Red	HEALCRY-1380-R	-	HEALCRY-1700-R	HEALCRY-1100-R	-	HEASIDE-1400
Yellow	HEALCRY-1380-Y	-	HEALCRY-1700-Y	HEALCRY-1100-Y	HEALCRY-3020-Y	HEASIDE-1500
Rainbow*	HEARNBW-2300	-	HEARNBW-2200	HEARNBW-2000	-	HEARNBW-2400

Laser Cryo-Tags® and Cryo-Babies® on Sheets formatted for laser printers - continued

Ideal Use	Microscope Slides	Racks & Boxes	Racks & Boxes	Racks & Boxes	General Purpose	General Purpose
Label Dimensions	1.0 x 1.0 in (25 x 25 mm)	2.50 x 0.50 in (64 x 13 mm)	2.625 x 1.0 in (67 x 25 mm)	1.125 x 1.75 in (29 x 45 mm)	8.5 x 11.0 in (216 x 279 mm)	1.50 x 0.75 in (38 x 19 mm)
Labels/Sheet	80	60	30	42	1	60
Labels/pk	1,600	1,200	600	840	20	1,200
UOM	20/pk	20/pk	20/pk	20/pk	20/pk	20/pk
Item No.						
White	HEALCRY-1400	HEALCRY-2500	HEALCRY-1258	HEALCRY-1300	HEALCRY-1000	HEALCRY-1200
Blue	-	-	HEALCRY-1258-B	-	-	HEALCRY-1200-B
Green	-	-	HEALCRY-1258-G	-	-	HEALCRY-1200-G
Orange	-	-	HEALCRY-1258-O	-	-	HEALCRY-1200-O
Red	-	-	HEALCRY-1258-R	-	-	HEALCRY-1200-R
Yellow	-	-	HEALCRY-1258-Y	HEALCRY-1300-Y	-	HEALCRY-1200-Y
Rainbow*	-	-	HEARNBW-2500	-	-	HEARNBW-2100

*Rainbow 4 sheets of each Blue, Green, Orange, Red, Yellow
Cryo-Babies® and Cryo-Tags® are registered trademarks of Diversified Biotech, Inc

Tough-Spots® Tube Labels

Pre-cut peel off labels for tube tops

- Adheres to most plastics, glass, and metals without cracking, peeling or degrading
- Flexible polyolefin material and acrylic adhesive



Tough-Spots® on a Roll






Ideal Use	0.2 mL Tubes	0.5 mL Tubes	1.5/2.0 mL Tubes	1.5/2.0 mL Tubes	15-50 mL Tubes	50 mL Tubes
Label Diameter	0.1875 in/4.8 mm	0.375 in/9.5 mm	0.4375 in/11 mm	0.50 in/13 mm	0.75 in/19 mm	1.0 in/25 mm
UOM	1,500/roll	1,000/roll	1,000/roll	1,000/roll	500/roll	500/roll
Item No.						
White	HEATS-TNY	HEAT-SPOTS	HEATSPOTS-44	HEAT-SPOTS-50	HEATSPOTS-75	HEAT-SPOTS-100
Blue	-	HEAT-SPOTS-B	-	HEAT-SPOTS-50-B	-	-
Coffee	-	HEAT-SPOTS-C	-	HEAT-SPOTS-50-C	-	-
Gold	-	HEAT-SPOTS-D	-	HEAT-SPOTS-50-D	-	-
Green	-	HEAT-SPOTS-G	-	HEAT-SPOTS-50-G	-	-
Lavender	-	HEAT-SPOTS-L	-	HEAT-SPOTS-50-L	-	-
Maroon	-	HEAT-SPOTS-M	-	HEAT-SPOTS-50-M	-	-
Orange	-	HEAT-SPOTS-O	-	HEAT-SPOTS-50-O	-	-
Pink	-	HEAT-SPOTS-P	-	HEAT-SPOTS-50-P	-	-
Red	-	HEAT-SPOTS-R	-	HEAT-SPOTS-50-R	-	-
Silver	-	HEAT-SPOTS-S	-	HEAT-SPOTS-50-S	-	-
Yellow	-	HEAT-SPOTS-Y	-	HEAT-SPOTS-50-Y	-	-



- Label sheets are top coated and heat-resistant to prevent jams in laser printers
- Printer formatting instruction included



Laser Tough-Spots® on Sheets formatted for laser printers

Ideal Use	0.2 mL Tubes	0.5 mL Tubes	1.5/2.0 mL Tubes	1.5/2.0 mL Tubes	15-50 mL Tubes
					
Label Diameter	0.375 in/9.5 mm	0.4375 in/11 mm	0.50 in/13 mm	0.75 in/19 mm	1.0 in/25 mm
Labels/sheet	192	234	192	108	63
Labels/pk	3,840	4,680	3,840	2,160	1,260
UOM	20/pk	20/pk	20/pk	20/pk	20/pk
Item No.					
White	HEASPOT-1000	HEASPOT-5000	HEASPOT-2000	HEASPOT-3000	HEASPOT-4000
Blue	HEASPOT-1100	-	HEASPOT-2100	-	-
Green	HEASPOT-1200	-	HEASPOT-2200	-	-
Orange	HEASPOT-1300	-	HEASPOT-2300	-	-
Red	HEASPOT-1400	-	HEASPOT-2400	-	-
Yellow	HEASPOT-1500	-	HEASPOT-2500	-	-
Rainbow*	HEARNBW-1000	-	HEARNBW-1100	-	-

*Rainbow 4 sheets of each Blue, Green, Orange, Red, Yellow
Tough-Spots® is a registered trademarks of Diversified Biotech, Inc

Tough-Tags® Tube Labels

Pre-cut labels for identification

- Adheres to most plastics, glass, and metals without cracking, peeling or degrading
- Flexible polyolefin material and acrylic adhesive



Tough-Tags® on a Roll

Ideal Use	0.2 mL Tubes	0.5 mL Tubes	1.5/2.0 mL Tubes	General Purpose	Microplates
Label Diameter	0.81 x 0.28 in/21 x 7 mm	0.94 x 0.50 in/24 x 13 mm	1.28 x 0.50 in/33 x 13 mm	1.50 x 0.75 in/38 x 19 mm	1.50 x 0.28 in/38 x 6 mm
UOM	1,500/roll	1,000/roll	1,000/roll	1,000/roll	1,000/roll
Item No.					
White	HEATT-TNY	HEATTSW-1000	HEATTLW-1000	HEATTGP-1000	HEATT-SWALL
Blue	-	HEATTSB-1000	HEATTLB-1000	-	-
Green	-	HEATTSG-1000	HEATTLG-1000	-	-
Orange	-	HEATTSO-1000	HEATTL0-1000	-	-
Red	-	HEATTSR-1000	HEATTLR-1000	-	-
Yellow	-	HEATTSY-1000	HEATTLY-1000	-	-
Clear	-	HEATTSC-1000	HEATTLc-1000	-	-

- Label sheets are top coated and heat-resistant to prevent jams in laser printers
- Printer formatting instruction included



Tough-Tags®
is a registered trademark
of Diversified Biotech, Inc

Laser Tough-Tags® on Sheets formatted for laser printers

Ideal Use	0.5 mL Tubes	1.5/2.0 mL Tubes	General Purpose	Microscope Slides
Label Dimensions	0.94 x 0.50 in/24 x 13 mm	1.28 x 0.50 in/33 x 13 mm	1.50 x 0.75 in/38 x 19 mm	0.94 x 0.94 in/22 x 22 mm
Labels/Sheet	119	85	60	96
Labels/pk	2,975	2,125	1,500	2,400
UOM	25/pk	25/pk	25/pk	25/pk
Item No.				
White	HEATTSW-2240	HEATTLW-2016	HEATTLW-1050	HEAMISL-1000
Blue	HEATTSW-2240-B	HEATTLW-2016-B	HEATTLW-1050-B	HEAMISL-1000-B
Green	HEATTSW-2240-G	HEATTLW-2016-G	HEATTLW-1050-G	HEAMISL-1000-G
Orange	HEATTSW-2240-O	HEATTLW-2016-O	HEATTLW-1050-O	HEAMISL-1000-O
Red	HEATTSW-2240-R	HEATTLW-2016-R	HEATTLW-1050-R	HEAMISL-1000-R
Yellow	HEATTSW-2240-Y	HEATTLW-2016-Y	HEATTLW-1050-Y	HEAMISL-1000-Y
Rainbow*	HEARNBW-3100	HEARNBW-3000	HEARNBW-3200	HEARNBW-3300

*Rainbow 4 sheets of each Blue, Green, Orange, Red, Yellow

Label Station

Organize your labels for easy accessibility


Portable label station holds up to six standard size boxes of Tough-Tags®, Cryo-Babies® and Tough-Spots®

Just drop the box in the label station and feed the paper through the slot

- Acrylic station includes a writing surface, a cutting edge for easy dispensing and non-skid rubber feet
- Counter-sunk holes provided for wall mounting

Item No.	Color	L x W x H in	L x W x H cm
HS23460	Blue	11 x 4.7 x 3.5	28 x 12 x 8.9

Tough-Tags®, Cryo-Babies® and Tough-Spots®
are registered trademarks of Diversified Biotech, Inc.



Labels shown not included,
order separately

Charts & References



Plastic Material References

Product Number Index

Product Alpha Index

Plastic Material References

G Good Resistance: continuous exposure to the substance does not cause damage within 30 days

F Fair Resistance: continuous exposure to the substance causes minor damage within 7-30 days

P Poor Resistance: not suitable for continuous exposure to the substance. Immediate damage may occur

Chemical Resistance of Plastics

	LDPE		HDPE		PP		PMP (TPX)		PVC		PC		ABS/PS		SAN		PMMA		PTFE		POM	
TEMPERATURE °C	20	50	20	50	20	50	20	50	20	50	20	50	20	50	20	50	20	50	20	50	20	50
Acetaldehyde	G	P	G	F	F	P	F	P	P	P	F	P	P	P	P	P	P	P	G	G	G	G
Acetic Acid (Glacial)	F	P	G	G	G	F	F	F	G	P	P	P	P	P	P	P	P	P	G	G	P	P
Acetic Anhydride	P	P	F	F	F	F	G	F	P	P	P	P	P	P	P	P	P	P	G	G	P	P
Acetone	F	P	P	P	F	P	G	G	P	P	P	P	P	P	P	P	P	P	G	G	G	G
Ammonium Chloride (10%)	G	G	G	G	G	G	G	G	G	F	F	F	G	G	G	G	F	F	G	G		
Ammonium Hydroxide (30%)	G	G	G	G	G	G	G	G	G	P	P	F	F	F	F	F	F	F	G	G	G	G
Amyl Acetate	F	P	G	F	F	P	G	F	P	P	P	P	P	P	P	P	G	G	G	G	G	G
Aniline (Phenylamine)	G	F	G	G	G	G	G	F	P	P	P	P	P	P	P	P	P	P	G	G	F	F
Aqua Regia	P	P	P	P	P	P	P	P	P	P	P	P	P	P	F	P	P	P	G	G	G	G
Benzaldehyde	G	F	F	P	P	F	G	F	P	P	P	P	P	P	P	P	P	P	G	G	G	G
Benzene	P	P	P	P	P	P	P	P	P	P	P	P	P	P	P	P	P	P	G	G	G	F
Benzoic Acid	G	G	G	G	G	F	G	G	G	F	G	F	F	F	G	G	G	G	G	G		
Boric Acid (10%)	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G
Butyl Acetate	F	F	F	F	F	F	F	F	G	P	P	P	P	P	P	P	P	P	G	G	G	F
Butyric Acid (Butanoic acid)	P	P	F	P	P	P	P	P	G	P	P	P	P	P	P	P	P	P	G	G		
Calcium Hydroxide (Saturated)	G	G	G	G	G	G	G	G	G	G	P	P	F	F	G	G	F	P	G	G	G	G
Carbon Disulphide	P	P	P	P	P	P	P	P	P	P	P	P	P	P	P	P	P	P	G	G	G	G
Carbon Tetrachloride	F	P	F	P	F	P	P	P	G	P	P	P	P	P	F	P	P	P	G	G	F	F
Chloroform	F	P	F	P	P	P	P	P	G	P	P	P	P	P	P	P	P	P	G	G	P	P
Citric Acid (1M)	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G		
Cresol	P	P	F	P	F	F	P	P	P	P	P	F	G	P	P	P	P	P	G	G		
Cyclohexane	F	P	F	P	F	P	P	P	P	P	G	F	P	P	F	P	F		G	G	G	G
Dibutyl Phthalate	F	P	F	P	F	P	F	F	P	P	F	P	P	P	P	P	P	P	G	G	G	G
Dichlorobenzene	F	P	F	P	F	P	P	P	P	P	P	P	P	P	P	P	P	P	G	G	P	P
Diethyl Ether	P	P	F	P	F	P	P	P	G	P	P	P	P	P	P	P	P	P	G	G	G	G
Diethylene Glycol	G	G	G	G	G	G	G	G	F	P	F	F	G	G	G	G	P	P	G	G	G	G
Dimethyl Formamide (DMF)	G	G	G	G	G	G	G	G	P	P	P	P	P	P	P	P	P	P	G	G	G	G
Dimethyl Sulfoxide (DMSO)	G	G	G	G	G	G	G	G	P	P	P	P		F	P	P	P	P	G	G		
Dioxane	F	F	F	F	P	P	F	F	P	P	P	P	P	P	P	P	P	P	G	G	F	F
Ethyl Acetate (Ethyl Ester)	G	G	G	G	F	P	F	P	P	P	P	F	P	P	P	P	P	P	G	G		
Ethyl Alcohol (Absolute Ethanol)	G	F	G	G	G	G	G	F	G	G	G	F	G	F	P	P	P	P	G	G	G	G
Ethyl Chloride (Chloroethane)	F	P	P	P	F	P	F	P	P	P	P	F	P	P	P	P	P	P	G	G		
Ethylene Chloride	F	P	P	P	F	P	F	P	P	P	P	F	P	P	P	P	P	P	G	G	P	P
Ethylene Oxide (Pure)	F	F	F	F	F	P	F	P	P	P	F	P	P	P	P	P	G	G	G	G	G	G
Ethylene Oxide (Gas)	F	F	G	G	G	G	F	F	F	F	G	G	G	F	F	F	G	G	G	G		
Formaldehyde (Formalin) 40%	G	F	G	F	G	G	G	G	G	G	G	G	F	F	F	P	G	F	G	G	G	G
Formic Acid (50%)	F	F	G	G	G	F	G	G	G	P	F	F	F	F	G	F			G	G		
Formic Acid (100%)	F	F	G	G	G	F	G	G	G	P	F	P	F	F	P	P			G	G	P	P
Glycerine (Glycerol)	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	F	F
Hexane	P	P	F	F	F	F	F	P	G	P	F	P	P	P	G	G	G		G	G	G	G
Hydrobromic Acid (69%)	G	G	G	F	G	F	G	G	G	G	F	F	F	F	F	F	P	P	G	G	P	P
Hydrochloric Acid (5%)	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	F		G	G	P	P
Hydrochloric Acid (35%)	G	G	G	G	G	F	G	F	G	G	F	P	G	G	G	F	P	P	G	G	P	P
Hydrofluoric Acid (48%)	G	G	G	G	G	F	G	F	F	F	F	P	P	P	G	F	P	P	G	G	P	P
Hydrogen Peroxide (30%)	G	G	G	G	G	F	G	F	G	G	G	G	G	F	G	G	F	P	G	G	G	P
Lactic Acid (85%)	G	G	G	G	G	F	G	F	G	F	G	F	G	G	G	G	G	G	G	G	G	P
Methyl Acetate	G	P	G	G	G	F	F	P	P	P	P	P	P	P	P	P			G	G		
Methyl Alcohol (Methanol)	G	F	G	G	G	G	G	F	G	G	F	F	F	F	P	P	P	P	G	G	G	G
Methyl Ethyl Ketone (Butanone)	P	P	P	P	G	G	G	P	P	P	P	P	P	P	P	P	P	P	G	G		
Methylene Chloride (Dichloro Methane)	P	P	G	P	G	P	G	P	P	P	P	P	P	P	P	P	P	P	G	G	P	P
Mineral Oil	F	G	G	G	G	F	G	F	G	G	G	G	G	G	G	G	G	G	G	G	G	G
Nitric Acid (10%)	G	G	F	F	F	F	F	F	F	G	G	F	F	G	F	F	G		G	G	P	P
Nitric Acid (70%)	G	P	G	P	P	P	G	P	P	P	P	P	P	P	P	P	P	P	G	G	P	P
Nitrobenzene	P	P	P	P	P	P	G	P	P	P	P	P	P	P	P	P	P	P	G	G	F	P
Oxalic Acid (10%)	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G
Perchloric Acid (70%)	F	P	F	P	F	P	F	P	F	P	P	F	F	F	F	F	G	G	F	F	P	P
Phenol (100%)	P	P	P	P	P	P	G	P	P	P	P	P	P	P	P	P	P	P	G	G	P	P
Phosphoric Acid (85%)	G	P	G	G	G	F	G	F	G	F	G	F	G	G	G	F	P	G	G	G	G	P

	LDPE		HDPE		PP		PMP (TPX)		PVC		PC		ABS/PS		SAN		PMMA		PTFE		POM	
TEMPERATURE °C	20	50	20	50	20	50	20	50	20	50	20	50	20	50	20	50	20	50	20	50	20	50
Picric Acid	P	P	P	P	P	P	G	G	P	P	P	P	F	F	G	G			G	G		
Potassium Hydroxide (30%)	G	G	G	G	G	G	G	G	G	G	P	P	F	F	G	G	G	G	G	G	G	G
Potassium Permanganate	G	G	G	G	G	F	G	G	G	F	G	G	G	F	F	F	F	P	G	G	F	F
Propylene Glycol	G	G	G	G	G	G	G	G	F	F	F	F	G	G	G	G			G	G	G	G
Pyridine	P	P	P	P	G	G	F	P	P	P	P	P	P	P	P	P	P	P	G	G	G	F
Salicylic Acid	G	G	G	G	G	G	G	G	F	F	G	F	G	F	G	F			G	G	P	P
Silver Nitrate	G	F	G	G	G	G	G	G	G	F	G	G	G	G	G	F			G	G	F	F
Sodium Hydroxide (50%)	F	F	G	G	G	G	G	G	G	G	P	P	G	F	G	F	G	P	G	G	G	G
Sodium Hypochlorite (15%)	G	F	G	F	F	P	G	G	G	G	F	F	G	F	G	G	G	G	G	G		
Sulphuric Acid (20%)	G	G	G	G	G	G	G	G	G	F	G	F	G	G	G	F	G	G	G	G		
Sulphuric Acid (60%)	G	F	G	F	G	F	G	F	G	F	G	F	G	F	P	P	P	P	G	G	P	P
Sulphuric Acid (98%)	F	F	F	F	F	G	F	F	F	P	P	P	F	F	P	P	P	P	G	G	P	P
Tetrahydrofuran (THF)	F	P	F	G	F	F	F	F	P	P	P	P	P	P	P	P	P	P	G	G	F	F
Toluene	F	P	P	P	P	P	F	F	P	P	P	P	P	P	P	P	P	P	G	G	P	P
Trichloroacetic Acid	F	P	F	P	F	F	G	G	F	P	F	P	F	P	P	P			G	G		
Trichloroethylene	P	P	P	P	P	P	G	P	P	P	P	P	P	P	P	P	P	P	G	G	P	P
Turpentine	F	P	F	P	F	P	F	P	F	F	F	P	P	P	F	P	F	P	G	G	G	G
Xylene	P	P	F	P	P	P	F	P	P	P	P	P	P	P	P	P	P	P	G	G	G	G
Zinc Chloride (10%)	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	F
Zinc Sulphate (10%)	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	G	F	P

Physical Properties of Plastics

	LDPE	HDPE	PP	PMP (TPX)	PC	PS	PMMA	PTFE
Max Usage Temp. °C	80	120†	135	145	135	90	50	260
Min Usage Temp. °C	-50	-100	-20‡	-180	-135	-40	-60	-200
Transparency	Translucent	Translucent	Translucent	Clear	Clear	Clear	Clear	Opaque
Flexability	Flexible	Rigid	Rigid	Rigid	Rigid	Rigid	Rigid	Rigid
Specific Gravity	0.92	0.95	0.9	0.83	1.2	1.05	1.2	2.2
Microwavable	Yes	No	Yes	Yes	Marginal	No	No	Yes
Sterilization								
Autoclavable	No	No	Yes	Yes	Yes	No	No	Yes
Gas	Yes	Yes	Yes	Yes	Yes	Yes	No	Yes
Dry Heat	No	No	No	Yes	No	No	No	Yes
Radiation (Gamma)	Yes	Yes	No	No	Yes	Yes	Yes	Yes
Disinfection	Yes	Yes	Yes	Yes	Yes	Some	Some	Yes

† Polymer may be malleable at temperatures above 80°C if the product is under structural stress

‡ Material may become brittle at low temperatures

Note:

The information provided in the Chemical Resistance of Plastics chart is intended to provide guidelines for reference only. Observe the tolerated temperature ranges of plastic when autoclaving laboratory plasticware. Remove any stoppers, fittings or caps from the plasticware prior to autoclaving. Plastic vessels should be autoclaved separately from their closures and other fittings. Autoclaving vessels with their closures in place can lead to deformation and destruction of the vessels. All statements are advisory only and imply no liability on the part of Heathrow Scientific.

Sterilization:

- Autoclaving (121°C, 16 psi for 15 minutes)—Clean and rinse items with distilled water before autoclaving. (Always completely disengage thread before autoclaving.) Certain chemicals which have no appreciable effect on resins at room temperature may cause deterioration at autoclaving temperatures unless removed with distilled water beforehand.
- Gas—Ethylene Oxide, formaldehyde, hydrogen peroxide.
- Dry heat (160°C, 120 minutes)
- Disinfectants—Benzalkonium chloride, formalin/formaldehyde, ethanol, etc.
- Radiation—gamma irradiation at 25 kGy (2.5 MRad) with unstabilized plastic.

Product Page Index *by Item No.*

120000	20	120116	55	120262	98	120569	9	120730	50
120003	59	120117	55	120265	120	120570	9	120731	50
120004	59	120118	55	120312	58	120575	46	120732	50
120005	59	120148	8	120313	58	120576	104	120733	50
120006	59	120155	11	120314	58	120577	104	120734	50
120007	59	120172	40	120315	58	120578	104	120735	50
120008	59	120177	114	120316	106	120579	104	120738	10
120010	62	120178	114	120318	13	120581	6	120739	10
120011	62	120184	79	120327	13	120582	6	120740	106
120012	62	120185	79	120328	13	120584	10	120741	106
120015	6,7	120186	79	120329	13	120585	10	120747	25
120018	92	120187	79	120330	13	120592	104	120748	25
120021	35	120188	79	120331	13	120593	104	120758	121
120026	47	120189	79	120332	13	120595	11	120759	84
120032	75,94	120190	79	120333	13	120598	14	120760	84
120033	75,94	120191	79	120356	93	120599	25	120761	84
120034	75,94	120195	58	120365	27	120604	17	120762	120
120035	75,94	120196	107	120366	98	120610	8	120763	84
120036	75,94	120197	107	120367	98	120611	8	120764	84
120037	75,94	120201	75,94	120370	98	120617	14	120765	84
120038	75,94	120202	75,94	120371	98	120618	104	120766	84
120039	75,94	120203	75,94	120372	98	120622	6	120767	84
120040	75,94	120204	75,94	120375	40	120623	6	120768	84
120041	75,94	120205	83	120376	98	120624	6	120769	84
120042	75,94	120206	83	120379	98	120625	59	120770	85
120043	75,94	120207	83	120383	98	120626	58	120771	85
120044	75,94	120209	13	120387	98	120627	55	120772	85
120052	36	120210	13	120388	98	120637	51	120791	28
120056	47	120212	13	120391	98	120638	51	120792	28
120057	47	120214	13	120392	98	120639	51	120793	28
120070	91	120215	18	120398	83	120640	51	120813	28
120071	91	120216	18	120399	83	120641	51	120814	14
120072	91	120217	18	120422	13	120642	51	120815	11
120073	91	120218	18	120428	98	120643	51	120817	78
120074	91	120219	52	120460	15	120644	51	120818	78
120075	91	120220	97	120479	17	120645	51	120827	14
120076	91	120221	97	120480	45	120646	51	120834	38,86
120077	36	120222	57	120481	45	120647	51	120835	38,86
120078	36	120223	57	120487	16	120648	51	120836	57
120085	85	120224	57	120491	8,15	120649	51	120845	54
120086	85	120225	57	120499	45	120650	51	120846	54
120087	85	120226	57	120506	25	120651	51	120847	54
120088	85	120227	57	120507	25	120652	51	120848	54
120089	84	120228	20	120519	97	120653	107	120851	61
120091	40,85	120241	73	120521	15	120654	107	120852	61
120097	18	120243	73	120524	15	120656	59	120853	61
120098	18	120245	111	120527	6,7	120657	59	120854	61
120099	18	120246	111	120538	71,93	120661	21	120855	61
120100	18	120247	111	120539	71,93	120667	121	120856	61
120101	18	120248	111	120540	71,93	120668	121	120857	61
120102	18	120249	111	120541	71,93	120669	16	120858	61
120103	18	120250	112	120547	12	120682	6	120880	57
120104	18	120251	112	120548	12	120683	38	120881	57
120105	18	120252	112	120549	12	120684	38	120882	57
120106	18	120253	112	120550	12	120685	38	120884	14
120107	118	120254	112	120551	12	120686	38		
120108	34	120255	112	120552	12	120687	38		
120110	34	120256	112	120553	12	120690	51	HEA-----	
120111	34	120257	112	120554	19	120710	57	HEACRY-0500	123
120112	38,86	120258	117	120556	108	120722	121	HEACRY-1000(x)	123
120113	38,86	120259	113	120557	108	120727	71,90	HEACRY-2000(x)	123
120114	38,86	120260	113	120567	14	120728	71,90	HEACSEAL-75(x)	122
120115	38,86	120261	35	120568	9	120729	71,90	HEACSEAL-XL(x)	122
								HEALCRY-1(x)00	124

HEALCRY-1100(x)	124	HS2345(x)	64	HS20162A	46	HS24319	82
HEALCRY-1200(x)	124	HS2860 A, B, C	97	HS20341 B	32	HS24320 B, P	82
HEALCRY-1258(x)	124	HS2860 CB, CR, CY	97	HS20611	49	HS24322 A, E	38,86
HEALCRY-1380(x)	124	HS2860(x)	97	HS20612	49	HS24323 A, E	38,86
HEALCRY-1700(x)	124	HS2862 (x),(xx)	100,101	HS20612M	49	HS24330A	81
HEALCRY-2380(x)	124	HS2862 KB, KC	101	HS20613	49	HS24490 D, M	27
HEALCRY-2500	124	HS2864	99	HS20613 C, E	48	HS24645	27
HEALCRY-3020(x)	124	HS3000	16	HS20615M	49	HS24670	27
HEAMINI-1700	124	HS4282	95	HS20620	48	HS2538A	31
HEAMISL-1000(x)	126	HS4283P	95	HS206203	48	HS27501 A, B, C, D	83
HEAPSEAL-108(x)	122	HS4284A	72	HS206204	48	HS27502 A, B, C, D	83
HEARNBW-1(x)00	125	HS4323	66	HS20622 A, B	44	HS27503 A, B, C, D	83
HEARNBW-2(x)00	124	HS4323K	66	HS20623	44	HS27504 A, B, C, D	83
HEARNBW-3(x)00	126	HS4325	66	HS20629 A, B	44	HS27505 A, B, C, D	83
HEASIDE-1(x)00	130	HS4422	66	HS20630(x)	44	HS27511 A, B, C, D	83
HEASPOT-1(x)00	125	HS4423	66	HS20631B	44	HS27512 A, B, C, D	83
HEASPOT-2(x)00	125	HS4426 B, R	66	HS20632 A, B, C	44	HS27513 A, B, C, D	83
HEASPOT-3000	125	HS4427 B, R, S	66	HS20821A	50	HS27514 A, B, C, D	83
HEASPOT-4000	125	HS4428	66	HS23051	80	HS27515 A, B, C, D	83
HEASPOT-5000	125	HS4429	66	HS23052	80	HS27551 A, B, C, D	83
HEAT-SPOTS-(x)	125	HS5000RFFIL	16/17	HS23053A	80	HS27552 A, B, C, D	83
HEAT-SPOTS-44	125	HS5000RFRFA	16	HS23072	84	HS28270	120
HEAT-SPOTS-50(x)	125	HS5000RFRP	16	HS23108	84	HS28271	120
HEAT-SPOTS-75	125	HS5360 A, B	21	HS23202 A, N	67	HS28272	120
HEAT-SPOTS-100	125	HS6594	21	HS23205	96	HS28280	120
HEATS-TNY	125	HS8151	63	HS23206	96	HS28622	100
HEATTGP-1000	126	HS8171 A, B	63	HS23210	96	HS28630	102
HEATT-SWALL	126	HS8610 A, B, D, E	26	HS23215	84	HS28714K	88
HEATT-TNY	126	HS8881 A, B, C, D	52,113	HS23224	85	HS28721G	88
HEATTL(x)-1000	126	HS8882 A, B, C, D	52,113	HS23229B	96	HS28724B	88
HEATTLW-1050(x)	126	HS9902	107	HS23243A	29	HS28729R	88
HEATTLW-2016(x)	126	HS9903	108	HS232461	85	HS29022(x)	77
HEATTS(x)-1000	126	HS9904	108	HS23263 A, B	73	HS29025(x)	73
HEATTSW-2240(x)	126	HS9913	107	HS23271A	78,95	HS29040 A, B	76
HS-----		HS9914	107	HS23300D	99	HS29041A	72
HS1000RFU	17	HS9930	107	HS23325A	99	HS29050A	76
HS1002 A, B, C	114	HS10044H	114	HS23381B	99	HS81121 A, B, C, D	62
HS1040 A, B, C	115	HS10050	7	HS23401	31	HS86655	63
HS1041	115	HS10060	67	HS23402	31	HS100011	8
HS1042	115	HS10320	114	HS23403	31	HS100012	8
HS1419 A, B, C	57	HS10321	114	HS23404	31	HS100502	6,7
HS1420 A, B, C	57	HS10322	114	HS23411	116	HS100503	7
HS1420 AA, BB, CC	57	HS10323	114	HS23412	116	HS100504	7
HS1421 A, B, C	57	HS14521A, B, C	55	HS2344A	70	HS100506	7
HS1422	57	HS14522	55	HS23452	64	HS159890B	106
HS1423 A, B, C	57	HS15092D	21	HS23453	32	HS206371 A, B, C	45
HS1423 AA, BB, CC	57	HS15092ZP	26	HS23457	117	HS206372 A, B, C	45
HS1424 A, B, C	57	HS15094	26	HS23458	117	HS206373 A, B, C	45
HS1424AA	57	HS15906	58	HS23460	126	HS234521	116
HS14251 A, B, C, D	57	HS15907	58	HS23461 A, B, C, D	70	HS234522	116
HS1425 A, B, C	57	HS15909	58	HS234631	48	HS234523	116
HS1425AA	57	HS15951A	109	HS234634	48	HS234524	121
HS1426 A, B, C	57	HS15979A	105	HS23471 A, B, C	117	HS234525 B, C, D, R	121
HS1426AA	57	HS15981 A, B, C	109	HS23491A	117	HS234526 A, B, C	120
HS1427 A, B, C	57	HS15982	108	HS23500	30	HS243071 B, R, W, Y	82
HS1427AA	57	HS15983 A, B, G	108	HS23502	30	HS243073 B, R, W, Y	82
HS2132A	89	HS15986	108	HS23503	30	HS243077 B, R, W, Y	82
HS2134A	89	HS15988H	106	HS23504	30	HS300502	16
HS2165 A, B, C	89	HS15989(x)	105	HS23534A	122		
HS2166	89	HS15990 A, B, C	106	HS23536A	122		
HS2200 A, B	34	HS15991 A, B, C	106	HS24306	81		
HS2201 A, B	34	HS15994(x)	105	HS24311 A, B	29,86		
HS2202	34	HS15996 A, B, C	105	HS24312 A, B	29,86		
		HS15998C	58	HS24318	82		

Product Alpha Index

A

ABS Dispenser for Parafilm® M	121
Acrylic Dispenser for Parafilm® M	121
Adapt-a-Rack® Tube Racks	79
Aluminum Weighing Boats & Dishes	55
Arctic Squares® Cryostorage Boxes	99
Arctic Squares® Upright Freezer Rack	99
Autoclavable Bags	114

B

Bag Sealer	120
Baskets, Wire	32
Biohazard Disposal Bags/Holder	114
Benchtop Counters	21
Benchtop Dispensing Bins	116
Blood Collection Tray	35
Box Top Dispenser for Parafilm® M	121
Bulbs	44

C

Cardboard Slide Tray/Clasps	107
Cargo Pipette Rack	49
Cell Spreaders	63
Centrifuges	6-9
Chillers, PCR	50
Clamp-ette® Pipette Holder	46
Clinical Tube Racks	38-39
Clinical 4-Way Tube Rack	38
Clocks	27
Coated Wire Tube Racks	84-85
Colony Counter, eCount	20
Conical Centrifuge Tubes	66
ContainerSEAL™	122
Cooling Containers	88-91
Counters	20-21
Cryo-Babies® Labels	124
CryoTags® Labels	124
Cryogenic Racks & Storage	96-99
Cryogenic Vial Boxes & Partitions	97
Cryogenic Vial Canes & Sleeves	96
Cryogenic Vial Rack	96
Cryogenic Vial Tube Storage Box	97
Cryogenic Vials	67
Cube Tube Rack	76

E

Earplug Dispenser	115
Eco smartFunnel™	52
eCount Colony Counter	20

D

Digital Orbital Shaker	15
Dispenser, Label	126
Dispenser, Parafilm®	121
Dispenser, Petri Dish	64

Dispensers, PPE	115-118
Droplet Blood Collection Tray, Lid	34
Droplet Corrugated Storage Tray	40
Droplet Sample Storage Tray	40
Drying Racks	28-29
Dual Solution Reservoir	50
DuraPorter® Transport Box	36-37

F

4-Way Tube Rack	77
Film, Sealing	120
Floating Tube Racks	89
Flip & Grip® Pipette Holder	47
Foam Tube Racks	95
Foil Rolls	122
Fold and Snap Tube Rack	82
Forceps	58
Freezer Boxes	97-99
Freezer Boxes, PP Thin-Film	98
Freezer Racks, Upright, Chest	100-101
Freezer Thermometer	102
Funnels	52

G

Glove Box Holders	117-118
Gusto® High Speed Mini Centrifuge	7

H

Hand Tally Counter	21
Heat Sealer	120

I

Ice Buckets	see Cooling Containers
Inoculating Needle/Loops	62
Inoculation Picks	62

L

Lab Alert Clocks & Timers®	27
Labels	123-126
Label Station	126
Lab Drying Rack	28
Lab Mats	24-25
Lab Wipe Holders	116
Laboratory Notebooks	26
Lead Rings	52, 113
Light Box	21
Low-Temp Microscope Slide Box	106
Low-Temp PCR Rack: 96-Well	71, 93

M

MagFuge® Centrifuge/Stirrer	6
Magnetic Lab Storage (MagLab)	30
Magnetic Induction Stirrer	10
Magnetic Stirrers	6, 10-11
Magnifier	21
Manual Pipette Controller	17

Mega Rack®	38, 86
Markers, Permanent	26
Microcentrifuge	See Centrifuges
Microcentrifuge Tube Racks & Storage	72-78
Microcentrifuge Tubes	66
Microtube Racks	73
Microtube Racks: 80-Well	73
Microtube Storage Boxes, Cardboard	97
Microtube Storage Boxes 50, 81, 100-Well	75
Mini Coolers	91
Mini Light Box	21
Mini Magnifier	21
Mini Stirrer	11
Mini Vortexer	14
Mixers	13-15

N

Notebooks	26
-----------	----

O

One-Well Connecting Tube Racks	80
One-Well Friction-Fit Tube Rack	80
One-Well Wire Tube Rack	80
OneRack® Tube Rack	83
OneRack® Multi Tube Rack	78
Organizers	30-32

P

Parafilm® M Sealing Film	120
Parafilm® Dispensers	121
PCR Cooler: 96-Well	71, 90
PCR Tube & Storage Racks	70-72, 90
PCR Tube Racks: 96-Well	70
PCR Tube Racks Reversible	70
Pearl® Pipettes	18
Peg Tube Racks	86
Petri Dish Dispenser	64
PetriSEAL™	122
Phlebotomy Trays	34-35
Picks	62
Pipette Controllers	17
Pipette Fillers	44
Pipette Racks	49
Pipette Stands	48
Pipette Stands & Holders	45-49
Pipette Starter Kit	19
Pipette Tip Disposal Box	48
Pipette Workstation	48
Pipettes	18, 45
Pop-Up Tube Racks	82
Portable Centrifuge Kit	9
Cryogenic Vial Boxes & Partitions	97
Premium Microscope Slide Boxes	105
Premium Plus Slide Box: 100-Place	104

R

Racks	
- Freezer Racks	<i>See Freezer Racks</i>
- Tube Racks	<i>See Tube Racks</i>
Razor Blade	113
Reagent Reservoirs	51
Reservoir Chiller	50
Reservoirs	50-51
Reversible Tube Racks: 96-Well	72
RF1000	17
RF3000®	16
Rota-Rack® Duo Tube Racks	72
Rota-Rack® Tube Racks	76

S

Safety Glasses Dispensers	115
Sample Collection Trays	34-35
Sample Storage Trays	40
Sample Transport Box	36-37
Sampling Utensils	55-58
Scoops	59-61
Screw-top Tubes	67
Sealing Tape	120, 122
Shakers	15
Sharpie® Pens	26
Sharps Containers	114
Silicone Lab Mat	25
Slide Boxes	104-106
Slide Holders	109
Slide Mailers	108
Slide Staining Tray Set	109
Slide Trays	107
smartBoats™	58
smartPans™	55
SmartPicks™	62
smartSlicer™	113
smartSpatula™	59
Snap Together Conical Tube Rack	81
Soft Cover Dispenser	115
Spatulas	58
Sprout® Plus Mini Centrifuge	8
Standard Microscope Slide Boxes	106
Stir Bar Retriever	12
Stir Bars	12
Stirrers	<i>See Magnetic Stirrers</i>
Stopwatches	27
Storage Rack/Microscope Slide Boxes	106
Storage Rack/Microtube 79-Well Racks	73

T

3-in-1 Dispenser	115
Test Tube Baskets	32
Thermometer	<i>See Freezer Thermometer</i>
Timers	27
Totemax® Blood Collection Tray	35
Tough-Spots® Labels	125
Tough-Tags® Labels	126

Tourniquets, X-tourn®	34
TracRack® Tube Rack/Box	93
Transfer Pipettes	45
Transformer Cube Tube Rack	78, 95
Tubby® Storage Container	32
Tube Racks	38-39, 70-86
Tubes	66-67
Tube Storage	70-86, 92-99
Tube Storage Boxes & Partitions, 15/50mL	97
Two Tier Microtube Rack: 20-Well	72

U

Universal Carousel Pipette Stand	45
Universal Pipette Holder	46
Universal Tube Rack	81
Urine Container Racks & Trays	40

V

Vacuum/Pressure Bulb	44
Vials	67
Vortex Mixers	13-14
Vortexer	13

W

Wash Bottles	111-112
Weighing Boats & Dishes	55-58
Weighing Papers	55
Wire Drying Rack	29
Wire Tube Racks	84-85
Work 2 Store® Tube Storage	92
Workstation Dispensing Bins	31
Workstation Lab Mat	25
Workstation Organizer	31
Workstation Storage Bin	116

